







Item # Product Name	★ STANDARD ★			★★ ENHANCED ★★	★★★ PREMIUM ★★★	
	15019	15076	15670	15770	15870	15970
	Therapeutic 5 Zone	Quick 'N Easy Comfort	Gravity 6 Mattress	Gravity 7 Mattress	Gravity 8 Mattress	Gravity 9 Mattress
						
Dimensions	36" x 80" x 6"	36" x 80" x 6"	36" x 80" x 6"	36" x 80" x 6"	36" x 80" x 6"	36" x 80" x 6"
Key Differences	Single layer 5 zone pressure redistribution	14" x 14" x 37.5" box footprint	Single layer 5 zone pressure redistribution with heel slope	Dedicated heel section	Channel cut heel slope	Visco foam top layer
Additional Sizes Available (Item#: WxLxH)	—	—	—	15876: 36" x 76" x 6" 15784: 36" x 84" x 6" 15777: 36" x 80" x 6" (w/raised rails) 15785: 36" x 84" x 6" (w/raised rails) 15877: 36" x 76" x 6" (w/raised rails)	15886: 36" x 76" x 6" 15884: 36" x 84" x 6" 15888: 36" x 80" x 6" (w/raised rails)	15996: 36" x 76" x 6" 15984: 36" x 84" x 6" G9-39WX-80LX6H: 39" x 80" x 6" 15970-4280: 42" x 80" x 6" 15970-4284: 42" x 84" x 6" 15977: 36" x 80" x 6" (w/raised rails) 15970-4280-RR: 42" x 80" x 6" 15970-4240-RR: 42" x 80" x 6"
Layers of Foam	1	1	1	3	3	4
Weight Capacity	350 lbs	275 lbs	350 lbs	350 lbs	400 lbs	450 lbs
Available with Raised Rail	—	—	—	Y - Standard	Y - Standard	Y - Standard
Dedicated Heel Protection Zone	—	—	—	Y	Y	Y
Heel Slope	—	—	Y	—	Y	Y
Non-Prorated Warranty	3 year foam, 1 year cover	1 year limited	3 year foam, 3 year cover	7 year foam, 3 year cover	10 year foam, 3 year cover	15 year foam, 3 year cover
HCPCS Code	E0184	E0272	E0184	E0184	E0184	E0184
NEW! Updated Features!	<ul style="list-style-type: none">• Softer and longer-lasting cover fabric• Breathable stretch cover with welded seams for better water resistance• Intuitive tie-down straps easily fasten to virtually any bed frame for a secure fit to reduce mattress migration and falls• Full length articulation cuts to extend mattress life cycle and improve conformity to all bed positions• Newly added handles make setup, delivery, and cleaning easier	<ul style="list-style-type: none">• Softer and longer-lasting cover fabric• Breathable stretch cover with welded seams for better water resistance• Upgraded handles make setup, delivery, and cleaning easier• Intuitive tie-down straps easily fasten to virtually any bed frame for a secure fit to reduce mattress migration and falls• Full length articulation cuts to extend mattress life cycle and improve conformity to all bed positions	<ul style="list-style-type: none">• Softer and longer-lasting cover fabric• Intuitive tie-down straps easily fasten to virtually any bed frame for a secure fit to reduce mattress migration and falls• Full length articulation cuts to extend mattress life cycle and improve conformity to all bed positions• Newly added handles make setup, delivery, and cleaning easier	<ul style="list-style-type: none">• Therapeutic die-cut foam spans the full width for maximum patient contact• Full length articulation cuts to extend mattress life cycle and improve conformity to all bed positions• High-density foam perimeters remain beneath for safe patient transfers• Visco foam heel section offloads pressure from vulnerable heels• Newly added handles make setup, delivery, and cleaning easier	<ul style="list-style-type: none">• Therapeutic die-cut foam spans the full width for maximum patient contact• Full length articulation cuts to extend mattress life cycle and improve conformity to all bed positions• High-density foam perimeters remain beneath for safe patient transfers• Heel slope with channel cut foam reduces pressure on patient heels• Newly added handles make setup, delivery, and cleaning easier	<ul style="list-style-type: none">• Head to toe hybrid visco memory foam top layer• Therapeutic die-cut foam spans the full width for maximum patient contact• Full length articulation cuts to extend mattress life cycle and improve conformity to all bed positions• High-density foam perimeters remain beneath for safe patient transfers• Heel slope with channel cut foam reduces pressure on patient heels• Newly added handles make setup, delivery, and cleaning easier



Item # Product Name	SPECIALTY	BARIATRIC	
	6500-DE-2-FB Dynamic Elite	15301 Bariatric Foam Mattress	1500SC-42-FB Shearcare 1500
			
Dimensions	36" x 80" x 6"	42" x 80" x 6.5"	42" x 80" x 6"
Key Differences	Gel infused visco foam top layer	Durable bariatric foam	Durable bariatric foam
Additional Sizes Available (Item#: WxLxH)	—	15301: 42" x 80" x 6.5" 15301-84: 42" x 84" x 6.5" 15310-84: 48" x 84" x 6.5" 15312: 54" x 80" x 6.5"	1500SC-48-FB: 80" x 48" x 6"
Layers of Foam	4		2
Weight Capacity	400 lbs	585 lbs (15301, 15301-84) 665 lbs (15310-84, 15310) 750 lbs (15312)	650 lbs (48" W-750 lbs)
Available with Raised Rail	Y - Custom	—	Y - Custom
Dedicated Heel Protection Zone	Y	—	Y
Heel Slope	Y	—	—
Non-Prorated Warranty	15 year foam, 3 year cover	1 year limited warranty	5 year foam, 3 year cover
HCPCS Code	E0184	E0184	E0184
NEW! Updated Features!	<ul style="list-style-type: none">• Therapeutic die-cut foam spans the full width for maximum patient contact• Full length articulation cuts to extend mattress life cycle and improve conformity to all bed positions• High-density foam perimeters remain beneath for safe patient transfers• Newly added handles make setup, delivery, and cleaning easier	<ul style="list-style-type: none">• Softer and longer-lasting cover fabric• Breathable stretch cover with welded seams for better water resistance• Intuitive tie-down straps easily fasten to virtually any bed frame for a secure fit to reduce mattress migration and falls• Full length articulation cuts to extend mattress life cycle and improve conformity to all bed positions• Newly added handles make setup, delivery, and cleaning easier	<ul style="list-style-type: none">• Softer and longer-lasting cover fabric• Breathable stretch cover with welded seams for better water resistance• Intuitive tie-down straps easily fasten to virtually any bed frame for a secure fit to reduce mattress migration and falls• Full length articulation cuts to extend mattress life cycle and improve conformity to all bed positions• Newly added handles make setup, delivery, and cleaning easier

LEGEND		
STANDARD ★ Entry level systems suitable to meet basic user needs for a good value.	ENHANCED ★★ Meets basic needs plus supplemental features to improve the user's experience.	PREMIUM ★★★ Includes multiple advanced features for the ultimate user experience.

Discover Why Our Therapeutic Support Surfaces are Leaders in Homecare & Long Term Care.

Drive DeVilbiss Healthcare’s newly improved foam mattresses offer enhanced comfort and safety for both patients and caregivers. With clinical features that meet a wide range of demands to prevent and treat pressure injuries, along with practical upgrades like durable handles and water-resistant covers with welded seams, our mattresses deliver both superior care and long-lasting performance.

