



# 2025

## Canes & Crutches Catalog

Enhancing quality  
of life is our  
driving force





# Enhance Independence Safely

For dependable mobility products you can trust, doctors, hospitals, and individuals alike turn to Drive DeVilbiss Healthcare. We offer a wide variety of canes and crutches to provide the precise, individual level of support you need to protect your independence and safely enhance your mobility. Whether you're looking for mobility assistance during a temporary period of recovery, or to better support your balance during everyday activities, you'll find a finely crafted solution from Drive.

## Table of Contents

Canes	3-16
Cane Accessories	17-24
Crutches	25-28
Crutch Accessories	28-30




Hurrycane®  
Freedom Edition™

Item #	Color
HCANE-BK-C2	Black
HCANE-BL-C2	Blue
HCANE-RD-C2	Red
HCANE-PR-C2	Purple
HCANE-PB-C2	Pastel Blue
HCANE-PP-C2	Pastel Pink
HCANE-PG-C2	Pastel Green

- Freestanding design
- Folds & unfolds in seconds
- 8 convenient height adjustments
- Whisperflex pivoting base design for stability & balance
- Steadigrip base for all-terrain traction
- Comfort-fit shock-absorbing handle
- Anodized aluminum shaft
- Wrist strap included

Shop Now



HurryCane® Go™

Item #	Color
HCANE-BK-G2	Black

- Stands on its own—Even on Carpet The wider SteadiGrip™ base with cross-hatched tips provides greater stability and traction to help your cane stand on any surface
- Improved Stability with a wider range of motion, the pivoting WhisperFlex™ design provides even better balance and stability
- Enhanced Traction with new cross-hatched tips provide maximum traction
- Adjustable height with 8 convenient height adjustments from 30.5"–37.5" give you a perfect fit

Shop Now



## HurryCane® Accessory Feet

### Item #

HCFEET-BK

- Large cross-hatched tips give you even greater traction and help the HurryCane® Freedom Edition stand on its own even on carpet!
- For use only with the new HurryCane® Freedom Edition™. Not compatible with previous editions with rounded tips.

[Shop Now](#)

## HurryFlex® Accessory Cane Tip

### Item #

HTIP-BK

- Transforms Single-Point Canes: Add this replacement cane tip to your single-point cane to turn it into a self-standing HurryCane.
- Freestanding Design: Self-standing cane tip lets your cane stand on its own so it's always by your side when you need it
- Enhanced Stability: Tripod cane base features three points of contact for all-terrain traction, plus a pivoting base for additional stability
- Lightweight and Compact: These rubber cane tips are 54 percent lighter and 80 percent smaller than large, quad cane tips, while offering incredible support
- Tool-Free Installation: The 3-point cane tip installs in seconds, no tools necessary

[Shop Now](#)



## Gel Grip, Aluminum Offset Canes, Height Adjustable


### Item #

RTL10372BK  
RTL10372BC  
RTL10372RC  
RTL10372FL  
RTL10336MA  
RTL10336MN  
RTL10336SF  
RTL10372GL  
RTL10372PL

### Color

Black  
Blue Crackle  
Red Crackle  
Floral  
Marble  
Mech. Navy  
Spring  
Green Leaves  
Plaid

- Ergonomically designed handle with soft Gel Grip provides comfort and security
- Manufactured with sturdy, extruded aluminum tubing
- Comes standard with premium Gel Grip and wrist strap
- Easy-to-use, one-button height adjustment with locking ring. Locking ring prevents rattling
- Handle height adjusts from 30" to 39"
- Weight Capacity 300 lbs. / Approx. user height 5'-6'5"

[Shop Now](#)


## Offset Handle Aluminum Canes

### Item #

RTL10306  
RTL10307  
RTL10303NV  
RTL10303AF  
RTL10303BF  
RTL10303GP  
RTL10303HR  
RTL10303LV  
RTL10303NA  
RTL10303PF  
RTL10303PN  
RTL10303SN  
RTL10303SW

### Color

Black  
Bronze  
Navy  
Amber Floral  
Burst Floral  
Green Plaid  
Vibrant Holiday Red  
Lavendar  
Navy Argyle  
Pink Floral  
Leafy Green Pine  
Frosted Snowflake  
Swirl

- Ergonomic soft Foam Grip provides comfort and security
- Manufactured with sturdy, extruded aluminum tubing
- Easy-to-use, one-button height adjustment with locking ring. Locking ring prevents rattling
- Handle height adjusts from 30" to 39"
- Weight Capacity 300 lbs. / Approx. user height 5'-6'5"
- 3/4" Cane Diameter

[Shop Now](#)


# Canes

## Offset Aluminum Canes with Tab-Loc Silencer

**Item #**

10381BLK-6

**Color**

Black

- Tab-Loc silencer provides added safety to ensure a locked height adjustment
- Features a padded grip
- The flat side of the grip faces up to reduce stress that the user's hand may experience
- Ergonomically designed padded handle
- Comes standard with wrist strap
- Manufactured with sturdy, extruded aluminum tubing
- Handle height adjusts from 28.75" to 37.75"

[Shop Now](#)

## Offset Handle Aluminum Canes

**Item #**

10309BKW-6

**Color**

Black Wave

10303-6

Silver

- Ergonomically designed handle with soft Foam Grip provides comfort and security
- Comes standard with foam rubber grip and wrist strap
- Manufactured with sturdy, extruded aluminum tubing
- Easy-to-use, one-button height adjustment with locking ring
- Locking ring prevents rattling
- Comes with Black wave, Blue wave, and Bronze wave
- Handle height adjusts from 30" to 39"

[Shop Now](#)

## Bariatric Offset-Handle Canes

**Item #**  
10318-6

**Color**  
Black

- Manufactured with steel tubing
- Easy-to-use, dual-button height adjustment with locking ring Locking ring prevents rattling
- Comfortable vinyl hand grip with wrist strap



[Shop Now](#)



## Comfort Grip Cane



**Item #**  
RTL10336RG  
RTL10336FG  
RTL10336GR  
RTL10336AN  
RTL10336PA  
RTL10336LM  
RTL10336MA  
RTL10336MN  
RTL10336SF

**Color**  
Rose Gold  
Forest Green  
Graphite  
Anchor  
Patriotic  
Lemons  
Marble  
Mechanical  
Spring Floral

- Manufactured with sturdy, extruded aluminum tubing
- Comes Standard with premium wrist strap
- Easy-to-use, one-button height adjustment with locking ring Locking ring prevents rattling
- Handle height adjusts from 30" to 39"
- Approx. user height 5'-6'5"



[Shop Now](#)



# Canes

## Aluminum Folding Canes, Height Adjustable


### Item #

RTL10304BZ  
RTL10304  
RTL10304BW  
RTL10304PN  
RTL10304DG  
RTL10304GD  
10304BF-1

### Color

Bronze  
Black  
Black/White Abstract  
Pink Abstract  
Puppy  
Garden  
Black Floral

- Cane folds into four convenient parts for easy storage
- Manufactured with sturdy, extruded aluminum tubing
- Easy-to-use, one-button height adjustment with locking ring Locking ring prevents rattling
- Handle height adjusts from 30" to 37"
- Approx. user height 5'-6'5"

[Shop Now](#)

## Bariatric Aluminum Folding Cane, Height Adjustable

### Item #

RTL10304HD

### Color

Black

- Cane folds into four convenient parts for an easy storage
- Attractive wood handle with brass collar
- Manufactured with sturdy, extruded aluminum tubing
- Comes with plastic clip to hang cane when folded
- Handle height adjusts from 33" to 37"

[Shop Now](#)



# Canes

## Aluminum Folding Canes with Gel Grip, Height Adjustable

### Item #

RTL10370RC

RTL10370BC

RTL10370BK

### Color

Red Crackle

Blue Crackle

Black

- New gel grip reduces hand stress and fatigue
- Aluminum folding cane, height adjustable
- Cane folds into four convenient parts for easy storage
- Manufactured with sturdy, extruded aluminum tubing
- Comes with plastic clip to hang cane when folded
- Available in Red Crackle, Blue Crackle, Black, and Hearts
- Handle height adjusts from 32" to 36"

[Shop Now](#)

## Folding Canes With Silicone Glow Grip Handle & Tip

### Item #

RTL10304CE

### Color

Celebration

- Handle and tip glow in the dark
- Handle height adjusts from 33" to 37"
- Manufactured with sturdy, anodized, extruded aluminum tubing
- Metal (lock nut) silencer
- Soft, silicone handle for added comfort

[Shop Now](#)

# Canes

## Flex N Go Cane

**Item #**  
RTL10305

**Color**  
Black



- Cane base with 3 feet, pivots in any direction providing added security and traction
- Soft grip handle relieves stress and pressure on fingers and wrist
- Newly designed tip is removable without tools, for easy replacement
- Height adjusts for use by most individuals
- 3 point base allow the cane to stand on its own
- Tip flexes and swivels to maintain contact with the ground when the cane is at an angle, providing superior balance
- Conveniently folds for easy storage

[Shop Now](#)



## Quad Canes, Small Base

**Item #**  
10301-4

**Color**  
Chrome

**Description**  
Vinyl Contoured Grip



- Four point base, combined with offset handle provides additional stability and support
- Easy-to-use, one-button height adjustment with locking bolt, allows quick, safe fitting
- No assembly required
- Handle height adjusts from 30" to 39"

[Shop Now](#)



# Canes

## Quad Canes, Small Base

Item #	Color	Description
RTL10310	Black	Foam Rubber Grip



- Four point base, combined with offset handle provides additional stability and support
- Easy-to-use, one-button height adjustment with locking bolt, allows quick, safe fitting
- No assembly required
- Handle height adjusts from 30" to 39"

[Shop Now](#)



## Bariatric Quad Canes

Item #	Color
10316-4	Chrome



- Lightweight, extra-strong base
- Manufactured with aluminum tubing, steel base
- Comfortable vinyl, contoured hand grips
- Four point base, combined with offset handle provides additional stability and support
- Aluminum handle, reinforced steel base
- Handle height adjusts from 30" to 39"

[Shop Now](#)





## Quad Canes, Large Base

Item #	Color	Description
10300-4	Chrome	Vinyl Contoured Grip

- Four point base, combined with offset handle provides additional stability and support
- Easy-to-use, one-button height adjustment with locking bolt, allows quick, safe fitting
- No assembly required
- Handle height adjusts from 30" to 39"



[Shop Now](#)

## Folding Quad Cane, Small Base

Item #	Color
RTL10310FCH	Charcoal

- Folds and disassembles for small footprint and easy storage
- Made of sturdy, extruded aluminum tubing
- Four-point base and offset handle for additional stability and support
- Contoured handle for comfortable grip
- Handle height adjusts with the push of a button from 30" to 39"
- Locking ring eliminates rattling and provides added safety



[Shop Now](#)

## Bariatric Quad Canes

Item #	Color	Description
10317-4	Chrome	Vinyl Contoured Grip

- Lightweight, extra-strong base
- Manufactured with aluminum tubing, steel base
- Comfortable vinyl, contoured hand grips
- Four point base, combined with offset handle provides additional stability and support
- Aluminum handle, reinforced steel base
- Handle height adjusts from 30" to 39"



[Shop Now](#)



## Round-Handle Cane with Foam Grip

Item #	Color
RTL10342	Chrome

- Manufactured with sturdy, extruded aluminum tubing
- Easy-to-use, one-button height adjustment with locking ring
- Locking ring prevents rattling
- Handle height adjusts from 28.75" to 37.75"



[Shop Now](#)



## Round-Handle Aluminum Cane

**Item #**


10302-6

**Color**

Chrome

- Comfortable, contoured handle with vinyl comfort grip
- Smaller handle than a standard cane to accommodate small hands
- Attractive finish resists scratches
- Manufactured with sturdy, extruded aluminum tubing
- Easy-to-use, one-button height adjustment with locking ring
- Locking ring prevents rattling
- Handle height adjusts from 28.75" to 37.75"

[Shop Now](#)



## Blind Folding Cane

**Item #**

10352-1

- High quality four section aluminum construction
- Shaft covered with white and red reflective tape for night visibility
- Cane length: 45.75"
- Comes standard with wrist strap and reinforced nylon tip

[Shop Now](#)





## Cane/Sling Seat

### Item #

RTL10360

- Provides a comfortable seat to rest on when open, and a sturdy quad cane when closed
- Manufactured with sturdy, extruded aluminum tubing
- Four sturdy legs with vinyl contoured tips
- Nylon sling bench, 15" x 5.75"
- Foam grip handles
- 34" tall when closed
- Seat to floor height: 20.5"
- Weighs 2.5 lbs



[Shop Now](#)

## Cane Seat

### Item #

RTL10365-ADJ

- Provides a comfortable seat to rest on when open, and a sturdy support cane when closed
- Manufactured with sturdy, extruded aluminum tubing
- Tripod design with vinyl contoured tipped legs
- 9" diameter seat
- 33" tall when closed
- Weighs 2.5 lb



[Shop Now](#)

# Walking Poles & Cane Accessories

## Drive Summit™ Walking Poles



### Item #

RTL10330MN

- Push-button locking system, designed with an easy-to-use push-button locking system
- 200 lb weight-bearing capacity per pole
- Lightweight yet sturdy aluminum shaft
- Ergonomic handle with rubber and foam grip
- Removable Wrist Strap
- Use a carry strap to attach the two poles together
- Bell-shaped tips offer greater stability, off-loading and traction
- Angled tips help users push off and propel themselves.
- Steel tip provides traction on nature trails, slippery or soft terrain and beaches
- Handle height adjusts from 37.75" to 47.75", Collapses to 30" or disassembles for storage and travel

[Shop Now](#)



## Cane Tips

### Item #

10383BK-8

- For use with manufacturers canes/crutches with a 1" diameter

[Shop Now](#)



# Cane Accessories

## Glow in the Dark Cane Tip ( $\frac{3}{4}$ " cane diameter)

**Item #**

RTL10324BB

- For use with Drive cane model #'s: 10302, 10303, 10304, 10306, 10307, 10309, 10317, 10331, 10335, 10339, 10342, 10345, 10347, 10350, 10351, 10370, 10372, 10381 series

[Shop Now](#)



## Cane Tip ( $\frac{3}{4}$ " cane diameter)

**Item #**

RTL10322BKB

- For use with Drive cane model #'s: 10302, 10303, 10304, 10306, 10307, 10309, 10317, 10331, 10335, 10339, 10342, 10345, 10347, 10350, 10351, 10370, 10372, 10381 series

[Shop Now](#)





# Cane Accessories

## Replacement Tips 1 1/8"

### Item #

RTL10386GB

- For use with manufacturers products with 1/18" diameter tubing (item # 10386G-8) and 7/8" tubing (Item # 10390G-8)
- Compatible replacement tip for the 11117N-1 Commode

[Shop Now](#)



## Replacement Tips 7/8"

### Item #

RTL10390GB

- For use with manufacturers products with 7/8" diameter tubing (item # 10386G-8) and 7/8" tubing (Item # 10390G-8)
- Compatible replacement tip for the 11117N-1 Commode

[Shop Now](#)



# Cane Accessories

## Small Base Quad Cane Tip (1/2" cane diameter)

**Item #**

RTL10320BK

- For use with Drive cane model #'s: 10310, 10311, 10312, 10315, 10384, 13078 series

[Shop Now](#)



## Large Base Quad Cane Tip (5/8" cane diameter)

**Item #**

RTL10321BKB

- For use with 10308, 10314, 10385 quad canes

[Shop Now](#)



# Cane Accessories

## Quad Cane Tip

### Item #

RTL10353

- Allows cane to stand upright on its own
- Provides better balance and stability with a 4-pronged base
- For use with  $\frac{3}{4}$ " diameter cane
- For use with Drive cane model #'s: 10302, 10303, 10304, 10306, 10307, 10309, 10317, 10331, 10335, 10339, 10342, 10345, 10347, 10350, 10351, 10370, 10372, and 10381 series

[Shop Now](#)

## Quad-Support Cane Tip

### Item #

RTL10351

- Replacement tip extends the life of your cane
- 54% lighter and 80% more compact than a large base quad cane
- Fits canes with  $\frac{3}{4}$ " diameter
- Provides superior grip and helps absorb impact shock throughout the walking motion
- Allows your cane to stand independently
- Engineered for maximum stability and traction over all varieties of terrain, Rubber
- Provides better balance and stability with a 4-pronged base

[Shop Now](#)



# Cane Accessories

## A-Frame Cane Holder

**Item #**

730-947

- Conveniently keeps your cane in place while you use your rollator
- Attaches to A-frame rollators only (eg: Elite/Adventure8/Comfort-Plus)

[Shop Now](#)



## Cane Holder and Strap

**Item #**

RTL102681

- Allows your cane to hang from any table in your home, restaurants and when traveling
- Easily attaches to most canes

[Shop Now](#)



# Cane Accessories

## Cane Strap

### Item #

730-949

- Attaches to any standard type cane
- Manufactured from durable nylon cord with a metal connecting ring and clip
- Provides the user with convenient control of the cane



[Shop Now](#)

## Cane Holder for Side-Fold

### Item #

730-946

- Allows you to safely and conveniently store your cane on your rollator, wherever you go
- Compatible with Airgo/Hugo side-folding rollators (except Sidekick)
- Tool-free installation



[Shop Now](#)

# Cane Accessories

## Cup Holder for Side-Fold

**Item #**

730-944

- Allows you to safely store and carry your beverage with you, wherever you go
- Swivels from left to right for maximum convenience
- Senior friendly tool-free installation
- Compatible with Airgo/Hugo side-folding rollators (except Sidekick)
- Compatible with Airgo eXcursion rollators of Revision 2.1 or higher (excluding 2.5)

[Shop Now](#)



## Replacement Sponge for Round Handle

**Item #**

730-340

- Comfortable replacement hypalon sponge cane handle

[Shop Now](#)



# Cane Accessories

## Apollo Rechargeable LED Safety Light

**Item #**

730-950

- Includes 3 super-bright LEDs
- Feature 3 usage modes: High, Low, and Strobe
- Rechargeable in 3.5 hours with included USB cable
- 20 lumen design for increased safety
- Compatible with all mobility aids

[Shop Now](#)



# Crutches

## Aluminum Crutches with Accessories

Item #	Description
10401-8	Youth
10400-8	Adult
10402-8	Tall Adult

- Lighter than wood, the standard aluminum crutch provides stability and durability
- Double extruded center tube provides additional strength to weight-bearing area
- Pushpin adjustment makes it easy-to-adjust crutch height in 1" increments
- Non-skid, jumbo-size vinyl contoured tips provide excellent traction
- Easy wing nut adjustments for hand grips
- Comfortable, durable accessories are latex free
- Weight Capacity 350 lb on aluminum crutch (excluding pediatric)

[Shop Now](#)



## Retail Packaged Aluminum Crutches with Accessories

Item #	Description
10401-1	Youth
RTL10400	Adult
RTL10402	Tall Adult

- Lighter than wood, the standard aluminum crutch provides stability and durability
- Double extruded center tube provides additional strength to weight-bearing area
- Pushpin adjustment makes it easy-to-adjust crutch height in 1" increments
- Non-skid, jumbo-size vinyl contoured tips provide excellent traction
- Easy wing nut adjustments for hand grips
- Comfortable, durable accessories are latex free
- Weight Capacity 350 lbs on aluminum crutch (excluding pediatric)

[Shop Now](#)





# Crutches

## EZ Adjust Aluminum Crutches with Euro-Style Clip and Accessories



Item #	Description
10431-8	Youth
10430-8	Adult
10432-8	Tall Adult

- Handgrip position adjusts easily and tool-free with Euro-style clip
- Height adjusts easily and tool-free with push pin
- Lighter than wood, the aluminum crutch provides stability and durability
- Non-skid, jumbo-size vinyl contoured tips provide excellent traction
- Hand grips are easily adjusted with Euro-style clips
- Comfortable, durable underarm pad and hand grip

[Shop Now](#)



## Universal Aluminum Crutch with Accessories



**Item #**  
RTL10433

- Universal adjustment for use by a child, adult or tall adult
- Euro-Style clips and pushpin adjustment make it easy-to-adjust height of hand grip
- Lighter than wood, the aluminum crutch provides stability and durability
- Crutch disassembles for convenient storage
- Non-skid, jumbo-size vinyl contoured tips provide excellent traction
- Comfortable, durable underarm pad and hand grip

[Shop Now](#)



# Crutches

## Bariatric Steel Crutches with Accessories

Item #	Description
10416-1	Youth
10406	Adult
10408	Tall Adult

- Durable, steel construction
- Push buttons easily adjust height in 1" increments
- Comfortable, contoured underarm design
- Easy wing nut adjustments for hand grips
- Large, non-skid jumbo-size vinyl contoured tips
- Weight Capacity 500 lbs

[Shop Now](#)



## Steel Forearm Crutches

Item #	Description
10403	Adult
10405	Tall Adult

- Vinyl-coated, ergonomically contoured arm cuffs molded for comfort and stability
- Leg and forearm sections adjust independently for optimal sizing
- Vinyl hand grips are comfortable and durable
- Extra-large tips for added stability
- Ortho K Handle is hard plastic

[Shop Now](#)



# Crutches & Accessories

## Bariatric Steel Forearm Crutch

**Item #**  
10403HD

**Description**  
Adult



- Vinyl-coated, ergonomically contoured arm cuffs molded for comfort and stability
- Leg and forearm sections adjust independently for optimal sizing
- Vinyl hand grips are comfortable and durable
- Extra-large tips for added stability

[Shop Now](#)



## Crutch Tips (7/8" crutch diameter)

Item #	Color	Description
RTL10439B	Grey	1 Pair
10439-8	Grey	8-Pairs
RTL10374BKB	Black	1 Pair

- For use with Drive crutch model #'s: 10430-8, 10431-8, 10432-8, 10400-8, 10401-8, 10402-8

[Shop Now](#)



# Crutches & Accessories

## Crutch Hand Grips (closed style)

**Item #**

RTL10453

- For use with Drive cane model #'s: 10400, 10401, 1402, 10416, 10430, 10431, 10432 series

[Shop Now](#)



## Crutch Underarm Cushions

**Item #**

RTL10446

- For use with: 10400, 10401, 10402, 10430, 10431, 10432 crutches

[Shop Now](#)



# Crutch Accessories

## Crutch Accessory Kit

### Item #

RTL10395

- Contains 1 pair Underarm Cushions (8 ½" x 7/8"), Hand Grips (4" x ¾") and Tips (7/8")

[Shop Now](#)

## Crutch Pillows

### Item #

RTL10355

- Adds increased comfort during long hours of crutch use
- Made from the finest materials that are resistant to body oils and perspiration
- Fast evaporation "moisture wicking" keeps your skin dry and comfortable
- Universal –fits all standard crutches
- Machine washable
- Latex free
- Crutch Pillow package includes: 2 underarm covers and 2 hand grip covers

[Shop Now](#)





CORPORATE OFFICE 🏠 99 Seaview Blvd., Port Washington, NY 11050 📞 877-374-8763 🌐 [www.drivemedical.com](http://www.drivemedical.com)

CANADA OFFICE 🏠 256 Aviva Park Drive, Vaughan, Ontario Canada L4L 9C7 📞 866-752-1363

©2024 Medical Depot, Inc. dba Drive DeVilbiss Healthcare, 99 Seaview Boulevard, Port Washington, NY 11050. All trademarks used in association with the sale of products of Drive DeVilbiss are trademarks owned by Medical Depot, Inc. All other trademarks, trade names, service marks and logos referenced herein belong to their respective companies. Refer to product instruction guide for warranty details.