



Shower Chairs | Good, Better, Best Options



★ GOOD ★									
Standard Shower Chairs									
Seat Type	Basic Round Shower Stool	Shower Stool	Diamond Spa Stool	Deluxe Aluminum Shower Bench without Back	Deluxe Aluminum Shower Bench without Back, Knock Down	Adjustable Bath Seat without Backrest, Not Assembled	Folding Shower Chair	Ergonomic Adjustable Bath Seat without Backrest	Round Shower Stool with Swivel Seat
Product Name	Basic Round Shower Stool	Shower Stool	Diamond Spa Stool	Deluxe Aluminum Shower Bench without Back	Deluxe Aluminum Shower Bench without Back, Knock Down	Adjustable Bath Seat without Backrest, Not Assembled	Folding Shower Chair	Ergonomic Adjustable Bath Seat without Backrest	Round Shower Stool with Swivel Seat
Item #	RTL12004KD	770-514	RTL12300	12203KD-4	RTL12203KDR	770-500	12486	770-535	RTL12061M
Product Image									
Sales Channel									
Key Differences	Standard shower stool with small footprint, ideal for smaller tubs/showers.	Attractive white stool with aqua rubber accents.	Standard shower stool with small footprint, ideal for smaller tubs/showers. Clean and modern design, featuring a caddy to store items.	Partial knock-down, only assembly required is the legs.	Full knock-down style to allow for small footprint retail package.	Attractive white bench with aqua rubber accents.	Folding shower bench allowing for easy storage.	Contoured white seat with aqua rubber accents. Seat includes cutouts on each side to hold shower head.	Rotating stool allowing for ease of movement getting in an out of tub/shower.
Assembly Required	✓	✓	✓	✓	✓	✓	✓	✓	✓
With Backrest	--	--	--	--	--	--	--	--	--
With Armrests	--	--	--	--	--	--	--	--	--
Contoured Design	--	--	--	--	--	--	--	✓	--
Frame Material	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum	Plastic
Seat Material	Plastic	Plastic	Plastic	Plastic	Plastic	Plastic	Plastic	Plastic	Plastic
Color	White	White	White	Grey	Grey	White	White	White	White
Weight Capacity	300 lbs	250 lbs	300 lbs	300 lbs	300 lbs	300 lbs	300 lbs	300 lbs	450 lbs
Handles	--	--	--	✓	✓	✓	✓	✓	--
Water Drainage Holes	--	--	✓	✓	✓	✓	✓	--	--
Non-Slip Rubber Tips	✓	✓	✓	✓	✓	✓	✓	✓	✓
Shower Spray Cutout	--	--	--	--	--	--	--	✓	--
Outside Leg Dimensions	12.5" x 12.5"	18" x 18"	18" x 18"	19.5" x 19.5"	19.5" x 19.5"	19.5" x 19.5"	20.5" x 19.25"	20" x 17"	15" x 15"
Seat Dimensions	13" x 13"	12.8" x 12.8"	12" x 12"	19.25" x 11.5"	19.25" x 11.5"	20" x 12"	19.75" x 12"	20" x 12"	14" x 14"
Width Between Arms	--	--	--	--	--	--	--	--	--
Adjustable Height	13.5" - 21"	15" - 20"	--	15.5" - 19.5"	15.5" - 19.5"	14" - 21"	17.5" - 18.5"	15" - 22"	16" - 22"
Product Weight	4.48 lbs	3.58 lbs	2.73 lbs	5.4 lbs	6.07 lbs	4.44 lbs	6.7 lbs	5.5 lbs	6.6 lbs
Shipping Weight	4.48 lbs	3.58 lbs, 21.48 lbs	15.43 lbs	21.6 lbs	6.07 lbs	4.44 lbs, 8.88 lbs	6.7 lbs	5.5 lbs, 11 lbs	6.6 lbs, 26.4 lbs
Case Pack	1/EA	1/EA, 6/CV	1/EA, 4/CV	4/CV	1/EA	1/EA, 2/CV	1/EA	1/EA, 2/CV	1/EA, 4/CV
Retail Packaged	✓	✓	✓	--	✓	✓	✓	✓	✓
Retail Box Dimensions (LxWxH)	15.75" x 3.25" x 14.25"	17" x 3.25" x 14.5"	13.19" x 3.35" x 18.9"	--	22" x 5.5" x 14.5"	13" x 4.75" x 22.5"	20" x 4" x 15.75"	13" x 4.5" x 22.25"	15" x 4.75" x 16.25"
Warranty	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime

LEGEND



ECOMMERCE
Available to order online



RETAIL
Available in large retail chain stores



PROVIDER
Available in independent DME provider locations

GOOD ★

Includes your most basic models and features that meet the customer minimum requirements and provide a good value.

BETTER ★★

Includes similar design and features from our Good models, our Better models offer supplementary features that give the consumer a better user experience.

BEST ★★★

Includes similar design and features from our Good and Better models, our Best models offers unique high-end features providing consumers with a top level experience.



Shower Chairs | Good, Better, Best Options



	★★ BETTER ★★					★★★★ BEST ★★★★★			★ GOOD ★	★★ BETTER ★★	★★★★ BEST ★★★★★
Seat Type	Standard Shower Chairs								Bariatric Shower Chairs		
Product Name	Deluxe Aluminum Bath Chair, Knock Down	Deluxe Aluminum Bath Chair, Knock Down	Adjustable Bath Seat with Backrest, Not Assembled	Ergonomic Adjustable Bath Seat with Backrest	Ergonomic Folding Bath Seat	Shower Chair with Back and Removable Padded Arms	Premium Series Shower Chair with Back and Arms	PreserveTech™ 360° Swivel Bath Chair	Deluxe Bariatric Shower Chair with Cross-Frame Brace, without Back, Knock Down	Deluxe Bariatric Shower Chair with Cross-Frame Brace Knock Down	Deluxe Bariatric Shower Chair with Cross-Frame Brace
Item #	12202KD-4	RTL12202KDR	770-510	770-537	770-525	12445KD-1	RTL12505	RTL12A001-GR	12022KD-1	12021KD-1	12029-2
Product Image											
Sales Channel											
Key Differences	Partial knock-down, only assembly required is the legs.	Full knock-down style to allow for small footprint retail package.	Attractive white stool with aqua rubber accents and large backrest.	Contoured white seat with aqua rubber accents. Seat includes cutouts on each side to hold shower head.	Contoured white seat with aqua rubber accents. Seat includes cutouts on each side to hold shower head. Folds for easy storage.	Standard bath seat with removable padded armrests.	Lightweight standard shower chair with removable legs and arms.	Compact bath chair designed to aid in preventing skin irritation and friction wounds featuring 360° swivel action in either direction and automatically locks at 90° intervals.	Featuring a cross-brace frame for added security.	Full knock-down style with back and cross-brace frame for added security.	Partial knock-down style with back and cross-brace frame for added security. Also featuring a 600lbs weight capacity.
Assembly Required	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓ (back only)
With Backrest	--	--	✓	✓	✓	✓	✓	✓	--	✓	✓
With Armrests	--	--	--	--	--	✓	✓	✓	--	--	--
Contoured Design	--	--	--	✓	✓	--	--	--	--	--	--
Frame Material	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum	Plastic	Aluminum	Aluminum	Aluminum	Aluminum
Seat Material	Plastic	Plastic	Plastic	Plastic	Plastic	Plastic	Plastic	Plastic, EVA	Plastic	Plastic	Plastic
Color	Grey	Grey	White	White	White	Grey	White	Grey	White	White	White
Weight Capacity	300 lbs	300 lbs	300 lbs	300 lbs	250 lbs	300 lbs	350 lbs	350 lbs	500 lbs	500 lbs	600 lbs
Handles	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Water Drainage Holes	✓	✓	✓	--	--	✓	✓	✓	✓	✓	✓
Non-Slip Rubber Tips	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Shower Spray Cutout	--	--	--	✓	✓	--	--	✓	--	--	--
Outside Leg Dimensions	19.5" x 19.5"	19.5" x 19.5"	--	20" x 17"	15.53" x 9.25"	17.75" x 19.25"	17" x 15.5"	17"	17.5" x 17.25"	17.5" x 17.25"	21.25" x 18.5"
Seat Dimensions	19.25" x 11.5"	19.25" x 11.5"	20" x 12"	20" x 12"	19.76" x 12.2"	16" x 14.5"	21.75" x 16.5"	16" x 18"	20" x 12"	20" x 12"	16.5" x 16.5"
Width Between Arms	--	--	--	--	--	18.5"	20.38"	18"	--	--	--
Adjustable Height	15.5" - 19.5"	15.5" - 19.5"	14" - 21"	15" - 22"	15" - 21"	16" - 21"	16" - 20.5"	15" - 22"	14" - 19"	14" - 19"	16.5" - 21.5"
Product Weight	7.45 lbs	9.07 lbs	7.04 lbs	7.04 lbs	8.86 lbs	7.72 lbs	11.5 lbs	8.8 lbs	6.6 lbs	9.1 lbs	11.8 lbs
Shipping Weight	29.8 lbs	9.07 lbs	7.04 lbs, 14.08 lbs	7.04 lbs, 14.08 lbs	8.86 lbs, 17.73 lbs	7.72 lbs	11.5 lbs, 46 lbs	10.3 lbs	6.6 lbs	9.1 lbs	23.6 lbs
Case Pack	4/CV	1/EA	1/EA, 2/CV	1/EA, 2/CV	1/EA, 2/CV	1/EA	1/EA, 4/CV	1/EA	1/EA	1/EA	2/CV
Retail Packaged	--	✓	✓	✓	✓	--	✓	✓	--	--	--
Retail Box Dimensions (LxWxH)	--	22" x 6" x 13"	13" x 5.5" x 22.75"	14" x 6" x 22"	14" x 20.75" x 6"	--	22.25" x 7" x 18"	20.8" x 7.09" x 16.15"	--	--	--
Warranty	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime

LEGEND



GOOD ★
Includes your most basic models and features that meet the customer minimum requirements and provide a good value.

BETTER ★★
Includes similar design and features from our Good models, our Better models offer supplementary features that give the consumer a better user experience.

BEST ★★★
Includes similar design and features from our Good and Better models, our Best models offers unique high-end features providing consumers with a top level experience.