



Raised Toilet Seats & Frames | Good, Better, Best Options



	★ GOOD ★							★★ BETTER ★★	
	Basic Raised Toilet Seat							Hinged Toilet Seat Risers	
Product Name	Economy Raised Toilet Seat	Locking Raised Toilet Seat	Raised Toilet Seat, 4"	Raised Toilet Seat w/ Lid, 4"	Raised Toilet Seat, 6"	Raised Toilet Seat w/ Lid, 2.5"	Raised Toilet Seat w/ Lid, 4"	Hinged Toilet Seat Riser, 3.5"	Hinged Toilet Seat Riser, 3.5'
Item #	770-610	12014	RTL12064	12065	12066	770-625	770-626, 770-629	RTL12607, RTL12608	770-630, 770-631
Product Image									
Sales Channel									
Key Differences	Lightweight, portable and universally sized to fit most toilet seats.	Wide opening in front and back for personal hygiene.	Lightweight and portable raised toilet seat. -Canada Only-	Lightweight and portable raised toilet seat with lid.	Lightweight and portable raised toilet seat.	Easy to clean, contoured raised toilet seat with lid providing 2.5" of additional height.	Easy to clean, contoured raised toilet seat with lid providing 4" of additional height.	Hinge feature allows for a more secure fit while facilitating cleaning and hygiene for all users. When placed in the raised position, provides an unobstructed area, for men.	Hinge feature allows for a more secure fit while facilitating cleaning and hygiene for all users. When placed in the raised position, provides an unobstructed area, for men. Featuring a 300lb weight capacity.
Toilet Seat Fit	Fits most toilets	Elongated	Standard	Standard	Standard	Standard	Standard, Elongated	Standard, Elongated	Standard, Elongated
Assembly Required	✓	✓	✓	✓	✓	✓	✓	Hardware Included	Hardware Included
Color	White	White	White	White	White	White	White	White	White
Material	Plastic	Plastic	Plastic	Plastic	Plastic	Plastic	Plastic	Plastic	Plastic
Weight Capacity	400 lbs.	300 lbs.	400 lbs.	400 lbs.	400 lbs.	300 lbs.	300 lbs.	250 lbs.	300 lbs.
Height	4"	4.5"	4"	4"	6"	2.5"	4"	3.5"	3.5"
Lid	--	--	--	✓	--	✓	✓	--	--
Removable Arms	--	--	--	--	--	--	--	--	--
Width Between Arms	--	--	--	--	--	--	--	--	--
Locking Mechanism	--	Knob (Worm Screw)	Tightening Clamps	Tightening Clamps	Tightening Clamps	Tightening Clamps	Tightening Clamps	Bolted Into Toilet Frame	Bolted Into Toilet Frame
Seat Dimensions (LxWxH)	15.75" x 15.5" x 4"	17" x 16.5" x 4.5"	16" x 14" x 4"	16" x 14" x 4"	16" x 14" x 6"	15.25" x 14" x 2.5"	15.25" x 14" x 4" (Standard) 17.5" x 14" x 4" (Elongated)	17.5" x 13.5" x 3.5" (Standard) 19.25" x 14" x 3.5" (Elongated)	17.25" x 13.75" x 3.5" (Standard) 19.25" x 13.75" x 3.5" (Elongated)
Product Weight	3.96 lbs	4 lbs.	3.5 lbs	3.5 lbs	3 lbs, 5 lbs	3.82 lbs	4.41 lbs, 5.07 lbs	5.1 lbs	3.3 lbs, 3.74 lbs
Shipping Weight	3.96 lbs	6 lbs	3.5 lbs	3.5 lbs	3 lbs, 5 lbs	3.82 lbs, 22.89 lbs	4.41 lbs, 17.64 lbs 5.07 lbs, 20.28lbs	5.1 lbs	3.3 lbs, 9 lbs 3.74 lbs, 11.22 lbs
Case Pack	1/EA	1/EA	1/EA	1/EA	1/EA	1/EA, 6/CV	1/EA, 4/CV	1/EA	1/EA, 3/CV
Retail Packaged	✓	✓	✓	✓	✓	✓	✓	✓	✓
Retail Package Dimensions	15.25" x 6.25" x 16"	18.25" x 7.5" x 17.75"	16.5" x 3.9" x 15"	16" x 5.75" x 14.5"	4" - 15.25" x 7" x 16.5"	8.25" x 18.5" x 16.25"	2" - 16.25" x 6" x 15.25" 4" - 16.5" x 7" x 15"	18.75" x 4.5" x 14.6"	14" x 4" x 14.5"
Warranty	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime

LEGEND



ECOMMERCE
Available to order online



RETAIL
Available in large retail chain stores



PROVIDER
Available in independent DME provider locations

GOOD ★

Includes your most basic models and features that meet the customer minimum requirements and provide a good value.

BETTER ★★

Includes similar design and features from our Good models, our Better models offer supplementary features that give the consumer a better user experience.

BEST ★★★

Includes similar design and features from our Good and Better models, our Best models offers unique high-end features providing consumers with a top level experience.



Raised Toilet Seats & Frames | Good, Better, Best Options



★★★★ BEST ★★★★★							
Raised Toilet Seat with Armrests							
Product Name	Premium Raised Toilet Seat, Standard	Premium Raised Toilet Seat, Elongated	2-in-1 Raised Toilet Seat	3-in-1 Raised Toilet Seat, 2"	3-in-1 Raised Toilet Seat, 4"	PreserveTech™ Secure Lock Raised Toilet Seat with Armrests	PreserveTech™ Universal Raised Toilet Seat
Item #	12402, 12403	12013	RTL12027RA	770-617	770-618	RTL12C003-WH	RTL12C002-WH
Product Image							
Sales Channel							
Key Differences	Featuring removable arms and allows the use of existing toilet seat and lid.	Raised toilet seat featuring extra-wide space between arms.	Featuring a body-contoured seat and removable arms.	Featuring a body-contoured seat and removable arms; use with no arms, one arm or two. Attractive white color with aqua accents.	Featuring a body-contoured seat and removable arms; use with no arms, one arm or two. Attractive white color with aqua accents.	Three locking points of security; located on the front, left and right side providing the highest level of stability and security.	Accommodates three different seat heights (2", 4" and 6") by maneuvering just two components.
Toilet Seat Fit	Standard, Elongated	Elongated	Fits most toilets	Fits most toilets	Fits most toilets	Standard, Elongated	Standard, Elongated
Assembly Required	✓	✓	✓	Hardware Included	Hardware Included	✓	✓
Color	White/Grey	White/Grey	White/Grey	White	White	White	White
Material	Plastic	Plastic	Plastic	Plastic	Plastic	Plastic	Plastic
Weight Capacity	300 lbs.	300 lbs.	300 lbs.	300 lbs.	300 lbs.	350 lbs	350 lbs
Height	3.5"	5"	5"	2"	4"	5"	2", 4", 6"
Lid	✓	--	--	--	--	--	✓
Removable Arms	✓	✓	✓	✓	✓	✓	✓
Width Between Arms (in)	19.75"	21"	18"	18.11"	18.19"	17.5"	18.5"
Locking Mechanism	Bolted Into Toilet Frame	Knob (Worm Screw)	Knob (Worm Screw)	Knob (Worm Screw)	Knob (Worm Screw)	Knob (Worm Screw)	Bolted Into Toilet Frame
Seat Dimensions (LxWxH)	17" x 13.75" x 3.5" (Standard) 19.25" x 13.75" x 3.5" (Elongated)	17" x 15.5" x 5"	17" x 16.5" x 5"	17.72" x 20.87"	17.72" x 20.87"	20.3" x 17.5" x 5"	16" x 19" x 6.5"
Product Weight	3.5 lbs.	4.5 lbs.	5 lbs	7.61 lbs	7.39 lbs	5.29 lbs	9.3 lbs
Shipping Weight	4.5 lbs	6.5 lbs	8.5 lbs	7.61 lbs, 15.21 lbs	7.39 lbs, 14.78lbs	6 lbs	11 lbs
Case Pack	1/EA	1/EA	1/EA	1/EA, 2/CV	1/EA, 2/CV	1/EA	1/EA
Retail Packaged	✓	✓	✓	✓	✓	✓	✓
Retail Package Dimensions	14.75" x 8" x 18.25"	18.5" x 8.5" x 16.25"	17.25" x 7.5" x 20.25"	21" x 7" x 18.25"	8.25" x 21" x 18.75"	20.9" x 7.5" x 18.1"	8.6" x 20.25" x 18.5"
Warranty	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime

LEGEND



ECOMMERCE
Available to order online



RETAIL
Available in large retail chain stores



PROVIDER
Available in independent DME provider locations

GOOD ★

Includes your most basic models and features that meet the customer minimum requirements and provide a good value.

BETTER ★★

Includes similar design and features from our Good models, our Better models offer supplementary features that give the consumer a better user experience.

BEST ★★★

Includes similar design and features from our Good and Better models, our Best models offers unique high-end features providing consumers with a top level experience.



	★ GOOD ★			★★ BETTER ★★		★★★ BEST ★★★
	Toilet Frames					
Product Name	Toilet Safety Frame	Adjustable Toilet Safety Frame	Toilet Safety Frame	Toilet Safety Frame	Toilet Safety Frame	Free Standing Toilet Safety Frame
Item #	12001-4	770-665	RTL12000	RTL12087	770-660	RTL12079
Product Image						
Sales Channel						
Key Differences	Adjustable toilet safety frame with plastic armrests.	Adjustable toilet safety frame with full length plastic armrests. Attractive white color with aqua accents.	Standard toilet safety frame with padded handles.	Adjustable toilet safety frame with padded armrests.	Standard toilet safety frame with full length padded armrests. Attractive white color with aqua accents.	Free standing toilet safety frame with padded armrests and magazine rack.
Assembly Required	✓	✓	✓	✓	✓	✓
Color	Silver/Grey	White/Blue	White/Grey	White/Grey	White/Blue	White/Grey
Padded Armrests	--	--	✓	✓	✓	✓
Material	Anodized aluminum	Steel	Anodized aluminum	Anodized aluminum	Steel	Anodized aluminum
Weight Capacity	300 lbs.	250 lbs.	300 lbs.	300 lbs.	300 lbs	300 lbs.
Width Between Arms	18.5" - 22.75"	18" - 21"	19" - 23.5"	16.5" - 19.5"	15.5" - 18"	21"
Adjustable Arm Height	22.5" - 29.5"	28" - 32"	28" - 32"	--	--	--
Product Weight	4 lbs	3.52 lbs	4.8 lbs	5.4 lbs	5.73 lbs	12.7 lbs
Shipping Weight	16 lbs	3.52 lbs	4.8 lbs, 20.28 lbs	5.4 lbs	5.73 lbs	12.7 lbs
Case Pack	4/CV	1/EA	1/EA, 4/CV	1/EA	1/EA	1/EA
Retail Packaged	--	✓	✓	✓	✓	✓
Retail Box Dimensions (LxWxH)	--	18" x 3" x 13"	19" x 3.25" x 16.75"	7.75" x 3.25" x 17"	17" x 3.25" x 7.75"	26" x 5.12" x 19.69"
Warranty	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime

LEGEND



ECOMMERCE
Available to order online



RETAIL
Available in large retail chain stores



PROVIDER
Available in independent DME provider locations

GOOD ★
Includes your most basic models and features that meet the customer minimum requirements and provide a good value.

BETTER ★★
Includes similar design and features from our Good models, our Better models offer supplementary features that give the consumer a better user experience.

BEST ★★★
Includes similar design and features from our Good and Better models, our Best models offers unique high-end features providing consumers with a top level experience.