



Transfer Benches | Good, Better, Best Options



	★ GOOD ★			★★ BETTER ★★		★★★ BEST ★★★			
Bench Type	Standard Transfer Bench						Specialty Transfer Bench		Bariatric Transfer Bench
Product Name	Transfer Tub Bench	Three Piece Transfer Bench	Transfer Bench	Padded Transfer Bench	PreserveTech™ Transfer Bench	Splash Defense™ Transfer Bench with Curtain Guard Protection	Combination Padded Transfer Bench/Commode	Folding Universal Sliding Transfer Bench	Bariatric Transfer Bench
Item #	12011KD-1	RTL12031KDR	770-406	12005KD-1	12011KD-1A	RTL12032KDR	12005KDC-1	RTL12075	12025KD-1
Product Image									
Sales Channel									
Key Differences	Standard transfer bench with reversible capabilities to accommodate any bathroom.	Standard transfer bench with reversible capabilities to accommodate any bathroom.	Attractive white transfer bench with aqua rubber accents and reversible capabilities to accommodate any bathroom.	Standard transfer bench featuring padded seat and back for added comfort.	Drive PreserveTech™ products are manufactured with an antimicrobial* agent that helps to minimize the growth of stain and odor-causing bacteria and mold on treated surfaces. *Antimicrobial properties are built in to protect the products. These products do not protect users or others against bacteria, viruses, germs or other disease organisms	Standard transfer bench with patented u-shape design enables the shower curtain to stay in place preventing water from spilling onto the bathroom floor.	Combines a transfer bench and commode into one, featuring a padded seat and back for added comfort.	Sliding Transfer Bench moves from left to right to provide easy entry into and out of the tub. Added features include a soap dish and shower spray holder.	Bariatric transfer bench with reversible capabilities to accommodate any bathroom.
Assembly Required	✓	✓	✓	✓	✓	✓	✓	✓	✓
Frame Material	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum
Seat Material	Plastic	Plastic	Plastic	PVC	Plastic	Plastic	PVC	Plastic	Plastic
Color	Grey	Grey	White	Grey	Grey	White	Grey	White	Grey
Weight Capacity	400 lbs	400 lbs	250 lbs	400 lbs	400 lbs	400 lbs	400 lbs	300 lbs	500 lbs
Outside Leg Dimensions	32" x 23"	32" x 23"	32" x 23"	29" x 19.75"	32" x 23"	33.5" x 23"	17.25" x 28.5"	18.3" x 42.4"	32.5" x 24"
Seat Dimensions	26" x 18.5" - 19.5"	26" x 18.5" - 19.5"	26.4" x 18.5"	24" x 16"	26" x 18.5" - 19.5"	28" x 16"	24" x 16"	22.5" x 18.5"	26.5" x 16.5"
Adjustable Height	1/2" increments	1/2" increments	1/2" increments	1/2" increments	1/2" increments	1/2" increments	1/2" increments	1/2" increments	1/2" increments
Adjustable Height Range	17.5" - 22.5"	17.5" - 22.5"	18" - 22"	17.75" - 21.75"	17.5" - 22.5"	17.5" - 22.5"	18.5" - 22.5"	21" - 25"	18.5" - 23.5"
Product Weight	13.6 lbs	13.96 lbs	11.68 lbs	13.5 lbs	10.8 lbs	14.1 lbs	15.4 lbs	14 lbs	16.5 lbs
Shipping Weight	13.6 lbs	13.96 lbs	11.68 lbs, 23.37 lbs	13.5 lbs	12.5 lbs	14.1 lbs	15.4 lbs	14 lbs, 28 lbs	16.5 lbs
Case Pack	1/EA	1/EA	1/EA, 2/CV	1/EA	1/CV	1/EA	1/EA	1/EA, 2/CV	1/EA
Retail Packaged	--	✓	✓	--	--	✓	--	✓	--
Retail Box Dimensions (LxWxH)	--	18.75" x 10.25" x 16.5"	17" x 10" x 19"	--	--	20" x 8.75" x 17"	--	30.25" x 9.5" x 19"	--
Warranty	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime

LEGEND



ECOMMERCE
Available to order online



RETAIL
Available in large retail chain stores



PROVIDER
Available in independent DME provider locations

GOOD ★

Includes your most basic models and features that meet the customer minimum requirements and provide a good value.

BETTER ★★

Includes similar design and features from our Good models, our Better models offer supplementary features that give the consumer a better user experience.

BEST ★★★

Includes similar design and features from our Good and Better models, our Best models offers unique high-end features providing consumers with a top level experience.