



	★ GOOD ★	★★ BETTER ★★			★★★ BEST ★★★
	Tub Rails				
Product Name	Plastic Tub Rail	Clamp-On Tub Rail	Adjustable Clamp-On Tub Rail	Adjustable Clamp-On Tub Rail	Multi-Adjust Bath Safety Rail
Item #	12044	12036	785-350	RTL12036-ADJ	785-340
Product Image					
Sales Channel					
Key Differences	<i>Plastic parallel tub rail featuring a white finish and rubber pads to prevent scratches on the tub and provide a secure fit.</i>	<i>Durable steel parallel tub rail featuring a white finish and rubber pads to prevent scratches on the tub and provide a secure fit.</i>	<i>Durable steel parallel tub rail featuring a white finish and rubber pads to prevent scratches on the tub and provide a secure fit. Attractive white color with aqua accents.</i>	<i>Durable steel parallel tub rail featuring a white finish and rubber pads to prevent scratches on the tub and provide a secure fit. Added rubber grip to provide extra security.</i>	<i>Features 4 different height settings and allows for both parallel and perpendicular orientation to the tub wall.</i>
Installation	Tool free	Tool free	Tool free	Tool free	Tool-free
Color	White	White	White/Blue	White/Black	White/Blue
Material	Plastic	Steel	Steel	Steel	Steel
Weight Capacity	300 lbs	300 lbs	300 lbs	300 lbs	300 lbs
Clamp Adjustments	3" - 6.25"	3" - 7"	3" - 7"	3" - 7"	3" - 7"
Product Height	16.5"	14.5"	15"	14.5" - 17"	13.5" - 16.5"
Product Weight	6.1 lbs	7.7 lbs	6.6 lbs	7.94 lbs	6.68 lbs
Shipping Weight	6.1 lbs, 36.6 lbs	7.7 lbs, 46.2 lbs	6.6 lbs, 39.6 lbs	7.94 lbs	6.68 lbs, 40.09 lbs
Case Pack	1/EA, 6/CV	1/EA, 6/CV	1/EA, 6/CV	1/EA	1/EA, 6/CV
Retail Packaged	✓	✓	✓	✓	✓
Retail Box Dimensions (LxWxH)	8.25" x 6.75" x 21.75"	8.25" x 6" x 19.5"	6.75" x 6" x 19.5"	8" x 6" x 20.75"	7.25" x 6.75" x 12.75"
Warranty	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	5 Year Limited

LEGEND



ECOMMERCE
Available to order online



RETAIL
Available in large retail chain stores



PROVIDER
Available in independent DME provider locations

GOOD ★

Includes your most basic models and features that meet the customer minimum requirements and provide a good value.

BETTER ★★

Includes similar design and features from our Good models, our Better models offer supplementary features that give the consumer a better user experience.

BEST ★★★

Includes similar design and features from our Good and Better models, our Best models offers unique high-end features providing consumers with a top level experience.