

NEW

Lightweight Bariatric Beds



Upgrade your bariatric patient homecare bed offering with Drive DeVilbiss' new lightweight Bariatric Homecare Beds. These beds offer an easier to transport lightweight bed frame while maintaining the same durability and safe working load for patients. **With frames up to 100 lbs.* lighter than current Bariatric Beds**, these beds make deliveries less cumbersome and maneuverability and assembly more manageable.



Lighter frame provides equal support



Split, slat deck frame offers ease of cleaning, additional frame strength, and "no sag" mattress support



Large glow in the dark hand pendant has easy-to-use controls



Assembly and disassembly is simple, tool free, and easy to complete by even 1 individual



Durable rich brown wood grain panel bed ends are attractive and scratch resistant

ITEM #	PRODUCT DESCRIPTION	PRODUCT WEIGHT*	UPC	UOM
15300LW	Lightweight, Full Electric 42" Bariatric Bed	~70 lbs lighter	822383008578	2 Boxes
15302LW	Lightweight, Full Electric 48" Bariatric Bed	~80 lbs lighter	822383007717	3 Boxes
15303LW	Lightweight, Full Electric 54" Bariatric Bed	~100 lbs lighter	822383007724	3 Boxes

PRODUCT SPECS	15300LW	15302LW	15303LW
Minimum/Maximum Height	15.75" / 24"	15.75" / 24"	15.75" / 24"
Product Dimensions	42" (W) x 88" (L)	48" (W) x 88" (L)	54" (W) x 88" (L)
Total Bed Weight	220 LBS	265 LBS	281 LBS
Weight Capacity	600 lbs	750 lbs	1000 lbs

Warranty	Limited Lifetime Welds & Frame 1 Year Parts & Components 5 Year Motors & Hand Pendants
----------	--

ACCESSORIES	
15302LWB	Backup Battery Pack
15021BV-LW	Half Rails
1500SC-42-FB	42" Therapeutic Foam Mattress
1500SC-48-FB	48" Therapeutic Foam Mattress
15301	42" Basic Foam Mattress
15310	48" Basic Foam Mattress
15312	54" Basic Foam Mattress

*Product weights in comparison to current Drive bariatric beds