


Delta[®] Pro Bed System




The innovative Delta[®] Pro Bed System is easy to assemble and maximizes inventory efficiency while maintaining competitive pricing. The first-of-its-kind bed system features a universal bed frame that is compatible with color-coded, interchangeable components that include bed ends and motor options to create six unique bed configurations. The components easily assemble and disassemble as needed to create a manual, semi-electric, full-electric or low bed depending on user needs.

See just how simple it is to create your ideal design!



BED FRAME



F

Bed Frame
Item # 15545



BED ENDS



E

Standard Bed End
Item # 15546



E

Low Bed End
Item # 15547



MOTOR OPTIONS



M

Dual Motor
Item # 15550MOHC



M

Dual Motor w/ Hi-Lo
Item # 15551



M M

Linear Motors
Item # 15553, 15552



M M M

Linear Motors w/ Hi Lo
Item # 15553, 15552, 15554



C

Manual Crank
Item # 15555

ITEM #	PRODUCT DESCRIPTION	COMPONENTS	UPC
15030N	Standard, Semi-Electric Dual Motor	F E M	822383006741
15033NC	Standard, Full-Electric Dual Motor + Hi/Lo	F E M	822383006758
15235N	Low, Full Electric Dual Motor + Hi/Lo	F E M	822383006765
15003PN	Standard, Manual	F E C	822383006772
15004PN	Standard, Semi-Electric Linear Motors	F E M M	822383006789
15005LPN	Low, Full Electric Linear Motors + Hi/Lo	F E M M M	822383006796

ITEM #	PRODUCT DESCRIPTION	BED FRAME	BED ENDS	MAIN MOTOR	FRAME ADJ.	HEIGHT ADJ.
15030N	Standard, Semi-Electric Dual Motor	15545	15546	15550MOHC		Manual
15033NC	Standard, Full-Electric Dual Motor + Hi/Lo	15545	15546		15551	
15235N	Low, Full Electric Dual Motor + Hi/Lo	15545	15547		15551	
15003PN	Standard, Manual	15545	15546	No Motor	15555	Manual
15004PN	Standard, Semi-Electric Linear Motors	15545	15546	15553	15552	Manual
15005LPN	Low, Full Electric Linear Motors + Hi/Lo	15545	15547	15553	15552	15554



F E M

Standard, Semi-Electric Dual Motor

Item # 15030N

Configuration components:
Universal Bed Frame, Standard Bed Ends,
Dual Motor



F E M

Standard, Full Electric Dual Motor with Hi/Lo

Item # 15033NC

Configuration components:
Universal Bed Frame, Standard Bed Ends,
Dual Motor with Hi/Lo



F E M

Low, Full Electric Dual Motor with Hi/Lo

Item # 15235N

Configuration components:
Universal Bed Frame, Low Bed Ends,
Dual Motor with Hi/Lo



F E C

Standard Manual

Item # 15003PN

Configuration components:
Universal Bed Frame, Standard Bed Ends,
Manual Crank



F E M M

Standard, Semi-Electric Linear Motors

Item # 15004PN

Configuration components:
Universal Bed Frame, Standard Bed Ends,
Linear Motors



F E M M M

Low, Full Electric Linear Motors with Hi/Lo

Item # 15005LPN

Configuration components:
Universal Bed Frame, Low Bed Ends,
Linear Motors plus Hi/Lo

PRODUCT SPECS

	STANDARD	LOW
Minimum/Maximum Height*	13.5" - 23"	10" - 23"
Product Dimensions	36" (W) x 86" (L)	36" (W) x 86" (L)
Weight Capacity	450 lb	450 lb
Warranty	Limited Lifetime Welds & Frame 1 Year Parts & Components 5 Year Motors	

ACCESSORIES

15561HRAIL	Half Rails
15561FRAIL	Telescoping Full Rails
15006	Innerspring Mattress
15670**	Gravity™ 6 Therapeutic Foam Mattress

*With Casters

**Visit drivemedical.com for additional therapeutic support mattress options