

PRIMECARE

Model P503

LONG TERM
Care



COMFORT. SAFETY. FLEXIBILITY.

The P503 Low Bed from Drive DeVilbiss is an affordable option for facilities wanting to provide residents and staff with a safe and comfortable bed while addressing tight capital budgets. The contemporary, new design will accommodate your residents with unmatched features and quality not seen in this category of beds.

STANDARD FEATURES


- Best-in-class range of motion, 7.8" - 30"
- Integrated 76" - 80" length extension
- Handheld, 8-function hand pendant with auto-contour
- Eight, 360-degree swivel casters that roll at any height, four casters lock
- Sturdy, quick-release bed ends
- Sealed bearings provide smooth, reliable and True Vertical Lift operation

OPTIONS

- Assist rails and assist bars that also fit the P703 and P903
- Six bed end styles and 12 finishes from which to choose

SPECIFICATIONS

Bed Width	35"
Bed Height Range	7.8" – 30"
Bed Length Range	76" – 80"
Weight	210 lbs.
Weight Capacity	Accommodate up to 450 lbs.
Warranty	Welds: 15 years, Frame: 10 years, Electronics: 2 years
Item #	P503



450 lbs.
Weight Capacity

450 lbs.
Weight Capacity



Bed Height Range
7.8" – 30"



Hand Controller



Floor Lock and
Wall Bumper

To Learn More/Request a Demo: Call our LTC Customer Service Hotline at 1-877-374-8763, Contact your Sales Representative, or Visit DriveMedical.com EXPAND THE POSSIBILITIES

drive

DeVilbiss
HEALTHCARE

Available in 35" & 42"

Brownlee Woods



Mill Creek



Quick Ship (QS) Bed Ends available in 3 finishes and 2 styles only:



OSLO
Oak



OSLO
Cherry



OSLO
Figured Mahogany

OSLO

35BWBE-A	35" BRNWOOD BED END, OSLO CHERRY, EA
35BWBE-B	35" BRNWOOD BED END, OSLO OAK, EA
35BWBE-FM	35" BRNWOOD BED END, OSLO FIG MAHOG, EA
42BWBE-A	42" BRNWOOD BED END, OSLO CHERRY, EA
42BWBE-B	42" BRNWOOD BED END, OSLO OAK, EA
42BWBE-FM	42" BRNWOOD BED END, OSLO FIG MAHOG, EA

35MCBE-A	35" MILL CRK BED END, OSLO CHERRY, EA
35MCBE-B	35" MILL CRK BED END, OSLO OAK, EA
35MCBE-FM	35" MILL CRK BED END, OSLO FIG MAHOG, EA
42MCBE-A	42" MILL CRK BED END, OSLO CHERRY, EA
42MCBE-B	42" MILL CRK BED END, OSLO OAK, EA
42MCBE-FM	42" MILL CRK BED END, OSLO FIG MAHOG, EA

For Long Term Care PrimeCare Standard Beds Only

Rotating Assist Bar


Compatible with the P503,
P703 and P903 Bed



Model #PLTCAB
Model #PLTCABL, left
Model #PLTCABR, right

Rotating Assist Rail

Compatible with the P503,
P703, and P903 Bed



Model #PLTCAR
Model #PLTCARL, left
Model #PLTCARR, right