


The Solite Pro from Drive Medical is a fully electrical profiled homecare bed that offers smooth, quiet operation and versatile positioning for maximum patient comfort. An easy-to-use hand pendant control allows patient or caregiver to adjust head and foot section position as well as the overall height of the bed.

The Solite Pro includes an electric backrest, knee break function, and improved actuators in the variable height mechanism that add to the strength of the bed, all while retaining an attractive, domestic look and feel.



Built in transport stand for storage and transport



Tool-free built-in bed extension extends from 80" to 84"



Available with 3 rail options - full, standard and low



Can be used in the low and reverse trendelenburg positions

## FEATURES

- Redesigned platform for improved patient comfort when profiling
- Individual sections strengthened for extended performance and reliability
- Full-length, standard height and half height side rails available for safety and ease of use (item #: 15006LEXRAIL, 15006LPSTDRAIL, 15006LSHRTRAIL)
- Improved integral bed length extension with safety features to prevent entrapment
- All individual sections of the bed weigh less than 55 lbs. (25 kg) for easier assembly
- Includes transport stand
- Improved auto-regressing backrest prevents pinching in the sacral area
- Infinitely variable foot-down tilt.
- Four locking casters for maneuverability and safety
- Easy-use handset controls all functions and includes lockout function
- Maximum user weight: 396.8 (180 kg)
- Maximum safe working load: 498.2 (226 kg)
- Warranty: Frame: 10-year | Mechanical: 5-year | Electrical: 2-year

## SPECIFICATIONS

ITEM #	PRODUCT DESCRIPTION	PRODUCT SPECIFICATIONS	PACKAGING SPECIFICATIONS	UPC
15006LP	Solite Pro Bed Package including: 15006BEDENDS (2X) 15006HEADDECK 15006FOOTDECK	See individual components below	See individual components below	822383014098
15006BEDENDS	Solite Pro Bed End Assembly (packed 1/ea)	<b>Dimensions</b> in: 39.5" (L) x 7.3" (W) x 34" (H) cm: 100.3 (L) x 18.5 (W) x 86.3 (H) <b>Weight</b> lbs: 46.3 kg: 21.0	<b>Dimensions</b> in: 49.6" (L) x 7.5" (W) x 37.8" (H) cm: 125.9 (L) x 19.5 (W) x 96.0 (H) <b>Weight</b> lbs: 55 kg: 24.9	822383014104
15006HEADDECK	Solite Pro Head Deck Assembly	<b>Dimensions</b> in: 45" (L) x 35.7" (W) x 6.85" (H) cm: 114.3 (L) x 90.6 (W) x 17.3 (H) <b>Weight</b> lbs: 50.7 kg: 22.9	<b>Dimensions</b> in: 35.4" (L) x 40.6" (W) x 7" (H) cm: 114.3 (L) x 90.6 (W) x 17.3 (H) <b>Weight</b> lbs: 62 kg: 28.1	822383014111
15006FOOTDECK	Solite Pro Foot Deck Assembly	<b>Dimensions</b> in: 40.8" (L) x 35.7" (W) x 7.4" (H) cm: 103.6 (L) x 90.6 (W) x 18.7 (H) <b>Weight</b> lbs: 44.1 kg: 20.0	<b>Dimensions</b> in: 35.4" (L) x 40.6" (W) x 8" (H) cm: 89.9 (L) x 103.1 (W) x 20.3 (H) <b>Weight</b> lbs: 55 kg: 24.9	822383014128


## ACCESSORIES (SOLD SEPARATELY)

ITEM #	PRODUCT DESCRIPTION	PRODUCT SPECIFICATIONS	PACKAGING SPECIFICATIONS	UPC
15006LEXRAIL	Grange Telescopic Side Rails	<b>Dimensions</b> in: 44.3"-63.4" (L) x 3.3" (W) x 21.7" (H) cm: 112.5-161.0 (L) x 8.3 (W) x 55.1 (H) <b>Weight</b> lbs: 35.3 kg: 16.0	<b>Dimensions</b> in: 46.4" (L) x 21.1" (W) x 5.1" (H) cm: 117.8 (L) x 53.5 (W) x 12.9 (H) <b>Weight</b> lbs: 40 kg: 18.1	822383014135
15006LPSTDRAIL	Elland Standard Height Grab Handle	<b>Dimensions</b> in: 10.6" (L) x 1.6" (W) x 29.5" (H) cm: 26.9 (L) x 135.8 (W) x 32.7 (H) <b>Weight</b> lbs: 6.6 kg: 2.9	<b>Dimensions</b> in: 28.8" (L) x 11.9" (W) x 4.5" (H) cm: 73.1 (L) x 30.2 (W) x 11.43 (H) <b>Weight</b> lbs: 15.4 kg: 7.0	822383014142
15006LSHRTRAIL	Elland Short Height Grab Handle	<b>Dimensions</b> in: 10.6" (L) x 3.1" (W) x 23.5" (H) cm: 26.9 (L) x 7.8 (W) x 59.69 (H) <b>Weight</b> lbs: 4.5 kg: 2.0	<b>Dimensions</b> in: 28.8" (L) x 11.9" (W) x 4.5" (H) cm: 68.3 (L) x 30.2 (W) x 11.43 (H) <b>Weight</b> lbs: 11.0 kg: 5.0	822383014159

## FREQUENTLY PURCHASED WITH...



**Therapeutic 5 Zone Mattress Replacement**  
Item # 15019



**Gravity 7 Mattress Replacement**  
item # 15770

Verified to IEC 60601-1 General Medical Electrical Equipment standard, IEC 60601-2-52 Medical Bed standard and IEC 606601-1-11 Medical Electrical Equipment standard specific to homecare requirements.