



Commode Type	★ GOOD ★				★★ BETTER ★★			★★★ BEST ★★★			★ GOOD ★		★★ BETTER ★★	
	Standard Commodes				Drop Arm Commodes			Specialty Commode		Bariatric Commode				
Item #	11148CE-4	11148-1	11148N-4	RTL11158KDR	11148-1A	11148NW-1	11105N-4	11125KD-1	11101W-2	11125PSKD-1	11114KD-1	11120SV-1F	11117N-1	11135-1
Product Name	Competitive Edge	Folding Steel Commode	Folding Steel Commode	Retail Packed	PreserveTech™ Folding Commode	Folding Deep Seat Commode	Deeper Seat Depth, Assembled	Deluxe Steel Drop Arm Commode	Steel Drop-Arm, Wheeled Commode	Drop Arm Steel Commode	Rolling Shower Chair and Commode	Wheeled Commode w/ Footrests	Bariatric Folding Commode	Deluxe Bariatric Drop-Arm Commode
Sales Channel														
Key Differences	Basic brown boxed commode, meets minimum customer needs requirement for feature set and coding for billing.	Standard brown boxed commode, similar to CE commode but with upgraded features including arm rests and increased bucket size. Comes with piggy-back labels.	Deeper seat, folding commode that is sold in a case of 4. Comes with piggy-back labels.	Only retail boxed commode, key features include increased maximum height, deeper seat depth and slightly wider width between arms.	Drive PreserveTech™ products are manufactured with an antimicrobial* agent that helps to minimize the growth of stain and odor-causing bacteria and mold on treated surfaces. *Antimicrobial properties are built in to protect the products. These products do not protect users or others against bacteria, viruses, germs or other disease organisms	Deep seat, folding commode in a brown box featuring a fresh take on color, an added toilet paper holder and a commode liner sample pack.	Standard steel brown boxed commode, including a deep seat and fixed frame. - Canada only -	Drop arm for ease of transfers, light weight commode with increased width between arms.	Drop arm for ease of transfers, comes standard with wheels, padded armrests and increased height range.	Drop arm for ease of transfers, comes standard with padded seat and deepest seat depth.	Wheeled shower commode for transportation with ease, including padded seat and upholstery back for added comfort.	Wheeled commode for transportation with ease, including padded seat, upholstery back, locking casters and removable footrests.	Standard bariatric commode with steel foldable frame.	Drop arm bariatric commode for ease of transfers, fixed frame and increased weight capacity.
Assembly	No Tools Required	No Tools Required	No Tools Required	No Tools Required	No Tools Required	Tools Required	No Tools Required	No Tools Required	No Tools Required	No Tools Required	No Tools Required	No Tools Required	No Tools Required	No Tools Required
Foldable	✓	✓	✓	✓	Y	✓	--	--	--	--	--	--	✓	--
Material	Steel	Steel	Steel	Steel	Steel	Steel	Steel	Steel	Steel	Steel	Aluminum	Steel	Steel	Steel
Finish	Powder Coated Matte	Powder Coated Matte	Powder Coated Matte	Powder Coated Matte	Powder Coated Matte	Powder Coated Gloss	Powder Coated Matte	Powder Coated Matte	Powder Coated Gloss	Powder Coated Matte	Matte	Powder Coated Gloss	Powder Coated Gloss	Powder Coated Gloss
Color(s)	Gray	Gray	Gray	Gray	Gray	White/Gray	Gray	Gray	Gray	Gray	Gray	Gray	Gray	Gray
Wheel Size	--	--	--	--	--	--	--	--	3"	--	3"	5"	--	--
Padded Seat	--	--	--	--	--	--	--	--	--	✓	✓	✓	--	--
Toilet Paper Holder	--	--	--	--	--	✓	--	--	--	--	--	--	--	--
Snap-On Seat & Lid	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Splash Guard	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	--	--	✓	--
Bucket w/ Lid	✓	✓	✓	✓	✓	✓	✓	--	--	--	✓	✓	✓	--
Bucket Capacity (qt)	7.5	7.5	7.5	7.5	7.5	7.5	12	7.5	7.5	12	12	7	12	12
Weight Capacity	350 lbs	350 lbs	350 lbs	350 lbs	350 lbs	350 lbs	350 lbs	300 lbs	300 lbs	300 lbs	300 lbs	350 lbs	650 lbs	1000 lbs
Seat Dimensions (WxD)	13" x 14.5"	13" x 14.5"	14" x 16.5"	14" x 16.5"	13" x 14.5"	14" x 16.5"	14" x 16.5"	13.75" - 16.5"	14" x 16.5"	15.25" x 16.75"	16" x 16"	17.5" x 16.5" (Cushion: 16" x 16")	13.75" x 16.5"	23.25" x 18.5"
Width Between Arms	19"	18"	18"	18"	18"	18"	18"	19.75"	18"	19"	18"	18.5"	24"	25.75"
Width Outside Arms	21"	22"	22"	22"	22"	22"	22"	23.5"	24"	23.75"	21.5"	22.25"	31.75"	26.25"
Height Adjustments	16.6" - 22.5"	16.5" - 22.5"	15.5" - 21.75"	17" - 23"	15.5" - 22.5"	15.5" - 21.75"	16" - 21.75"	20" - 24"	18" - 24"	19" - 23"	21"	--	15.5" - 22"	23.25" - 18.5"
Product Weight	15 lbs	11.5 lbs	13 lbs	13 lbs	11.5 lbs	13 lbs	12 lbs	18 lbs	22.5 lbs	21 lbs	14 lbs	28 lbs	19 lbs	21 lbs
Shipping Weight	46 lbs	12.9 lbs	59 lbs (case of 4)	17.5 lbs	13 lbs	14.75 lbs	17.5 lbs	20 lbs	48 lbs (case of 2)	23 lbs	18 lbs	31 lbs	25 lbs	26 lbs
Case Pack	4/CS	1/CS	4/CS	1/BX	1/CS	1/CS	4/CS	1/CS	2/CS	1/CS	1/CS	1/CS	1/CS	1/CS
Retail Box	--	--	--	✓	--	--	--	--	--	--	--	--	--	--
Retail Box Dimensions (LxWxH)	--	--	--	14.5" x 8" x 22.5"	--	--	--	--	--	--	--	--	--	--
HPCP	0163	0163	0163	0163	E0163	0163	0163	0163	0163	0163	0163	0163	0163	0163
Warranty	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime

COMMUNE ACCESSORIES							
Type	Commode Bucket	Splash Guard	Commode Bucket	Oblong Oversized Toilet Seat	Oblong Oversized Toilet Seat	Round Toilet Seat with Lid	Commode Liners
Item #	11106	11107	11108	11160-1	11161N-1	11150-1	RTL12095
For use with Commode Model #'s	11148N, 11125KD, 11101W, NRS185006, NRS185007 and other leading manufacturers' commodes	11105N-4, 11112-2, 11117N, 11114KD, 11125KD, 11148N-4	11114KD, 11125PSKD, 11135, 11105N, 11117N	11125KD, 11148N-4, 11105N-4 (has replaced 11116-1)	11117N-1	11105N-4, 11148-1, 11148-4, RTL11148KDR, 11148CE-4	All (except Bariatric Units) - features powder absorbent, tie handles, can hold up to 2 quarts of liquid
Size	6.5 qt	--	12 qt	16.5" Seat Depth	16.5" Seat Depth	14.5" Seat Depth	TBD
Metal Handle	✓	--	✓	--	--	✓	--
Cover/Lid	✓	--	✓	✓	✓	✓	--
UOM	12/CS	24/CS	12/CS	1/BX	1/BX	1/BX	7/BX
HPCP	E0167	--	E0167	--	--	--	--

LEGEND

ECOMMERCE
Available to order online

RETAIL
Available in large retail chain stores

PROVIDER
Available in independent DME provider locations

GOOD ★
Includes your most basic models and features that meet the customer minimum requirements and provide a good value.

BETTER ★★
Includes similar design and features from our Good models, our Better models offer supplementary features that give the consumer a better user experience.

BEST ★★★
Includes similar design and features from our Good and Better models, our Best models offers unique high-end features providing consumers with a top level experience.