

12440KD-1 BATH SEAT

- 1. Arms (pr)- 12445KDARMS
- 2. Seat- 12440KD-S
- 3. Leg Extension (4 ea)- 12440KD-LE
- 4. Leg (1 ea)- 12440KD-L
- 5. Tip (4 ea)- 12002T



SPECIFICATIONS

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight
12440KD-1	18.5"	16"	17"	13"	21.5"	16.5"	6.75 lbs.

[Diagrams](#)

[Home](#)