

Kanga Adult Folding Tilt-in-Space Wheelchair

WHEELED MOBILITY



The Kanga™ Adult Wheelchair is ideal for anyone who requires adjustability in the tilt-in-space market, whether in an adult rehab, capped rental or long term care setting.

ACCESSORIES



70° Footrests, Adjustable Angle Footplates



Adult Elevating Legrests



12" Rear Wheel Assembly



Contoured Chest Harness



Padded Center Pull Pelvic Belt



Folds for Easy Transport



Adjustable Width Crossbrace

BENEFITS

Adjustable Design

FEATURES

- Folding frame tilt-in-space wheelchair with 45° of tilt
- Frame width and depth adjusts without added components
- Adjustable height back canes and backrest angle
- Height-adjustable & removable "T" style armrests
- Width-adjustable, solid seat pan included
- 20" rear wheels included, 12" rear wheels optional

Built-In Growth

- Width-adjustable frame in 16, 18 or 20"

Quick & Easy Transport

- Removable swing-away foot hangers included
- Transit brackets and removable anti-tippers included
- Frame folds when not in use

Standard Accessories

- Seat Cushion/Backrest/Headrest/Headrest Hardware (sold separately), Angle Adjustable Footrest, Elevating Legrest, 12" Rear Wheel, Poziform™ Pelvic Belt, Chest Harness, and Ankle Positioners

Inspired
by drive

For more information, call us today!

877-224-0946 | www.inspiredbydrive.com