

# Marlin™ Bath Lift

BATHING



## ACCESSORIES



Head and Trunk Support Pads



Contoured Chest Harness



Swivel and Transfer Aid



Rapid Dry Positioning Belt



The compact and easy-to-use Marlin™ Bath Lift features a reclining back and rechargeable battery-operated lift mechanism that is easy to transport or store when not in use.

### BENEFITS

### FEATURES

#### User Accommodation

- Seat height adjusts from 2.5–19" and back reclines from 10° off vertical to 50° at its lowest position
- Multiple positioning accessories available for postural support

#### Compact & Versatile

- Weighs 20.5 lb. and fits most bathtubs
- Supports up to 300 lb. and includes side flaps for easy lateral transfers

#### Easy Operation

- Waterproof, simple-to-use hand control safely and gradually lowers user into bathtub
- Safety mechanism prevents lift from lowering if battery charge is not sufficient to raise it back up

#### Travel-Friendly

- Tool-free assembly and installation
- Disassembles for easy storage or transport

For more information, call us today!

877-224-0946 | [www.inspiredbydrive.com](http://www.inspiredbydrive.com)