

Moxie™ GT Gait Trainer

WALKING AID



ACCESSORIES




Hand Loops



Forearm Platforms



Handgrips



Guide Bar



Ankle Prompts



Soft Seat Harness



Trunk Support



Hip Positioner & Pad

The Moxie™ GT offers a tremendous amount of versatility to meet the needs of most users. Height-adjustable and available in three sizes, the Moxie GT can be used in either a posterior or anterior position.

AVAILABLE COLORS



Warrior Blue
(Small & Medium Only)



Magic Green
(Small & Medium Only)



Sword Gray
(Large Only)

BENEFITS

Customizable Frame Options

Safe Ambulation

FEATURES

- Lightweight and durable aluminum and steel frame construction
- Medium and large frames fold easily for storage and transportation
- Three frame sizes with 1" height adjustments allow for range of growth in each size
- Large, 5.5" rubber multi-function front and rear casters smoothly navigate indoors and outdoors
- Wheel locks, swivel locks, one-directional controls and drag control for added safety
- Accommodates anterior or posterior configurations



For more information, call us today!

877-224-0946 | www.inspiredbydrive.com