



Item #	GOOD		BETTER	
	SEAT2GO		FIRST CLASS SCHOOL CHAIRS	
	CRS 2000	CRS 3000	FC 2000N	FC 4000N
Product Name	Seat2Go, Small	Seat2Go, Medium	First Class School Chair, Small	First Class School Chair, Large
Sales Channel				
Colors	Blue Exterior w/ Black Seat	Blue Exterior w/ Black Seat	Blue w/ Red Legs	Blue w/ Red Legs
Material	Plastic, Wood, Cloth	Plastic, Wood, Cloth	Plastic, Steel, Aluminum	Plastic, Steel, Aluminum
Hip Width	8"	12"	12"	16"
Seat Depth	9"	12"	9.5"-14.5"	13"-17"
Head Height	12"	14"	15"	17"
Product Weight	5 lb	7 lb	16.4 lb	22.8 lb
Case Pack	1/cs	1/cs	1/cs	1/cs
HCPC	E0143	E0143	E1399	E1399
Warranty	Limited 1 Year	Limited 1 Year	Limited 1 Year	Limited 1 Year

Item #	ACCESSORIES FOR SEAT2GO		ACCESSORIES FOR FIRST CLASS SCHOOL CHAIRS		
	CRS 8000	CRS 8025	FC 2029N	FC 4029N	FC 8000N
Product Name	Headrest for Seat2Go & Back2Go	Abductor for Seat2Go Small/Medium	Legs w/ Casters for FC 2000N	Legs w/ Casters for FC 4000N	Multi-Axis Headrest
Sales Channel					
Product Weight	0.7 lb	0.4 lb	5.6 lb	6.3 lb	2.9 lb
Case Pack	1/cs	1/cs	1 set/cs	1 set/cs	1/cs
HCPC	E1399	E1399	E1399	E1399	E1399
Warranty	Limited 1 Year	Limited 1 Year	Limited 1 Year	Limited 1 Year	Limited 1 Year

Item #	ACCESSORIES FOR FIRST CLASS SCHOOL CHAIRS					
	FC 8001	FC 2024	FC 4024	FC 8025N	FC 8025 LN	FC 8027N
Product Name	Push Handles	Dining Tray for FC 2000N	Dining Tray for FC 4000N	Support Kit, Small, for FC 2000N	Support Kit, Large, for FC 4000N	Hip Guide
Sales Channel						
Product Weight	0.7 lb	6.9 lb	8.1 lb	4.5 lb	4.6 lb	4.8 lb
Case Pack	1 pair/cs	1/cs	1/cs	1 set/cs	1 set/cs	1 pair/cs
HCPC	E1399	E1399	E1399	E1399	E1399	E1399
Warranty	Limited 1 Year	Limited 1 Year	Limited 1 Year	Limited 1 Year	Limited 1 Year	Limited 1 Year

Item #	ACCESSORIES FOR FIRST CLASS SCHOOL CHAIRS			
	FC 2028	FC 4028	FC 2030N	FC 4030N
Product Name	Footrest for FC 2000	Footrest for FC 4000	Anti Tippers for FC 2000N	Anti Tippers for FC 4000N
Sales Channel				
Product Weight	3 lb	3.8 lb	0.9 lb	1.1 lb
Case Pack	1 pair/cs	1 pair/cs	1 pair/cs	1 pair/cs
HCPC	E1399	E1399	E1399	E1399
Warranty	Limited 1 Year	Limited 1 Year	Limited 1 Year	Limited 1 Year

LEGEND

PROVIDER
Available in independent DME provider locations

ECOMMERCE
Available to order online



BEST P PODS								
Item #	APP-2000S-MP, APP-2000S-PB	APP-2000S-DX	APP-3000M-MP, APP-3000M-PB	APP-3000M-DX	APP-4000L-MP, APP-4000L-WB, APP-4000L-JP	APP-4000L-DX	APP-5000XL-MP, APP-5000XL-WB, APP-5000XL-JP	APP-5000XL-DX
Product Name	P Pod, Small		P Pod, Medium		P Pod, Large		P Pod, Extra Large	
Sales Channel								
Colors	Mermaid Pink, Pool Blue		Starfish Red		Mermaid Pink, Whale Blue, Jellyfish Purple		Starfish Red	
Water Proof, Wipeable Dartex	--		✓		--		✓	
Material	Fabric, Foam, Plastic		Fabric, Foam, Plastic		Fabric, Foam, Plastic		Fabric, Foam, Plastic	
Hip Width	10.25"		12.25"		14.25"		16.5"	
Trunk Width	10"		10.5"		12.5"		13.5"	
Seat Depth	11.5"		13"		13.75"		17"	
Axilla Height	11.5"		12.5"		15"		16.5"	
Shoulder Height	16.5"		17.25"		20"		21.75"	
Head Height	22"		24.5"		26"		28.75"	
Product Weight	28 lb		28 lb		32 lb		32 lb	
Case Pack	1/cs		1/cs		1/cs		1/cs	
Warranty	Limited 1 Year		Limited 1 Year		Limited 1 Year		Limited 1 Year	

ACCESSORIES FOR P POD									
Item #	APP-A4500M	APP-A4500L	APP-A4501M	APP-A9001M	APP-A9002M	APP-A9002L	APP-A7300M	APP-A5330S, APP-A5330M, APP-A5330L, APP-A5330XL	APP-A3800S, APP-A3800M, APP-A3800L, APP-A3800XL
Product Name	Freestanding Tray, Medium	Freestanding Tray, Large	Lap Tray, Medium	Activity Frame, Medium	Mobility Base, Medium	Mobility Base, Large	Foot Bolster, Medium	Lateral Bolster, Small, Medium, Large, Extra Large	Chest Harness, Small, Medium, Large, Extra Large
Sales Channel									
Product Weight	23.3 lb	24 lb	3 lb	12 lb	34 lb	44 lb	5 lb	1 lb	1 lb
Case Pack	1/cs	1/cs	1/cs	1/cs	1/cs	1/cs	1/cs	1/cs	1/cs
HCPC	E1399	E1399	E1399	E1399	E1399	E1399	E1399	E1399	E1399
Warranty	Limited 1 Year	Limited 1 Year	Limited 1 Year	Limited 1 Year	Limited 1 Year	Limited 1 Year	Limited 1 Year	Limited 1 Year	Limited 1 Year

LEGEND

PROVIDER
Available in independent DME provider locations

ECOMMERCE
Available to order online