

# Spirit Plus™

## CAR SEAT



The Spirit Plus™ Car Seat provides positioning support for users with mild to complex positioning needs. The Spirit Plus is the only Car Seat on the market today offering swing-away padded trunk supports, hip supports and abductor.

## ACCESSORIES



Incontinence Covers



3-inch Extensor Thrust Wedge



Buckle Guard



Retainer Clip Guard

## AVAILABLE COLORS



Convertible Pink (CP)



Go Kart Green (GG)



Rally Blue (RB)



Roadster Red (RR)



Speedway Gray (SG)

## BENEFITS

### User Need Accommodation

## FEATURES

- Two sizes with seat depth extender (12–16") allow for built-in growth
- Height & weight capacities support users up to 66" in height and 130 lb. in weight

### Optimal Support

- Recline bar adds 25° of extra tilt for added comfort
- Trunk & hip supports and abductor help maintain proper positioning and postural support

### User Safety & Security

- Durable, machine-washable padded cover
- Memory foam back/seat padding reduces sores and fatigue
- Open shell design with low-profile sides and swing-away supports for easy transfers

### Accessory & Expansion Options

- Optional incontinence covers, anti-thrust wedge and buckle guard available



For more information, call us today!

877-224-0946 | [www.inspiredbydrive.com](http://www.inspiredbydrive.com)