

Spirit™ Spica

CAR SEAT



The Inspired by Drive Spirit™ Spica Car Seat features a low and wide profile with positioning wedges to safely transport users with spica casts. The Spirit Spica seat with wipeable cover accommodates a wide range of user heights and weights to meet individual needs.



ACCESSORIES



Positioning Wedge
(2 included)



Memory foam
support



Wipeable,
low maintenance
upholstery



BENEFITS

Built-In Safety & Security

FEATURES

- 5-point adjustable seat harness with included tether strap provides added safety
- Easily installed in most vehicles
- Two seat depth options with a wide weight and height range makes this seat a great choice for multiple clients

Adjusts to User

- Seat depth options of 12" and 14"
- Accommodates users up to 66" tall and has a 25-130 lb weight capacity
- Adjustability and seat configurations can be obtained for different sized clients and Spica casts

Comfort & Convenience

- High quality memory foam on the back and seating surface provides a comfortable environment
- Padded, wipeable Stealth Black upholstery
- 2 wipeable foam seat wedges included standard
- 1 year limited warranty

For more information, call us today!

877-224-0946 | www.inspiredbydrive.com