


Versa™ Transfer System

BATHING



The lightweight and multi-functional Versa™ Transfer System is most suitable for users who require minimal positioning support. Constructed of high-grade aluminum, stainless steel and padded polyurethane, the Versa can be utilized as a bath transfer system, roll-in shower chair or bedside commode. The Versa features a polyurethane seat with opening and a tub base that connects to the rolling base for safe transfers.

ACCESSORIES



Fixed-Length Bridge



Adjustable-Length Bridge



Rapid Dry Positioning Belt



Lap Belt



BENEFITS

Versatile Design

FEATURES

- Compact, 18" x 18", tub base transfer frame fits inside most bathtubs
- Multiple Bridge Options to fit your bathroom
 - Standard removable 13" bridge
 - Optional fixed bridge in 9"
 - Optional adjustable bridges from 13½"–17¾" or 17½"–21¼" lengths

Multi-Functional

- Allows incontinent users to bathe and shower comfortably
- Opening for easy toileting or bathing access

User Accommodation

- Use for adolescents, teens or adults
- Supports up to 300 lb.

Built-In Safety

- Locking system limits movement and slipping

Available Accessories

- Headrest, commode pail, footrests and locking casters come standard
- Fixed and Adjustable Bridges, Lap Belt and Rapid Dry Positioning Belt options available

For more information, call us today!

877-224-0946 | www.inspiredbydrive.com