



# Crutches | Good, Better, Best Options



Product Description	★ GOOD ★			★★ BETTER ★★		BARIATRIC	★ GOOD ★	BARIATRIC	★★★ SPECIALTY ★★★
	Standard Underarm Crutches						Standard Forearm Crutches		Platform Crutches
Item Numbers	Standard Crutch	Standard Crutch	Standard Crutch	EZ Adjust Crutch with Euro-Style Clip	Universal EZ Adjust Crutch with Euro-Style Clip	Bariatric Standard Crutch	Forearm Crutch	Bariatric Forearm Crutch	M+D Comfort Crutch™
	10401-8 (Youth) 10400-8 (Adult) 10402-8 (Tall Adult)	10401-1 (Youth) RTL10400 (Adult) RTL10402 (Tall Adult)	721-790 (Youth) 721-785 (Adult) 721-780 (Tall Adult)	10431-8 (Youth) 10430-8 (Adult) 10432-8 (Tall Adult)	RTL10433	10406 (Adult) 10408 (Tall Adult)	10403 (Adult) 10405 (Tall Adult)	10403HD	MD20041 (Assembled) MD30002 (Non-Assembled)
Sales Channel									
Product Material:	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum	Steel	Steel	Steel	Aluminum
Handle Height	Youth	23" - 37"	17.25" - 31.25"	22.5" - 36.5"		--	--		
	Adult	23" - 37"	24.75" - 38.5"	24.5" - 38.5"	24.5" - 47"	24.5" - 38.5"	28" - 37"	28" - 37"	39" - 52"
	Tall Adult	29" - 43"	29.25" - 43.25"	31" - 45"		31.5" - 45.5"	33" - 42"		
Underarm Height	Youth	37" - 46"	54" - 62"	36" - 44"		--	--	--	--
	Adult	45" - 53"	62" - 70"	45" - 53"	38.5" - 62.5"	46" - 54"	--	--	--
	Tall Adult	53" - 61"	70" - 78"	52.5" - 60.5"		55" - 63"	--	--	--
User Height	Youth	4'6" - 5'2"	4'6" - 5'2"	4'6" - 5'2"		--	--		
	Adult	5'2" - 5'10"	5'2" - 5'10"	5'2" - 5'10"	4'6" - 6'6"	5'0" - 6'2"	5'0" - 6'2"	5'0" - 6'2"	4'11" - 6'8"
	Tall Adult	5'10" - 6'6"	5'10" - 6'6"	5'10" - 6'6"		5'10" - 6'6"	5'10" - 6'6"		
Weight Capacity	350 lbs.	350 lbs.	300 lbs.	300 lbs.	300 lbs.	500 lbs.	300 lbs.	500 lbs.	300 lbs.
Color	Silver	Silver	Silver	Silver	Silver	Grey	Silver	Silver	Silver/Black
Piggy Back Labels	✓	--	--	--	--	--	--	--	--
Unit of Measure	8 Pairs/Case	1 Pair, 8 Pairs/Case	1 Pair, 8 Pairs/Case	8 Pairs/Case	1 Pair, 8 Pairs/Case	1 Pair/Case	1 Pair/Case	1 Pair/Case	MD20041: 2 Pairs/Box MD30002: 1 Pair/Box, 3 Pairs/Box
HCPC	E0114	E0114	E0114	E0114	E0114	E0116	E0110   E0111	E0110   E0111	--
Warranty	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited One Year from Receiving Date

Product Description	Accessories						
	Crutch Handgrips (Closed Style)	Crutch Underarm Cushions	Crutch Tip	Crutch Tip	Crutch Tip	Crutch Pillows (2 underarm; 2 handgrip covers)	Crutch Accessory Kit (2 underarm; 2 handgrip; 2 tips)
Item Numbers	RTL10453	RTL10446	RTL10439B	RTL10374BKB	10439-8	RTL10355	RTL10385/RTL10395
Sales Channel							
Diameter	--	--	7/8"	7/8"	7/8"	--	Underarm Cushions: 8 1/2" x 7/8" Handgrips: 4" x 3/4", Tips: 7/8"
Color	Grey	Grey	Grey	Black	Grey	Grey	Tan (RTL10385) Grey (RTL10395)
Unit of Measure	1 Pair/Case, 8 Pairs/Case	1 Pair/Case, 8 Pairs/Case	2/Blister, 10/Case	2/Blister, 10/Case	8/Case	1/Box, 8 Boxes/Case	
HCPC	A4636	A4635	A4637	A4637	A4637	--	A4635   A4636   A4637
For Use With	10430-8 10431-8 10432-8 10400-8 10401-8 10402-8	10430-8 10431-8 10432-8 10400-8 10401-8 10402-8	10430-8 10431-8 10432-8 10400-8 10401-8 10402-8	10430-8 10431-8 10432-8 10400-8 10401-8 10402-8	10430-8 10431-8 10432-8 10400-8 10401-8 10402-8	10403 10430-8 10431-8 10432-8 10400-8 10401-8 10402-8	Universal - Standard Crutches  Universal - Standard Crutches

### LEGEND

**E-COMMERCE**  
Available to order online

**RETAIL**  
Available in large retail chain stores

**PROVIDER**  
Available in independent DME provider locations

**GOOD ★**  
Good models are the base within the class, comparable to the competition of the same level. These meet the basic requirements with less refinements and minimal extras vs. the "BETTER" and "BEST" classifications. They tend to be entry level in the market with good value.

**BETTER ★★**  
Mid-level models within the class, offer more features and benefits to give the end user a better user experience. These tend to be lighter weight than the good models and have more options available. They may have a higher user weight limit as well.

**BEST ★★★**  
Best models combine similar design features from our "Good" and "Better" models but offer additional features that pushes them to the top ranking. They may be very light or could come with more options such as anti tipppers or flip back arms included. They tend to be made from lighter materials without compromising strength to give the end user the experience they are looking for.