



Power Wheelchair Type	★★ BETTER ★★			★★★★ BEST ★★★★★		★★★★ BEST ★★★★★		★★★★ BEST ★★★★★		
	Standard Power Wheelchair			Lightweight Power Wheelchair - Rear Wheel	Heavy Duty Power Wheelchair		Folding Power Wheelchair			
Item #	TITANAXS-18CS	TITANAXS-20CS	TITANAXS-22CS	TITANLTE-18FS	2850HD-22	2850HD-24	CPN18FBAN	CPN20FBAN	CPN22FBAN	
Product Name	Titan AXS	Titan AXS	Titan AXS	Titan LTE	Trident HD	Trident HD	Cirrus Plus LT	Cirrus Plus LT	Cirrus Plus LT	
Image										
Sales Channel										
Key Differences	<i>K0823 coded mid-wheel drive power wheelchair, interchangeable color panels and 18" captain seat.</i>	<i>K0823 coded mid-wheel drive power wheelchair, interchangeable color panels and 20" captain seat.</i>	<i>K0823 coded mid-wheel drive power wheelchair, interchangeable color panels and 22" captain seat.</i>	<i>K0821 coded rear wheel drive power wheelchair with 18" stadium seat. Disassembles for transport & storage.</i>	<i>K0825 coded front wheel drive heavy duty wheelchair with a 450 lb. weight capacity and 22" captain seat.</i>	<i>K0825 coded front wheel drive heavy duty wheelchair with a 450 lb. weight capacity and 24" captain seat.</i>	<i>K0822 coded rear wheel drive power wheelchair. Its foldable, lightweight frame is ideal for portability, and has an 18" sling seat.</i>	<i>K0822 coded rear wheel drive power wheelchair. Its foldable, lightweight frame is ideal for portability, and has a 20" sling seat.</i>	<i>K0822 coded rear wheel drive power wheelchair. Its foldable, lightweight frame is ideal for portability, and has a 22" sling seat.</i>	
Drive Type	Mid-Wheel	Mid-Wheel	Mid-Wheel	Rear Wheel	Front Wheel	Front Wheel	Rear Wheel	Rear Wheel	Rear Wheel	
Weight Cap	300 lb	300 lb	300 lb	300 lb	450 lb	450 lb	300 lb	300 lb	330 lb	
Top Speed	4 mph	4 mph	4 mph	4.2 mph	4 mph	4 mph	4 mph	4 mph	4 mph	
Maximum Range	19 miles	19 miles	19 miles	8 miles	15 miles	15 miles	12 miles	12 miles	12 miles	
Turning Radius	20"	20"	20"	27"	29"	29"	30"	32"	33"	
Climbing Angle	6°	6°	6°	6°	6°	6°	6°	6°	6°	
Ground Clearance	2.5"	2.5"	2.5"	2.5"	2.5"	2.5"	4"	4"	4"	
Dimensions (L x W x H)	39.25" x 23.5" x 44"	39.25" x 23.5" x 44"	39.25" x 23.5" x 44"	34" x 25.75" x 36"	48" x 25.75" x 52"	48" x 25.75" x 52"	41.7" x 26.2" x 35.4"	41.7" x 28" x 35.4"	41.7" x 30" x 35.4"	
Colors / Shrouds	Red & Blue	Red & Blue	Red & Blue	Red	Black	Black	Black	Black	Black	
Seat Type	Captain	Captain	Captain	Stadium	Captain	Captain	Sling	Sling	Sling	
Seat Dimensions	18" x 18"	20" x 18"	22" x 20"	18" x 18"	22"	24"	18" x 16.5"	20" x 16.5"	22" x 16.5"	
Controller	Dynamic Linx 40A	Dynamic Linx 40A	Dynamic Linx 40A	Dynamic Linx 40A	PG VR2 70A	PG VR2 70A	Micon 50A	Micon 50A	Micon 50A	
Motor(s)	24V x 320W x 4600 rpm	24V x 320W x 4600 rpm	24V x 320W x 4600 rpm	200W x 3900RPM	24V x 450W x 2700 rpm	24V x 450W x 2700 rpm	24V x 340W x 4500 rpm	24V x 340W x 4500 rpm	24V x 340W x 4500 rpm	
Batteries	12V x 33AH (2)	12V x 33AH (2)	12V x 33AH (2)	12V x 20AH (2)	12V x 50AH (2)	12V x 50AH (2)	12V x 33AH (2)	12V x 33AH (2)	12V x 33AH (2)	
Charger	4A	4A	4A	3A	5A	5A	4A	4A	4A	
Product Weight	171 lb	171 lb	171 lb	122 lb	212 lb	212 lb	146 lb	146 lb	156 lb	
Drive Wheels	10" x 3"	10" x 3"	10" x 3"	9" x 3"	14" x 3"	14" x 3"	12.5" x 2.2"	12.5" x 2.2"	12.5" x 2.2"	
Caster Wheels	6" x 2"	6" x 2"	6" x 2"	8" x 2"	10" x 3"	10" x 3"	7.5" x 1.75"	7.5" x 1.75"	7.5" x 1.75"	
HCP	K0823	K0823	K0823	K0821	K0825	K0825	K0822	K0822	K0822	
Warranty	Limited Lifetime Warranty on Frame, 24 Month Limited Warranty on Electronics, 12 Month Limited Warranty on Batteries	Limited Lifetime Warranty on Frame, 24 Month Limited Warranty on Electronics, 12 Month Limited Warranty on Batteries	Limited Lifetime Warranty on Frame, 24 Month Limited Warranty on Electronics, 12 Month Limited Warranty on Batteries	Limited Lifetime Warranty on Frame, 24 Month Limited Warranty on Electronics, 12 Month Limited Warranty on Batteries	Limited Lifetime Warranty on Frame, 24 Month Limited Warranty on Electronics, 12 Month Limited Warranty on Batteries	Limited Lifetime Warranty on Frame, 24 Month Limited Warranty on Electronics, 12 Month Limited Warranty on Batteries	Limited Lifetime Warranty on Frame, 24 Month Limited Warranty on Electronics, 12 Month Limited Warranty on Batteries	Limited Lifetime Warranty on Frame, 24 Month Limited Warranty on Electronics, 12 Month Limited Warranty on Batteries	Limited Lifetime Warranty on Frame, 24 Month Limited Warranty on Electronics, 12 Month Limited Warranty on Batteries	Limited Lifetime Warranty on Frame, 24 Month Limited Warranty on Electronics, 12 Month Limited Warranty on Batteries

LEGEND



PROVIDER
Available in independent DME provider locations



ECOMMERCE
Available to order online

GOOD ★

Good models are the base within the class, comparable to the competition of the same level. These meet the basic requirements with less refinements and minimal extras vs. the "Better" and "Best" classifications. They tend to be entry level in the market with good value.

BETTER ★★

Mid-level models within the class, offer more features and benefits to give the consumer a better user experience. These tend to have larger, more padded seats, longer travel ranges, possibly a tighter turning radius and a generally more refined feel.

BEST ★★★

Best models combine similar design features from our "Good" and "Better" models but offer additional features that pushes them to the top ranking. They may separate for easy transport, fold, or have higher weight capacity. Sometimes these features come at a cost to other benefits such as tighter turning radius or seat style. They tend to include the most desirable features for the customer and work to give them the experience they are looking for.