

# NITRO

## Rollator



RTL10266

RTL10266WT

RTL10266BK



### FEATURES

- Attractive, Euro-style design.
- Brake cable inside frame for added safety.
- Handle height easily adjusts with unique push button.
- Back support height easily adjusts with tool-free thumb screw.
- Caster fork design enhances turning radius.
- Large 10" front casters allow optimal steering and rolling comfort.
- Easily folds with one hand to ultra-compact size for storage. Handle is located on seat.
- Plastic clip holds the Nitro closed when frame is folded.
- Cross-brace design allows for side-to-side folding and added stability.
- Seat is durable and comfortable.
- Handsome and removable zippered storage bag. Unique bag attachment keeps bag securely in place when rollator is open or folded.
- Lightweight, aluminum frame.
- Comes in retail box.
- Weight capacity: 300 lbs.

For more information, contact Customer Service today.  
t: 800-371-2266 • f: 516-998-4601 • [www.drivemedical.com](http://www.drivemedical.com)

*for the road ahead™*  
**drive**