



Item Number(s)	★ GOOD ★		★★★ BEST ★★★	★ GOOD ★			★★ BETTER ★★			
	3 Wheel			Standard 4-Wheel						
Product Name	RTL199	10289RD, 10289BL	RTL10266TWHL	R800RD, R800BL	R800KD-RD, R800KD-BL	RTL800KD-BK	10257RD-1, 10257BL-1	R726RD, R726BL, R726BK	RTL10261RD, RTL10261BL	R728RD, R728BL
Product Name	Winnie Lite Supreme/Go Lite 3 Wheel Aluminum Rollator	Deluxe 3 Wheel Aluminum Rollator, 7.5" Casters	Nitro 3-Wheel Rollator	Steel Rollator, 6" Casters	Steel Rollator, Knockdown, 6" Casters	Travel Rollator, 6" Casters	Steel Rollator, 7.5" Casters	Aluminum Rollator, 6" Casters	Aluminum Rollator, 6" Casters	Aluminum Rollator, 7.5" Casters
Sales Channel										
Key Differences	Easy, one hand folding 3 wheel aluminum rollator.	Easy, one hand folding 3 wheel aluminum rollator with tray and folding wire basket.	Modern, patented design offering a premium look. Hidden brake cable for added safety. Comes fully assembled, and easily folds with one hand.	Standard steel rollator with under-seat fabric pouch for storage.	Knockdown standard steel rollator with under seat fabric pouch for storage. Requires tool-free assembly of all wheels and handbrakes.	Retail packaged opening price point rollator with under-seat fabric pouch for storage. Requires tool-free assembly of wheels and handbrakes.	Standard steel rollator with under-seat fabric pouch for storage.	Lightweight aluminum rollator with under-seat fabric pouch for storage.	Adjustable height aluminum rollator with removeable padded backrest and under-seat storage.	Aluminum rollator with compact design. Removeable padded backrest and under-seat storage.
Assembly Required	✓	✓	--	✓	✓	✓	✓	✓	✓	✓
Foldable	✓	✓	✓	✓	✓	✓	✓	--	--	--
Material	Aluminum	Aluminum	Aluminum	Steel	Steel	Steel	Steel	Aluminum	Aluminum	Aluminum
Storage	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Wheel Size	7.5"	7.5"	9" Front / 8" Back	6"	6"	6"	7.5"	6"	6"	7.5"
Handle Height	30.75" - 37.5"	31" - 38"	32" - 37"	31.25" - 36.25"	30.75" - 36"	30" - 36"	31" - 35"	32" - 37"	29.5" - 38"	33"-38"
Unit Dimensions (W x D)	27" x 22.75"	25" x 24"	22" x 29"	24" x 24.75"	24" x 23.25"	22" x 24"	23.5" x 25.5"	24.5" x 27.5"	24" x 27"	24" x 28"
Seat Dimensions (W x D)	--	--	--	12" x 12"	12" x 12"	14" x 12"	12" x 12"	14" x 12"	13" x 12"	14" x 12"
Seat Height	--	--	--	20.25"	20.25"	20"	21"	20"	18"-22"	23"
Width Between Arms	14"	18.5"	14"	19.25"	19.25"	18"	17"	17"	17"	17"
Product Weight	11 lbs	14 lbs	12.5 lbs	16.5 lbs	17 lbs	17 lbs	18.5 lbs	13 lbs	14 lbs	14 lbs
Shipping Weight	15 lbs	16 lbs	16.5 lbs	19.25 lbs	20 lbs	19 lbs	22.5 lbs	16 lbs	16 lbs	17 lbs
Weight Capacity	300 lbs	300 lbs	300 lbs	300 lbs	350 lbs	350 lbs	300 lbs	300 lbs	300 lbs	300 lbs
Available Colors	Tan Plaid	Flame Red, Flame Blue	Red	Red, Blue	Red, Blue	Black	Red, Blue	Red, Blue, Black	Red, Blue	Red, Blue
Case Pack	1/CS	1/CS	1/CS	1/CS	1/CS	1/CS	1/CS	1/CS	1/CS	1/CS
Retail Packaged	✓	--	--	--	--	✓	--	--	✓	--
Retail Box Dimensions	20" x 9.5" x 26.75"	--	--	--	--	22" x 7" x 18"	--	--	20.25" x 6" x 21.5"	--
HCPC	E0143	E0143	E0143	E0143 + E0156	E0143 + E0156	--	E0143 + E0156	E0143 + E0156	E0143 + E0156	E0143 + E0156
Warranty	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime

LEGEND



E-COMMERCE
Available to order online



RETAIL
Available in large retail chain stores



PROVIDER
Available in independent DME provider locations

GOOD ★

Includes your most basic models and features that meet the customer minimum requirements and provide a good value.

BETTER ★★

Includes similar design and features from our Good models, our Better models offer supplementary features that give the consumer a better user experience.

BEST ★★★

Includes similar design and features from our Good and Better models, our Best models offers unique high-end features providing consumers with a top level experience.



Rollators | Good, Better, Best Options



	★ GOOD ★			★★ BETTER ★★				★★★ BEST ★★★
	Nitro							
Item Number(s)	RTL10266, RTL10266BK, RTL10266WT	RTL10266BK-H	RTL10266-T, RTL10266BK-T	RTL10266CH-HS, RTL10266BL-HS	102662RD, 102662BL, 102662BK	102662RD-T, 102662BL-T, 102662BK-T	102662RD-H, 102662BL-H, 102662BK-H	RTL10266CF
Product Name	Nitro Aluminum Rollator Standard Height	Nitro Aluminum Rollator Hemi Height	Nitro Aluminum Rollator Tall Height	Nitro Deluxe Rollator	Nitro Sprint Rollator Standard Height	Nitro Sprint Rollator Tall Height	Nitro Sprint Rollator Hemi Height	Nitro Elite Carbon Fiber Rollator
Sales Channel								
Key Differences	Modern, patented design, offering a premium look. Hidden brake cable for added safety.	Modern, patented design, offering a premium look. Hidden brake cable for added safety.	Modern, patented design, offering a premium look. Hidden brake cable for added safety.	Modern, patented design, offering a premium look. Hidden brake cable for added safety. Featuring a deluxe storage bag and durable, comfortable padded seat.	Modern, patented design, offering a premium look. Hidden brake cable for added safety. Features windowed backrest, patented slow down brake, and accessory ports	Modern, patented design, offering a premium look. Hidden brake cable for added safety. Features windowed backrest, patented slow down brake, and accessory ports	Modern, patented design, offering a premium look. Hidden brake cable for added safety. Features windowed backrest, patented slow down brake, and accessory ports	Modern, patented design, offering a premium look. Hidden brake cable for added safety. Made of super lightweight carbon fiber and featuring a deluxe storage bag.
Assembly Required	--	--	--	--	--	--	--	--
Foldable	✓	✓	✓	✓	✓	✓	✓	✓
Material	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum	Carbon Fiber
Storage	✓	✓	✓	✓	✓	✓	✓	✓
Wheel Size	10" Front	10" Front	10" Front	10" Front	10" Front	10" Front	10" Front	10" Front
Handle Height	33.5" - 38.25"	28" - 31"	36" - 41"	33.5" - 38.5"	30" - 36.5"	32.0" - 39.4"	29.1" - 33.4"	32.5" - 37.5"
Unit Dimensions (W x D)	23" x 27.75"	23" x 27.75"	23" x 27.75"	24" x 28"	24" x 27.5"	24" x 27.5"	24" x 27.5"	23.25" x 38"
Seat Dimensions (W x D)	18" x 10"	18" x 8.5"	18" x 8.5"	18" x 10"	18" x 10"	18" x 10"	18" x 10"	17.75" x 8.5"
Seat Height	20.5"	18"	23.5"	20"	21"	23"	18"	20.75"
Width Between Arms	18"	18"	18"	18"	19"	19"	19"	18"
Product Weight	18 lbs	17.5 lbs	18 lbs	22 lbs	19 lbs	19.5 lbs	19 lbs	12 lbs (Excluding Bag)
Shipping Weight	21.5 lbs	21 lbs	22 lbs	25 lbs	23.1 lbs	23.9 lbs	23 lbs	19 lbs
Weight Capacity	300 lbs	300 lbs	300 lbs	300 lbs	350 lbs	350 lbs	350 lbs	300 lbs
Available Colors	Red, Black, White	Black	Red, Black	Champagne, Blue	Red, Blue, Black	Red, Blue, Black	Red, Blue, Black	Black Carbon Fiber
Case Pack	1/CS	1/CS	1/CS	1/CS	1/CS	1/CS	1/CS	1/CS
Retail Packaged	--	--	--	✓	--	--	--	✓
Retail Box Dimensions	--	--	--	28.5" x 10.5" x 40"	--	--	--	28.75" x 10.25" x 34"
HCPC	E0143 + E0156	E0143 + E0156	E0143 + E0156	E0143 + E0156	E0143 + E0156	E0143 + E0156	E0143 + E0156	E0143 + E0156
Warranty	Frame: Limited Lifetime Brake Handle & Housing: 5 Years	Limited Lifetime	Frame: Limited Lifetime Brake Handle & Housing: 5 Years	Frame: Limited Lifetime Brake Handle & Housing: 5 Years	Frame: Limited Lifetime, Brake Handle & Housing: 5 Years	Frame: Limited Lifetime, Brake Handle & Housing: 5 Years	Frame: Limited Lifetime, Brake Handle & Housing: 5 Years	Frame: Limited Lifetime Brake Handle & Housing: 5 Years

LEGEND



E-COMMERCE
Available to order online



RETAIL
Available in large retail chain stores



PROVIDER
Available in independent DME provider locations

GOOD ★

Includes your most basic models and features that meet the customer minimum requirements and provide a good value.

BETTER ★★

Includes similar design and features from our Good models, our Better models offer supplementary features that give the consumer a better user experience.

BEST ★★★

Includes similar design and features from our Good and Better models, our Best models offers unique high-end features providing consumers with a top level experience.



Rollators | Good, Better, Best Options



	★ GOOD ★	★★ BETTER ★★	★★★ BEST ★★★	★★ BETTER ★★	★★★ BEST ★★★	UPRIGHT
	Bariatric			Combo		Elevate
Item Number(s)	10215RD-1 10215BL-1	10266HD-BK	102662BK-HD	795BU, 795B, 795BK, 795RG, 795SL, 795PR	RTL10266DT	UPWKRBL10, UPWKRGY10
Product Name	Go-Lite Bariatric Rollator, 8" Casters	Nitro HD Rollator	Nitro Sprint HD Rollator	Duet Rollator/Transport Chair	Nitro Rollator/Transport Chair	Elevate Upright Walker Rollator
Sales Channel						
Key Differences	Standard bariatric steel rollator with padded seat and storage basket.	Modern, patented design, offering a premium look. Hidden brake cable for added safety.	Modern, patented design, offering a premium look. Hidden brake cable for added safety. Features windowed backrest, patented slow down brake, and accessory ports	Combines the features of a rollator and transport chair into one unit. Engage locking flip-up armrests and flip-down footrest to use as a transport chair.	Modern, patented design, offering a premium look. Hidden brake cable for added safety.	Modern design, large backrest, comfortable seat
Assembly Required	✓	--	--	✓	✓	--
Foldable	--	✓	✓	✓	✓	✓
Material	Steel	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum
Storage	✓	✓	✓	✓	✓	✓
Wheel Size	8"	10" Front	10" Front	8"	10" Front	10" / 8"
Handle Height	35.25" - 39.5"	36" - 41"	32.7" - 37.8"	31.5" - 37"	33" - 37"	39" - 44"
Unit Dimensions (W x D)	30.5" x 26.75"	28.25" x 28"	29" x 28"	25.5" x 26"	23.5" x 40"	25.8" x 31.5"
Seat Dimensions (W x D)	18" x 13"	23" x 10"	22.5" x 10"	13.5" x 13"	17.5" x 13"	18" x 10"
Seat Height	22"	23.5"	22"	21"	20.5"	21"
Width Between Arms	23"	23"	23"	17"	18"	29"
Product Weight	26 lbs	20.5 lbs	20 lbs	20 lbs	22 lbs	23 lbs
Shipping Weight	30.5 lbs	23.5 lbs	24.7 lbs	25 lbs	25 lbs	25 lbs
Weight Capacity	500 lbs	450 lbs	500 lbs	300 lbs	300 lbs	300 lbs
Available Colors	Red, Blue	Black	Black	Burgundy, Blue, Black, Rose Gold, Silver, Purple	Red	Blue / Grey
Case Pack	1/CS	1/CS	1/CS	1/CS	1/CS	1/CS
Retail Packaged	--	--	No	--	--	--
Retail Box Dimensions	--	--	N/A	--	--	--
HCPC	E0143 + E0156	E0143 + E0156	E0143 + E0156	E0143 + E0156	E0143 + E0156	--
Warranty	Limited Lifetime	Frame: Limited Lifetime Brake Handle & Housing: 5 Years	Frame: Limited Lifetime, Brake Handle & Housing: 5 Years	Limited Lifetime	Limited Lifetime	Limited lifetime

LEGEND



ECOMMERCE
Available to order online



RETAIL
Available in large retail chain stores



PROVIDER
Available in independent DME provider locations

GOOD ★

Includes your most basic models and features that meet the customer minimum requirements and provide a good value.

BETTER ★★

Includes similar design and features from our Good models, our Better models offer supplementary features that give the consumer a better user experience.

BEST ★★★

Includes similar design and features from our Good and Better models, our Best models offers unique high-end features providing consumers with a top level experience.



	★★★★ BEST ★★★★★			★★★★ BEST ★★★★★	★★★★ BEST ★★★★★		
	4-Wheel			Bariatric	Combo		
Item Number(s)	700-920, 700-920PR, 700-921, 700-922	700-927	700-917, 700-917BK	700-928, 700-929	700-935, 700-935BK, 700-935CR	700-938	700-932, 700-932PB
Product Name	eXcursion™ X20 Rollator Standard Height	eXcursion™ X18 Rollator Hemi Height	eXcursion™ X23 Rollator Tall Height	eXcursion™ XWD Rollator, Wide	Fusion™ Side-Folding Rollator & Transport Chair Standard Height	Fusion™ Side-Folding Rollator & Transport Chair Hemi Height	Fusion™ Side-Folding Rollator & Transport Chair Tall Height
Image							
Sales Channel							
Key Differences	Heavy-duty, ultra-lightweight sturdy frame. Patented handgrips that absorb vibrations while placing your hand at a natural angle. Padded seat and contoured backrest.	Heavy-duty, ultra-lightweight sturdy frame. Patented handgrips that absorb vibrations while placing your hand at a natural angle. Padded seat and contoured backrest.	Heavy-duty, ultra-lightweight sturdy frame. Patented handgrips that absorb vibrations while placing your hand at a natural angle. Padded seat and contoured backrest.	Heavy-duty, ultra-lightweight sturdy frame. Patented handgrips that absorb vibrations while placing your hand at a natural angle. Padded seat and contoured backrest.	Easily converts from a rollator to a transport chair with Swing n' Lock™ footrests. Safety belt is designed exclusively for the Fusion™ and installs securely onto the seat rails via the existing upholstery screws.	Easily converts from a rollator to a transport chair with Swing n' Lock™ footrests. Safety belt is designed exclusively for the Fusion™ and installs securely onto the seat rails via the existing upholstery screws.	Easily converts from a rollator to a transport chair with Swing n' Lock™ footrests. Safety belt is designed exclusively for the Fusion™ and installs securely onto the seat rails via the existing upholstery screws.
Assembly Required	✓	✓	✓	✓	✓	✓	✓
Foldable	✓	✓	✓	✓	✓	✓	✓
Material	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum
Storage	✓	✓	✓	✓	✓	✓	✓
Wheel Size	8"	8"	8"	8"	8"	8"	8"
Handle Height	30.5" - 37"	29.5" - 35"	33.5" - 39.5"	30" x 36.5"	30" - 36.5"	29.5" - 35"	33.5" - 39"
Unit Dimensions (W x D)	25" x 28.5"	25" x 28.5"	25" x 28.5"	27" x 28.5"	27" x 28.5"	27" x 28.5"	27" x 28.5"
Seat Dimensions (W x D)	17.5" x 10"	17.5" x 10"	17.5" x 10"	19.5" x 10"	19.5" x 10"	19.5" x 10"	19.5" x 10"
Seat Height	20"	18"	23"	20"	20"	18"	23"
Width Between Arms	17.5"	17.5"	17.5"	19.5"	19.5"	19.5"	19.5"
Product Weight	15.5 lbs	15.5 lbs	17 lbs	16.5 lbs	17 lbs	17 lbs	19 lbs
Shipping Weight	21 lbs	21 lbs	24 lbs	22 lbs	27 lbs	27 lbs	28 lbs
Weight Capacity	300 lbs	300 lbs	300 lbs	400 lbs	300 lbs	300 lbs	300 lbs
Available Colors	Pearl Black, Pacific Blue, Cranberry, Purple	Cranberry	Pacific Blue, Black	Pacific Blue, Pearl Black	Pearl Black, Pacific Blue, Cranberry	Cranberry	Pacific Blue, Pearl Black
Case Pack	1/CS	1/CS	1/CS	1/CS	1/CS	1/CS	1/CS
Retail Packaged	--	--	--	--	--	--	--
Retail Box Dimensions	--	--	--	--	--	--	--
HCPC	E0143 + E0156	E0143 + E0156	E0143 + E0156	E0143 + E0156	E0143 + E0156	E0143 + E0156	E0143 + E0156
Warranty	Limited lifetime	Limited lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime

LEGEND



ECOMMERCE
Available to order online



RETAIL
Available in large retail chain stores



PROVIDER
Available in independent DME provider locations

GOOD ★
Includes your most basic models and features that meet the customer minimum requirements and provide a good value.

BETTER ★★
Includes similar design and features from our Good models, our Better models offer supplementary features that give the consumer a better user experience.

BEST ★★★
Includes similar design and features from our Good and Better models, our Best models offers unique high-end features providing consumers with a top level experience.



	★ GOOD ★		★★ BETTER ★★			★★★ BEST ★★★	
	4-Wheel		4-Wheel			Combo	
Item Number(s)	700-913	700-957	700-959, 700-961	700-959E	700-979	700-977	793-943
Product Name	Fit Rolling Walker with Padded Seat	Portable Rolling Walker with Seat	Elite Rolling Walker with Seat	Elite Rolling Walker with Seat eCom Edition	Sidekick™ Side-Folding Rolling Walker with a Seat	Explore Rollator	Navigator™ Side-Folding Rolling Walker & Transport Chair
Key Differences	Lightweight design combined with adjustable seat & handle heights for a perfect fit. Patented ergonomic hand grips place your hands at a natural angle for comfort and better brake access.	Featuring patented handgrips that absorb vibrations while placing your hand at a natural angle. Padded seat & backrest with all-weather covering. Convenient under seat storage.	Featuring patented handgrips that absorb vibrations while placing your hand at a natural angle. Convenient under seat storage bag, with additional saddle bag for extra storage.	Featuring patented handgrips that absorb vibrations while placing your hand at a natural angle. Convenient under seat storage bag, with additional saddle bag for extra storage.	"Sleek design, heavy-duty yet ultra-light steel frame. Built in seat and contoured backrest provide a comfortable place to power up whenever you need a rest."	Patented height-adjustable cushion-grip handles that position your hands at a natural angle.	Lightweight design combined with unparalleled stability. Easily converts from a rollator to a full transport chair.
Assembly Required	✓	✓	✓	✓	✓	✓	✓
Foldable	✓	✓	✓	✓	✓	✓	✓
Material	Steel	Aluminum	Aluminum	Aluminum	Steel	Aluminum	Aluminum
Storage	✓	✓	✓	✓	✓	✓	✓
Wheel Size	6"	8"	8"	8"	8"	8"	8"
Handle Height	31" - 38"	32" - 38"	29.5" - 38"	29.5" - 38"	30.5" - 35"	30.5" - 37.0"	30" - 36.5"
Unit Dimensions (W x D)	24.5" x 27"	24.5" x 27"	24.5" x 27.5"	24.5" x 27.5"	23.2" x 28.5"	25" x 28.5"	27" x 28.5"
Seat Dimensions (W x D)	14" x 12"	14" x 12"	14" x 12"	14" x 12"	18" x 9"	17.5" x 12"	19.5" x 12"
Seat Height	20" - 24"	21.5"	19.5" - 22"	19.5" - 22"	20.5"	20"	20"
Width Between Arms	17"	17.5"	17"	17"	15.8"	17.5"	19.5"
Product Weight	15 lbs	18 lbs	18 lbs	18 lbs	17 lbs	15.5 lbs	17 lbs
Shipping Weight	17 lbs	21 lbs	21 lbs	21 lbs	21 lbs	25 lbs	27 lbs
Weight Capacity	250 lbs	300 lbs	300 lbs	300 lbs	250 lbs	300 lbs	300 lbs
Available Colors	Pacific Blue	Pacific Blue	Pacific Blue, Garnet Red	Pacific Blue	Blueberry	Pearl Black	Pacific Blue
Case Pack	1/CS	1/CS	1/CS	1/CS	1/CS	1/CS	1/CS
Retail Packaged	✓	✓	✓	--	✓	✓	--
Retail Box Dimensions	21.5" x 6" x 20.25"	22" x 7" x 19.5"	22" x 7" x 21"	--	13.5" x 8.25" x 22.5"	29" x 10.25" x 30.75"	--
HCPC	E0143 + E0156	E0143 + E0156	E0143 + E0156	E0143 + E0156	E0143 + E0156	E0143 + E0156	E0143 + E0156
Warranty	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime

	ACCESSORIES		
	781-020	A01-831	700-864
Item Number(s)	781-020	A01-831	700-864
Product Name	Hugo Cup Holder	Hugo Seatbelt with End Clips	Seatbelt Accessory for TC's
Key Differences	Universal cup holder fits your favorite beverage container. 3" diameter. Tool-free installation. Securely fits tubing diameters from 5/8" to 1".	Seatbelt with end-clips	Safety Belt for Transport Chair
Available Colors	Black	Black	Black
Case Pack	1/CS	1/CS	1/CS
Retail Packaged	✓		
Retail Box Dimensions	4" x 4" x 3.25"		

LEGEND

ECOM-
Available to order online

PROVID-
Available in independent DME provider locations

RE-
Available in large retail chain stores

GOOD ★
Includes your most basic models and features that meet the customer minimum requirements and provide a good value.

BETTER ★★
Includes similar design and features from our Good models, our Better models offer supplementary features that give the consumer a better user experience.

BEST ★★★
Includes similar design and features from our Good and Better models, our Best models offers unique high-end features providing consumers with a top level experience.