



Wheelchairs | Good, Better, Best Options



Product Name	★ GOOD ★				★★ BETTER ★★		★★★ BEST ★★★					
	Blue Streak	Silver Sport I	Silver Sport II	Silver Sport Recliner	Cruiser III	Cruiser X4	Viper Jr.	Viper Plus GT	Pediatric Viper Plus Recliner	Viper Plus Recliner	Lynx	
Item Numbers	BLS16FBD-SF BLS16FBD-ELR BLS18FBD-SF BLS18FBD-ELR BLS20FBD-SF BLS20FBD-ELR	SSP118FA-SF	SSP216DDA-SF, SSP216DDA-ELR, SSP216DFA-SF, SSP216DFA-ELR, SSP218FA-SF, SSP218FA-ELR, SSP218DDA-SF, SSP218DDA-ELR, SSP218DFA-SF, SSP218DFA-ELR, SSP220DDA-SF, SSP220DDA-ELR, SSP220DFA-SF, SSP220DFA-ELR	SSP16RBDDAV SSP16RBDFAV SSP18RBDDAV SSP18RBDFAV SSP20RBDDAV SSP20RBDFAV	K316DDA-SF, K316DDA-ELR, K316DFA-SF, K316DFA-ELR, K316ADDA-ELR, K316ADFA-ELR, K316ADFA-SF, K316ADDA-SF, K318DDA-SF, K318DDA-ELR, K318DFA-SF, K318DFA-ELR, K318ADDA-SF, K318ADFA-SF, K318ADFA-ELR, K318ADDA-ELR, K320DDA-SF, K320DDA-ELR, K320DFA-SF, K320DFA-ELR, K320ADDA-SF, K320ADDA-ELR, K320ADFA-ELR, K320ADFA-SF	CX416ADDA-ELR, CX416ADDA-SF, CX416ADFA-ELR, CX416ADFA-SF, CX418ADDA-ELR, CX418ADDA-SF, CX418ADFA-ELR, CX418ADFA-SF, CX420ADDA-ELR, CX420ADDA-SF, CX420ADFA-ELR, CX420ADFA-SF	L412DDA-SF L412DDA-ELR L414DDA-SF L414DDA-ELR	PLA416FBUARAD-SF PLA416FBUARAD-ELR PLA418FBUARAD-SF PLA418FBUARAD-ELR PLA420FBUARAD-SF PLA420FBUARAD-ELR PLA422FBUARAD-SF PLA422FBUARAD-ELR	PL412RBDDA PL414RBDDA	PLA416RBDDA PLA418RBDDA PLA418RBDDA PLA418RBDFA PLA420RBDDA PLA420RBDFA	K516FBADDA-SF K516FBADDA-ELR K518FBADDA-SF K518FBADDA-ELR K520FBADDA-SF K520FBADDA-ELR	
Sales Channel												
Key Differences	Blue powder-coated frame with flip-back, padded desk-length arms and nylon upholstery that is durable, lightweight and attractive while easy to clean. Comes with swing-away footrests or elevating leg rests and calf strap.	Powder-coated silver vein finish, this wheelchair features urethane tires, push to lock wheels and nylon upholstery that is lightweight and easy to clean. Along with folding back-canes for transport and storage.	Powder-coated silver vein finish, this wheelchair features 8" front casters that are adjustable in three positions, urethane tires, push to lock brakes and embossed vinyl upholstery that is durable, lightweight and easy to clean.	Powder-coated silver vein finish, this wheelchair has a hydraulic reclining mechanism allowing for infinite adjustment up to 180 degrees. Arms have support extensions for reclining position. 8" front casters and precision sealed wheel bearings to ensure long lasting reliability. Comes standard with anti tipplers and swing-away elevating leg rests.	Powder-coated silver vein finish, this wheelchair has built-in seat rail extensions and extendable upholstery easily adjust seat depth from 16" to 18". Composite, mag style wheels are light weight and maintenance free. Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings.	Weighs 34 lbs. or less when properly equipped in the basic configuration (Excluding front rigging). Aluminum and Carbon steel blended frame with combination powder coated and silver vein finish is attractive, chip resistant and easy-to-maintain. 3 position quick release rear axles mounted on composite mag-style wheels with composite hand rims are lightweight and maintenance free.	Weighs under 32 lbs. (excluding front rigging). Carbon steel black powder-coated frame double coated for an attractive, chip-proof, maintainable finish. Composite mag-style wheels with aluminum hand rims are lightweight and maintenance free. Dual axle provides easy transition of seat height to hemi level.	Weighs under 35 lbs. in basic configuration (excluding front rigging). All-Aluminum frame, Dual Axle allows for multiple seat-to-floor height positions. Adjustable length and height armrests allow for conversion between desk, half, and full arm. Adjustable angle back easily adjusts from 5° - 20°. Adjustable height back adjusts from 17" - 19" in 1" increments, adaptable for individual patient needs. Adjustable angle caster forks have three height adjustments and angle adjustability. Tool-free adjustable swing-away footrests also feature easy-to-access release handles. Urethane rear tires and casters offer superior performance, smoother ride and are lightweight. Adjustable anti-tipper with flip-up wheels Never needs to be removed.	Double-coated chrome carbon steel frame for maximum durability. State-of-the-art hydraulic reclining mechanism. Backrest reclines from 90° to 180°. Padded flip-back detachable desk armrests. Composite lightweight 24" mag-style rear wheels with urethane tires and 6" front rubber tires. Offset rear wheel position provides additional stability and evenly distributes users weight in all reclining positions.	Carbon steel frame with double-coated chrome for an attractive, chip-proof, maintainable finish. Hydraulic reclining mechanism allows for adjustments up to 180°. Spreader bar acts as a back stabilizer and a convenient push handle for the attendant. Offset rear wheel position adds additional stability and equally distributes the patients weight while in the various reclining positions. Built-in seat rail extension and extendable upholstery easily adjusts seat depth from 16" or 18".	Ultra-light, weighs less than 29 lbs. Multi-position back can adjust 2". Quick-release rear wheels, semi-pneumatic tires on a 24" spoke wheel provide a smooth, true roll. 2" of seat depth adjustment. Front casters are height adjustable in three positions. Arm rest are flip back and height adjustable 9"-12" from seat.	
Back Height	16"	16"	16"	21" - 31" Reclining	16"	16" - 18"	16" - 18"	17" - 19"	26" Reclining	23" - 33" Reclining	17" - 19"	
Seat Width	16", 18", 20"	18"	16", 18", 20"	16", 18", 20"	16", 18", 20"	16", 18", 20"	12", 14"	16", 18", 20", 22"	12", 14"	16", 18", 20"	16", 18", 20"	
Seat Depth	16"	16"	16"	16"	16" - 18"	16"	12", 14"	16" - 18"	12", 14"	16", 18"	16" - 18"	
Seat To Floor Height	19.5"	19.5"	17.5" - 19.5"	20.5"	17.5" - 19.5"	17", 18", 19"	13.5" - 15.5" or 17.5" - 19.5"	17.5" - 19.5"	15" - 16"	17.5" - 19.5"	17.5" - 19.5"	
Arms	Arm Style	Flip Back Padded	Fixed Padded	Detachable Padded Fixed Padded	Reclining Back Padded	Flip Back Padded	Flip back Padded	Flip Back, Slope Padded	Flip Back, Removable	Flip Back Padded	Flip Back Padded	Flip Back, Height Adjustable, Detachable
	Non removeable	✓	✓	-- ✓	--	--	--	--	--	--	--	--
	Full	--	✓	✓	✓	--	--	✓	--	--	✓	--
	Desk	✓	--	✓	✓	--	--	✓	--	✓	✓	✓
	Full Adjustable	--	--	--	--	✓	✓	--	✓	--	--	--
Desk Adjustable	--	--	--	--	✓	✓	--	✓	--	--	✓	
Front Riggings	Plastic	Plastic	Plastic	Metal	Composite	Composite	Composite	Composite	Composite	Composite	Composite	Aluminum
Axle Type	Single	Single	Dual	Single	Dual	Triple	Dual	Dual	Dual	Dual	Dual	Dual
Swing-away Footrest	✓	✓	✓	--	✓	✓	✓	✓	✓	--	--	✓
Swing-away Elevating Leg Rest	✓	--	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Seat Rail Extensions	--	--	--	--	✓	✓	--	✓	--	✓	✓	✓
Quick Release Axle Wheels	--	--	--	--	✓	✓	--	✓	--	--	--	✓
Built in Heel Loops	--	✓	✓	--	✓	✓	✓	✓	✓	--	--	✓
Anti Tipplers	--	--	--	✓	--	--	✓	✓	✓	✓	✓	✓
Product Weight	41 lbs	37 lbs	40 lbs / 41 lbs / 42 lbs	59 lbs / 60 lbs / 61 lbs	35 lbs / 36 lbs	32 lbs / 34 lbs	30 lbs / 31 lbs	34.5 lbs / 35.25 lbs / 36 lbs	32 lbs	49 lbs	28 lbs / 28.5 lbs / 29 lbs	
Weight Capacity	250 lbs	250 lbs	250 lbs / 300 lbs / 350 lbs	300 lbs	300 lbs / 350 lbs	300 lbs	250 lbs	300 lbs / 300 lbs / 300 lbs / 350 lbs	250 lbs	300 lbs	275 lbs	
Color	Blue	Silver Vein	Silver Vein	Silver Vein	Silver Vein	Silver Vein	Black	Deep Green	Chrome	Chrome	Black	
Frame style	Steel	Steel	Steel	Steel	Steel	Steel/ Aluminum	Steel	Aluminum	Steel	Steel	Aluminum	
HCPC	K0001-1	K0001-1, K0001	K0001 / K0002-1, K0001 OR K0002	K0001 / E1226, K0001+E1226	K0003, K0003 + E0990 OR K0195	K0004-1, K0004	K0004, E1236+K0195 OR E0990	K0004, K0004+E0973+K0195 OR E0990+E0971(ELR)+E0951 (SF)+E2201(20"-22")	E1236 E1226-1, E1236, E1226	K0004, E1226-1, K0004 + E1226 + K0195 OR E0990	K0005-1, K0005+E0951+K0195	
Warranty	3 Year	3 Year	3 Year	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	



	★ HEAVY DUTY ★				★★ SPECIALTY ★★		
Product Name	Sentra EC HD	Sentra HD	Sentra EC HD Extra Wide	Sentra Recliner	Rebel Wheelchair	Poly-Fly	
Item Numbers	STD20ECDDAHD-SF, STD20ECDDAHD-ELR, STD20ECDFAHF-SF, STD20ECDFAHF-ELR, STD22ECDDA-SF, STD22ECDDA-ELR, STD22ECDFAS-SF, STD22ECDFAS-ELR, STD24ECDDA-SF, STD24ECDDA-ELR, STD24ECDFAS-SF, STD24ECDFAS-ELR	STD20DDA-SF, STD20DDA-ELR, STD22DDA-SF, STD22DDA-ELR, STD24DDA-SF, STD24DDA-ELR,	STD26ECDDA, STD26ECDFAS, STD28ECDDA, STD28ECDFAS, STD30ECDDA, STD30ECDFAS	STD22RBDDA		RTLREB18DDA-SF	TR18 TR20
Images							
Sales Channel							
Key Differences	Weight Capacity: 450 lbs. Carbon Steel frame with silver vein finish. Double embossed vinyl upholstery is durable, lightweight, attractive and easy-to-clean. Standard dual axle provides easy transition of seat to hemi-level. Urethane tires mounted on composite wheels provide durability and low maintenance. 8" front casters, adjustable in three positions.	Weight Capacity: 500 lbs. Carbon steel frame with silver vein finish. Reinforced steel gussets at all weight bearing points provide additional strength. Reinforced side frame, caster journals and dual cross bars. Dual-position front fork. Dual-position front frame.	Weight Capacity: 700 lbs. Carbon steel frame with, chip-proof, maintainable finish. Reinforced steel gussets at all weight bearing points provide additional strength. Durable reinforced nylon upholstery. Steel, spoke wheels are strong and maintenance free. Heavy-duty, flat-free 8" front casters are adjustable in two positions. Extra-heavy-duty front forks.	Weight Capacity: 450 lbs. This wheelchair has a hydraulic reclining mechanism allowing for infinite adjustment up to 180 degrees. Arms have support extension for reclining position. 8" front casters and precision sealed wheel bearings to ensure long lasting reliability. Comes standard with anti tipplers and swing-away elevating leg rests. Front caster forks are adjustable in 2 positions.	Packaged in a retail box, this lightweight wheelchair is easy to navigate in tight areas. It has a fold down back and quick release wheels for easy storage and transport.	When it comes to unique and different, Ploy-Fly is hard to beat. It can be effortlessly transformed from wheelchair to transport chair. The arm are easily removable for easy patient transfer.	
Back Height	16", 18"	18"	18"	23" - 33" Reclining	17"	16"	
Seat Width	20", 22", 24"	20", 22", 24"	26", 28", 30"	22"	18"	18", 20"	
Seat Depth	18"	18"	20"	18"	16"	16"	
Seat To Floor Height	17.5" - 19.5"	17.5" - 19.5"	19.5"	21"	19.5"	19.5"	
Arms	Arm Style	Detachable Padded	Detachable Padded	Detachable Padded	Fixed Padded	Detachable Padded	
	Non removable	--	--	--	✓	--	
	Full	✓	--	✓	--	--	
	Desk	✓	✓	✓	✓	✓	
	Full Adjustable	--	--	--	--	--	
Desk Adjustable	--	--	--	--	--	--	
Front Riggings	Aluminum	Aluminum	Aluminum	Aluminum	Plastic	Composite	
Axle Type	Dual	Dual	Single	Single	Single	Single	
Swing-away Footrest	✓	✓	--	--	✓	✓	
Swing-away Elevating Leg Rest	✓	✓	--	--	--	--	
Seat Rail Extensions	--	--	--	--	--	--	
Quick Release Axle Wheels	--	--	--	--	✓	--	
Built in Heel Loops	--	--	--	--	--	✓	
Anti Tipplers	--	--	--	--	--	--	
Product Weight	56 lbs / 57 lbs / 58 lbs	54 lbs / 56 lbs / 66 lbs	62 lbs / 64 lbs / 66 lbs	62 lbs	37 lbs	34 lbs / 36 lbs	
Weight Capacity	450 lbs	500 lbs	700 lbs	450 lbs	300 lbs	250 lbs	
Color	Silver Vein	Silver Vein	Silver	Silver	Flame Red	Silver	
Frame style	Steel	Steel	Steel	Steel	Steel	Aluminum	
HCPC	K0007, K0007+K0195 OR E0990	K0007, K0007+K0195 OR E0990	K0007-1, K0007	K0001 / E1226-1, K0001 OR K0002+E1226+E0990 OR K0195	N/A	E1038 or K0004+E0951+E0971+E0978	
Warranty	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	

LEGEND

ECOMMERCE
Available to order online

RETAIL
Available in large retail chain stores

PROVIDER
Available in independent DME provider locations

GOOD ★
Good models are the base within the class, comparable to the competition of the same level. These meet the basic requirements with less refinements and minimal extras vs. the "BETTER" and "BEST" classifications. They tend to be entry level in the market with good value.

BETTER ★★
Mid-level models within the class, offer more features and benefits to give the end user a better user experience. These tend to be lighter weight than the good models and have more options available. They may have a higher user weight limit as well.

BEST ★★★
Best models combine similar design features from our "Good" and "Better" models but offer additional features that pushes them to the top ranking. They may be very light or could come with more options such as anti tipplers or flip back arms included. They tend to be made from lighter materials without compromising strength to give the end user the experience they are looking for.