

NEW! **Rose Comfort Max**
Tilt and Recline Chair with Casters



- ✓ **ANGLE ADJUSTABLE**
OFFERS UP TO 38° OF TILT-IN-SPACE AND 90° OF RECLINE
- ✓ **PADDED HEADREST**
PROVIDES CUSTOMIZABLE HEAD AND NECK SUPPORT
- ✓ **BUILT-IN VERSATILITY**
ARMRESTS AND FOOTREST ADJUST TO INDIVIDUAL HEIGHT, LENGTH AND ANGLE PREFERENCES
- ✓ **ADJUSTABLE SEAT**
SEAT-TO-FLOOR HEIGHT STARTING AT 17"
- ✓ **LOCKING CASTERS**
REAR-DIRECTIONAL FOR EASY MANEUVERABILITY
- ✓ **COLOR-CODED LEVERS**
RED LEVER ON RIGHT SIDE CONTROLS RECLINE. BLACK LEVER ON LEFT SIDE CONTROLS TILT.



The new Drive DeVilbiss Long Term Care Tilt and Recline Chair with Casters offers exceptional positioning and comfort over traditional wheeled geriatric chairs. The chair design helps minimize carton size to optimize shelf and storage space.

ITEM #	PRODUCT DESCRIPTION	UOM
CLN400-20	Tilt and Recline Chair with Casters	1/ea
CLN400-HPP	Huntington's Padding Package includes Long Shoulder Bolsters, Positioning Thigh Belt, Padded Footbox	1/ea

PRODUCT SPECS

Overall Width	25.8"
Arm Pad Length	14.0"
Arm Pad Width	2.87"
Seat Depth	17.7"
Back Height	32.6"
Width Between Two Shoulder Rest	13.0"
Shoulder Pillow Length	17.7"
Leg Pad Height	12.4"
Leg Pad Length	19.3"
Foot Pad Width	12.2"
Width Between Shoulder Tube	22.4"
Weight Capacity	400 lb.
Limited Warranty*	Lifetime (10 Year) on Frame

ENHANCED POSITIONING



Huntington Pads provide comfortable positioning, increased sitting tolerance, and weight redistribution to help prevent pressure injuries.



Trendelenburg capable. Lays flat for added comfort and fewer transfers in and out of a bed or standard chair.



Multiple angle adjustments create a supportive and customizable seating option for varying needs.



Rear-directional locking casters provide easy maneuverability when in use.

* Refer to the manual for more information.