



# Patient Lifts | Good, Better, Best Options



Item Number	★ GOOD ★		★★ BETTER ★★		★★★ BEST ★★★				★★★ BEST ★★★	
	13023	13023SV	13240	13244	FLNP500	FLNP600	FLP500	FLP600	STSM450	STSP450
Lift Type	Hydraulic Floor		Powered Floor		Powered Floor				Powered Sit-To-Stand	
Base	Manual Base	Manual Base	Manual Base	Manual Base	Manual Base	Manual Base	Powered Base	Powered Base	Manual Base	Power Base
Weight Capacity	450 lbs	450 lbs	450 lbs	600 lbs	500 lbs	600 lbs	500 lbs	600 lbs	450 lbs	450 lbs
Material	Steel	Steel	Steel	Steel	Steel	Steel	Steel	Steel	Steel	Steel
Foldable Design	--	--	--	--	--	--	--	--	--	--
Floor Clearance	5.25"	5.25"	4.75"	4.75"	5.25"	5.25"	5.25"	5.25"	4.45"	4.45"
Weight	85 lbs	85 lbs	108 lbs	108 lbs	115 lbs	115 lbs	115 lbs	115 lbs	148 lbs	148 lbs
Battery Quantity	0	0	1	1	1	1	1	1	1	1
Charger	--	--	Built-in	Built-in	Built-in & charger stand	Built-in & charger stand	Built-in & charger stand	Built-in & charger stand	Built-in & charger stand	Built-in & charger stand
Visual Battery Fuel	--	--	--	--	✓	✓	✓	✓	✓	✓
Audio Battery Fuel	--	--	✓	✓	✓	✓	✓	✓	✓	✓
Base Length	40.5"	40.5"	45.5"	45.5"	51"	51"	51"	51"	42"	42"
Base Open (O)	22"	22"	24"	24"	43"	43"	43"	43"	40.5"	40.5"
Base Close (O)	42.25"	42.25"	41"	41"	29"	29"	29"	29"	25"	25"
Boom Minimum	17"	17"	27.5"	27.5"	26"	26"	26"	26"	32"	32"
Boom Maximum	61"	61"	76"	76"	76"	76"	76"	76"	62"	62"
Dual Controls	--	--	--	--	✓	✓	✓	✓	✓	✓
Scale	Option - 13046-HD	Option - 13046-HD	Option - 13046-HD	Option - 13046-HD	Option - FLSCALE	Option - FLSCALE	Option - FLSCALE	Option - FLSCALE	--	--
SMART	--	--	--	--	✓	✓	✓	✓	✓	✓
Overload Protection	--	--	--	--	✓	✓	✓	✓	✓	✓
Anti-Entrapment	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Maintenance Track	--	--	--	--	✓	✓	✓	✓	✓	✓
Lift Frame	Lifetime Limited	Lifetime Limited	Lifetime Limited	Lifetime Limited	Lifetime Limited	Lifetime Limited	Lifetime Limited	Lifetime Limited	Lifetime Limited	Lifetime Limited
Hyd Pump and Lift	1 year	1 year	--	--	--	--	--	--	--	--
Electronics	--	--	1 year	1 year	5 years	5 years	5 years	5 years	5 years	5 years
Battery	--	--	90 days	90 days	1 year	1 year	1 year	1 year	1 yr	1 year

### LEGEND



ECOMMERCE  
Available to order online



RETAIL  
Available in large retail chain stores



PROVIDER  
Available in independent DME provider locations

#### GOOD ★

Includes your most basic models and features that meet the customer minimum requirements and provide a good value.

#### BETTER ★★

Includes similar design and features from our Good models, our Better models offer supplementary features that give the consumer a better user experience.

#### BEST ★★★

Includes similar design and features from our Good and Better models, our Best models offers unique high-end features providing consumers with a top level experience.



		★ GOOD ★		★★ BETTER ★★		★★★ BEST ★★★				★★★ BEST ★★★	
Item Number		13023	13023SV	13240	13244	FLNP500	FLNP600	FLP500	FLP600	STSM450	STSP450
Lift Type		Hydraulic Floor		Powered Floor		Powered Floor				Powered Sit-To-Stand	
	Full Body Sling w/ Commode	13221L	✓	✓	✓	✓	✓	✓	✓	--	--
		13221M	✓	✓	✓	✓	✓	✓	✓	--	--
		13221XL	✓	✓	✓	✓	✓	✓	✓	--	--
		13221XXL	✓	✓	✓	✓	✓	✓	✓	--	--
	Full Body Sling Mesh	13223M	✓	✓	✓	✓	✓	✓	✓	--	--
		13223L	✓	✓	✓	✓	✓	✓	✓	--	--
	Full Body Sling Solid	13222M	✓	✓	✓	✓	✓	✓	✓	--	--
		13222L	✓	✓	✓	✓	✓	✓	✓	--	--
		13224XL	✓	✓	✓	✓	✓	✓	✓	--	--
	Divided Leg Sling w/ Headrest	13262L	✓	✓	✓	✓	✓	✓	✓	--	--
		13262M	✓	✓	✓	✓	✓	✓	✓	--	--
	Seat/Back Sling Canvas	13012	✓	✓	✓	✓	✓	✓	✓	--	--
		13025	✓	✓	✓	✓	✓	✓	✓	--	--
		13026	✓	✓	✓	✓	✓	✓	✓	--	--
		13060	✓	✓	✓	✓	✓	✓	✓	--	--
		13061	✓	✓	✓	✓	✓	✓	✓	--	--
	Straps for Patient Slings	13232	✓	✓	✓	✓	✓	✓	✓	--	--
	U-Sling w/ Head Support	13220S	✓	✓	✓	✓	✓	✓	✓	--	--
		13220M	✓	✓	✓	✓	✓	✓	✓	--	--
		13220L	✓	✓	✓	✓	✓	✓	✓	--	--
		13220XL	✓	✓	✓	✓	✓	✓	✓	--	--
	Stand-Up Sling	13263A	--	--	--	--	--	--	--	✓	✓
		13263B	--	--	--	--	--	--	--	✓	✓
	Sit-to-Stand Sling	13263E	--	--	--	--	--	--	--	✓	✓
		13263F	--	--	--	--	--	--	--	✓	✓
	Sit-to-Stand-Sling-Strap	13263H	--	--	--	--	--	--	--	✓	✓
	Transport Sling	13264A	--	--	--	--	--	--	--	✓	✓