



DURABLE
WELDED FRAME

ADJUSTABLE
HANDLE LENGTH

5 ADJUSTMENT
OPTIONS

WHEELED BASE
FOR STORAGE

WIDE, STURDY
BASE

The new free-standing bariatric trapeze, P750TPZ, is a floor mounted device and optional accessory to the bariatric long term care bed, P750. This heavy-duty trapeze is intended to support a portion of an individual's weight to aid in the repositioning of their body while in bed or assist in the transfer from bed to a chair.



For use with P750 Bed



5 adjustment
options



Universal base fits
under most beds

ITEM

UPC

UOM

P750TPZ

822383015835

1/cs

FEATURES

- Allows for independent resident positioning
- Assists with transfers into and out of the bed
- 600 lb. weight capacity
- Durable welded frame construction
- Adjustable handle with strap available with 5 adjustment options (2" increments)
- Adjustable sling mounting position on boom
- Warranty: 2 years, 10-year weld

PRODUCT SPECS

Total Height	74.25"
Total Width	36.5"
Total Depth	59.25"
Width Between Legs	29.25"
Base Clearance Height	1.5"
Boom Length	45"
Product Weight	127 lbs.
Shipping Weight	138 lbs.

*Refer to instruction guide or website for warranty details.