



LONG TERM  
*Care*

# RESIDENT FLOOR LIFT SCALE

## STANDARD FEATURES

- Recommended for use on DRIVE LEVANTAR (FLNP500/FLP500) or GRAVIS (FLNP600/FLP600) floor lifts
- Style and color compliment LEVANTAR and GRAVIS floor lifts
- 660 lbs/300 kg max weight capacity
- Easy to install
- Will not fit on old style 13240, 13242, 13244 and 13245 DRIVE Floor lifts

## SPECIFICATIONS

| Item #             | FLSCALE   |
|--------------------|---|
| Weight Range       | 660 lbs. / 300 kg   |
| Display Resolution | 0.2 lb / 0.1 kg   |
| Accuracy           | 0.2 lb / 0.1 kg   |
| Display Type       | Numerical LCD<br>0.7 mm / 17.8 mm                         |
| Size               | 3.2 in x 1.8 in x 5.5 in<br>(81.3mm x 45.7 mm x 139.7 mm) |
| Unit Weight:       | 1.37 lbs / 620 g  |
| Power Supply       | 3-AA batteries  |
| Battery Life       | Approximately 150 hours                                   |

To Learn More/Request a Demo: Call our LTC Customer Service Hotline at 1-877-374-8763, Contact your Sales Representative, or Visit [DriveMedical.com](http://DriveMedical.com) TODAY!

**drive**

**DeVilbiss**  
HEALTHCARE