

Blood Pressure Monitor

Upper Arm Model



Brand		Drive DeVilbiss	Omron 5-series
Model		BP3400	BP742N
Type		Upper-Arm Type	Upper-Arm Type
Measurement Method		Oscillometric Method	Oscillometric Method
Measurement Range	Pressure	0 to 300mmHg	0 to 299mmHg
	Pulse	30 to 180 Beat/Minute	40 to 180 Beat/Minute
Measurement Range	Pressure	±3mmHg	±3mmHg
	Pulse	±5%	±5%
Pressurization		Automatic Pressurization	Automatic Pressurization
LCD Size		2.17" x 3.31" / 55.1 x 84.1mm	N/A
Memory Bank (Users)		2x60	2x50
Data/Time		Month+Day+Hour+Minute	Year+Month+Day+Hour+Minute
Irregular Heartbeat Detection		Yes	Yes
WHO Classification Indicator		Yes	No; has indicator but not WHO
Averaging Function		Yes	Yes
Low Battery Detection		Yes	Yes
Automatic Power-off		Yes	Yes
Power Source		4 Alkaline Batteries Size AA included; AC Adapter sold separately	4 Alkaline Batteries Size AA Not Included; AC Adapter sold separately
Battery life		About 2 months (test 3 minutes per day, 30 days/per month)	Approximately 1000 uses
Unit Dimensions		5.83" x 3.94" x 2.2" / 14.8 x 10.0 x 5.6 cm	4.25" x 3.125 x 5.5" / 10.8 x 8 x 14 cm
Included Cuff Size		8.66" to 14.17" / 22 to 36 cm	9" to 17" / 22 to 42 cm
Storage Case		Yes	No



Visit Us at www.DriveMedical.com



Call Us Today! 877.224.0946



Contact Your Sales Representative

Blood Pressure Monitor

Upper Arm Model



Brand		Drive DeVilbiss	Omron 3-series
Model		BP3600	BP710N
Type		Upper-Arm Type	Upper-Arm Type
Measurement Method		Oscillometric Method	Oscillometric Method
Measurement Range	Pressure	0 to 300mmHg	0 to 299mmHg
	Pulse	30 to 180 Beat/Minute	40 to 180 Beat/Minute
Measurement Range	Pressure	±3mmHg	±3mmHg
	Pulse	±5%	±5%
Pressurization		Automatic Pressurization	Automatic Pressurization
LCD Size		1.81" x 2.45" / 46.0 x 62.3 mm	N/A
Memory Bank (Users)		120	Up to 14
Data/Time		Month+Day+Hour+Minute	Year+Month+Day+Hour+Minute
Irregular Heartbeat Detection		Yes	Yes
WHO Classification Indicator		Yes	No; has indicator but not WHO
Averaging Function		Yes	Yes
Low Battery Detection		Yes	Yes
Automatic Power-off		Yes	Yes
Power Source		4 Alkaline Batteries Size AAA included; AC Adapter sold separately	4 Alkaline Batteries Size AA Not Included; AC Adapter sold separately
Battery life		About 2 months (test 3 minutes per day, 30 days/per month)	Approximately 1000 uses
Unit Dimensions		5.51" x 3.86" x 1.89" / 14.0 x 9.8 x 4.8 cm	4" x 3.125" x 5.125" / 10.2 x 7.9 x 13cm
Included Cuff Size		8.66" to 14.17" / 22 to 36 cm	9" to 17" / 22 to 42 cm
Storage Case		Yes	No



Visit Us at www.DriveMedical.com



Call Us Today! 877.224.0946



Contact Your Sales Representative

Blood Pressure Monitor

Wrist Model



Brand		Drive DeVilbiss	Omron
Model		BP3200	BP629N
Type		Wrist Type	Wrist Type
Measurement Method		Oscillometric Method	Oscillometric Method
Measurement Range	Pressure	0 to 300mmHg	0 to 299mmHg
	Pulse	30 to 180 Beat/Minute	40 to 180 Beat/Minute
Measurement Range	Pressure	±3mmHg	±3mmHg
	Pulse	±5%	±5%
Pressurization		Automatic Pressurization	Automatic Pressurization
LCD Size		1.93" x 1.50" / 49.0 x 38.0 mm	N/A
Memory Bank (Users)		2x60	60
Data/Time		Year+Month+Day+Hour+Minute	Year+Month+Day+Hour+Minute
Irregular Heartbeat Detection		Yes	Yes
WHO Classification Indicator		Yes	No; has indicator but not WHO
Averaging Function		Yes	Yes
Low Battery Detection		Yes	Yes
Automatic Power-off		Yes	Yes
Power Source		2 Alkaline Batteries Size AAA included	2 Alkaline Batteries Size AAA Not Included
Battery life		About 2 months (test 3 minutes per day, 30 days/per month)	Approximately 300 uses
Unit Dimensions		3.31" x 2.52" x 1.14" / 8.4 x 6.4 x 2.9 cm	3" x 2.375" x 7.875" / 7.6 x 6 x 20 cm
Included Cuff Size		5.31" to 8.46" / 13.5 to 21.5 cm	5.25" to 8.5" / 13.5 to 21.5 cm
Storage Case		Yes	Yes



Visit Us at www.DriveMedical.com



Call Us Today! 877.224.0946



Contact Your Sales Representative