



99 Seaview Blvd.,
Port Washington, NY 11050



www.drivemedical.com
Customer Service | 877.224.0946

ELECTRIC BREAST PUMP

COMPARISON CHART



Model	Pure Expressions RTLBP2000	Pure Expressions RTLBP2000	Pure Expressions RTLBP1000
Item Description	Double Electric Breast Pump	Economy Double Electric Breast Pump	Single Electric Breast Pump
Adjustable Vacuum	Yes	Yes	Yes
BPA Free	Yes	Yes	Yes
Latex Free	Yes	Yes	Yes
Double Electric Pumping	Yes	Yes	No
Single Electric Pumping	Yes	Yes	Yes
Closed System	Yes	No	Yes
Expression Phases	2	2	2
Vacuum Chamber	2	1	1
Vacuum Strength	30 mm Hg - 270 mm Hg	30 mm Hg - 285 mm Hg	30 mm Hg - 270 mm Hg
Suction Level	9	5	9
Independent Suction Strength for each breast	YES	No	No
Clock/ Timer	Yes (30 min countdown)	Yes (30 min countdown)	Yes (30min countdown)
Memory Function	Yes (settings last used)	Yes (settings last used)	Yes (settings last used)
LCD Display	Yes	Yes	Yes
Portability	Yes	Yes	Yes
Pump Dimensions	4" x 6.5" x 8"	3.5" x 5" x 7.2"	2.2" x 4.2" x 5"
Pump Weight	0.82 lbs (13.1oz)	0.92 lbs (14.74 oz)	0.54 lbs (8.64 oz)
Weight with box	3 lbs 4oz	2.70 lbs	1.92 lbs
AC Adapter	Yes	Yes	Yes
Battery Option	External battery pack (4 AA batteries)	External battery pack (6 AA batteries)	Internal Battery Pack (4 AA batteries)
Breast Shield Sizes*	(1) 27mm, 24mm	(1) 24mm	(1) 27mm, 24mm
Breastmilk Containers	2 (4 oz)	2 (4 oz)	1 (4 oz)
Carrying Bag	Yes- Carry bag w/zipper	No	Yes- drawstring bag
Insulated Cooler Bag*	No	No	No
Tubing Length	36	36	36

*Sold Separately