



**135**  
YEARS  
EST. 1888

## Your Guide to Nebulizers

---

Enhancing quality  
of life is our  
driving force

**drive**

---

 **DeVilbiss**  
HEALTHCARE

## Taking Respiratory Therapy to New Places Expand Your Care & Achieve Optimal Patient Outcomes

Drive DeVilbiss Healthcare relies on US based manufacturing, research and development, and quality control, as well as over 100 years of combined experience on Team Respiratory to bring products to our Provider Partners which will exceed the market demands for clinical efficacy, reliability, and price.

Aerosol peak season typically begins in October and runs through May. The season follows the onset of temperature changes throughout much of the country, along with the arrival of the Influenza (Flu), Respiratory Syncytial Virus (RSV) and Coronavirus (COVID-19) seasons, which run concurrently and contribute to exacerbations, pneumonia, bronchitis and bronchiolitis.

### How do nebulizers work?

- Compressed air is used to break up a liquid medication into small aerosol droplets (mist) that can be inhaled through a mouthpiece or mask
- Inhaled mist only penetrates into the lungs if the nebulizer creates mist droplets that have a small diameter (1-5 microns). Otherwise, the mist droplets get caught in the mouth and trachea and the medication does not get to where it needs to go
- Compressor nebulizers can utilize disposable and reusable neb kits. Reusable neb kits have been proven to provide smaller and more uniform particle size during nebulization. Reusable neb kits are not only superior clinically but also a wise business model

### Disposable vs. Reusable Nebs

**HCPCS Code:** A7003 - Disposable Neb Kit

**HCPCS Code:** A7005 - Reusable Neb Kit

# Standard Nebulizers

## Pacifica Elite

Item # 18070 With Disposable Neb Kit

Item # 18071 With Disposable and Reusable Neb Kit

- This retail-packaged compressor neb is a workhorse with higher liter flows, shorter treatment times and quieter operation than similar economy compressor nebs.
- Ultra compact design makes it extremely portable
- Powerful piston pump is efficient, dependable and quiet
- Delivers consistent particle size of <5 microns and fast treatment times

[Learn More](#)



## Power Neb Ultra

Item # 18080 With Disposable Neb Kit

Item # 18081 With Disposable and Reusable Neb Kit

- A competitive edge price point for providers as one of the smallest compressor nebs at a mere 3.2 lbs
- Built-in neb kit holder for hands-free use
- Powerful piston pump is efficient, dependable and quiet
- Delivers consistent particle size of <5 microns and fast treatment times

[Learn More](#)



# Standard Nebulizers

## PulmoNeb® LT Compressor Nebulizer

Item # 3655LT With Disposable Neb Kit

Item # 3655LTR With Disposable and Reusable Neb Kit

- This affordable yet advanced system was developed to exceed market standards and provide superior clinical results and performance
- Small footprint and lightweight design make it easy-to-store and transport
- Ball bearing connecting rod design equips the unit for years of dependable performance
- Top-of unit connector facilitates easy attachment/removal of tubing

[Learn More](#)



## Compartment Style Compressor Nebulizer

Item # MQ6000 With Disposable Neb Kit

- A heavy duty compressor neb with a storage compartment for neb accessories, perfect for home and for long term care facilities.
- Includes disposable neb kit and tubing
- Attractive, sleek design
- Durable, quiet system
- Piston-driven compressor
- Convenient storage compartment

[Learn More](#)



# Standard Nebulizers

## MEDNEB+ Compressor Nebulizer

Item # MQ5600 With Disposable Neb Kit

Item # MQ5700 With Disposable with Reusable Neb Kit

Item # MQ5700B With Disposable and Reusable Neb Kit and Carry Bag

- Includes Airial™ compressor, disposable and optionable reusable nebulizer kit, 7' kink resistant tubing, mouthpiece, and extra filters
- Compact, desk-top unit
- Piston-driven system
- Available with a carry bag, to maximize portability

[Learn More](#)



## Compact Compressor Nebulizer

Item # MQ5800 With JetNeb Reusable Neb Kit

Item # MQ5900 With Disposable and JetNeb Reusable Neb Kit

Item # MQ5900P With Disposable and JetNeb Reusable Neb Kit and Pediatric Mask

- Our premium compressor nebulizer with integrated carry handle comes standard with a breath enhanced JetNeb Plus reusable neb kit to optimize treatments.
- Durable, lightweight, quiet compressor is perfect for a wide range of patients

[Learn More](#)



# Standard Nebulizers

## PulmoMate® Compressor Nebulizer

Item # 4650D With Disposable Neb Kit

- The sleek case design of the PulmoMate® Compressor Nebulizer System fits easily into luggage, backpacks and diaper bags for easy transport. The DeVilbiss PulmoMate features an updated compressor/motor for long-term durability and performance.
- Consistent performance over the life of the unit.
- Metal frame motor construction.
- Integrated handle for easy carrying.
- Large built-in compartment for storing nebulizers and tubing.
- Intuitive design is easy-to-use and maintain.

[Learn More](#)



## Pulmo-Aide® Compressor Nebulizer

Item # 5650D With Disposable Neb Kit

- The Pulmo-Aide Compressor Nebulizer System is the only approved compressor for use with TOBI (Tobramycin).
- One of the longest lasting compressors on the market with an average product life of 10 years.
- Unmatched quality and durability – Great for nursing homes and long-term users.
- Metal frame motor construction.
- Large built-in compartment for storing nebulizers and tubing.

[Learn More](#)



## Pulmo-Aide® Compact Compressor Nebulizer

Item # 3655D With Disposable Neb Kit

- The Pulmo-Aide Compact Compressor Nebulizer System was built to last with high-quality metal mechanical parts, ensuring quality therapy for an extended period of time
- Small footprint
- Small, compact design
- Convenient carrying handle

[Learn More](#)

## DeVilbiss Traveler® Portable Compressor Nebulizer System

Item # 6910P-DR With Battery

Item # 6910D-DR Without Battery

- Delivers a consistent and concentrated treatment of small particles that penetrate deeply into the lungs. Includes disposable and reusable neb kits, pediatric and adult masks, AC and DC adapters, and carry bag.
- Easy, one-touch operation
- Fast treatment time (approximately 10 minutes) with DeVilbiss 900 Series Reusable Nebulizer encourages patient compliance

[Learn More](#)



# Pediatric Character Nebulizers

## Fish Pediatric Compressor Nebulizer

Item # 18090-FS With Disposable Neb Kit

Item # 18091-FS With Disposable and Reusable Neb Kit

- Effective treatments in a child friendly design with a fish bowl style carry bag included for all of your accessories.
- Easy-to-replace filter
- Powerful piston pump is efficient, dependable and quiet
- Complete with disposable neb kit, mouth-piece, pediatric aerosol mask, carry bag, replacement filters and instructions

[Learn More](#)



## Beagle Pediatric Compressor Nebulizer

Item # 18090-BE With Disposable Neb Kit

Item # 18091-BE With Disposable and Reusable Neb Kit

- The dependable and efficient Beagle Nebulizer is designed with pediatric patients in mind, and includes a doggy house carry bag and doggy stickers.
- Complete with disposable neb kit, mouthpiece, pediatric aerosol mask, replacement filter and instructions
- Powerful piston pump is efficient, dependable and quiet

[Learn More](#)





## Penguin Pediatric Compressor Nebulizer

- Item # MQ6002 With Disposable Neb Kit
- Item # MQ6002R With Disposable and Reusable Neb Kit
- Item # MQ6006 With Disposable Neb Kit, No Carry Bag

- Powerful piston driven compressor system in a child friendly design and available with an igloo carry bag.

[Learn More](#)



## Panda Pediatric Compressor Nebulizer

- Item # MQ6003 With Disposable Neb Kit 1/ea
- Item # MQ6003R With Disposable and Reusable Neb Kit
- Item # MQ6004 With Disposable Neb Kit, No Carry Bag
- Item # MQ6005 With Disposable Neb Kit
- Item # MQ6005R With Disposable and Reusable Neb Kit

- Powerful piston pump for fast and effective treatment times, in a child friendly design. Available with a carry bag.

[Learn More](#)



# Holding Chambers

## Airial™ Holding Chambers

- Item # MQ8000 Airial™ Holding Chamber
- Item # MQ8000S Airial™ Holding Chamber with Small Mask
- Item # MQ8000M Airial™ Holding Chamber with Medium Mask
- Item # MQ8000L Airial™ Holding Chamber with Large Mask

- For use with most CFC and HFA metered dose inhalers to momentarily hold aerosolized medication minimizing the need for hand-breath coordination.
- Helps to eliminate waste of medication
- Minimizes need for hand-breath coordination
- Maximizes drug delivery
- Easy-to-use and clean

[Learn More](#)



## Airial Spotz® the Dog Mask with Meter Dose Inhaler Chamber

Item # MQ8100

- Pediatric holding chamber with Spotz® the Dog mask for use with most CFC and HFA metered dose inhalers to momentarily hold aerosolized medication, minimizing the need for hand-breath coordination.
- Helps to eliminate waste of medication
- Minimizes need for hand-breath coordination
- Includes Spotz® pediatric mask for use with chamber
- Easy-to-use and clean

[Learn More](#)



## Universal Neb Carry Bag

Item # 18031      9" (H) x 12" (W) x 7" (D)

[Learn More](#)



## Pediatric Aerosol Pacifier

Item # MQ0385      Pacifier Only

[Learn More](#)



## Puppy Aerosol Mask

Item # DL1050      Puppy Aerosol Mask

[Learn More](#)



# Aerosol Accessories

## Adult & Pediatric Aerosol Mask

Item # MASK 001A

Adult

Item # MASK 001P

Pediatric

[Learn More](#)



MASK 001A

MASK 001P

## Disposable Neb Kit

Item # NEB KIT 500

With Mouthpiece

Item # NEB KIT 600

With Adult Aerosol Mask

Item # NEB KIT 700

With Pediatric Aerosol Mask

[Learn More](#)



NEB KIT 500

NEB KIT 600

NEB KIT 700

## JetNeb Plus Reusable Neb Kit

Item # 18260

Reusable Neb Kit

[Learn More](#)



## Reusable Neb Kit

Item # 18254

Reusable Neb Kit

[Learn More](#)





CORPORATE OFFICE  99 Seaview Blvd., Port Washington, NY 11050  877-374-8763  [www.drivemedical.com](http://www.drivemedical.com)

CANADA OFFICE  375 New Huntington Rd., Unit 1, Vaughan, Ontario L4H 5B5  866-752-1363

©2023 Medical Depot, Inc. dba Drive DeVilbiss Healthcare, 99 Seaview Boulevard, Port Washington, NY 11050. All trademarks used in association with the sale of products of Drive DeVilbiss are trademarks owned by Medical Depot, Inc. All other trademarks, trade names, service marks and logos referenced herein belong to their respective companies. Refer to product instruction guide for warranty details.