



Cushions | Good, Better, Best Options



Product Name	SEAT CUSHIONS						BACK CUSHIONS		COMBO KITS
	GENERAL USE				SKIN PROTECTION	ADJUSTABLE SKIN PROTECTION	GENERAL USE		GENERAL USE COMBO KITS
	★ GOOD ★	★★ BETTER ★★		★★★ BEST ★★★	★ GOOD ★	★★★ BEST ★★★	★ GOOD ★	★★ BETTER ★★	
	Premier One Foam Cushions	Molded Wheelchair Cushions	Gel "E" Wheelchair Cushions	Gel Cushions	Gel-U-Seat Wheelchair Cushions	Balanced-Aire Adjustable Cushions	Premier One Foam Back Cushions	General Use Adjustable Back Cushions	Seat & Back Cushion Kits
Key Differences	Aerated foam helps keep the user cool & dry	Designed for low risk patients seeking comfort and support	Viscous gel bladder is encased in polyurethane foam shell and provides optimal pressure redistribution, support and comfort	Unique dual-chambered gel bladder keeps the pelvis in contact with pressure redistributing cooling gel	Dual chambered gel bladder provides proper immersion and envelopment for pressure redistribution"	Flexible interconnected, independent air cells allow for deep immersion into the cushion	Durable Masonite board prevents hammocking in sling-back wheelchair. Easily removed as per patient preference	5 Adjustment straps featured for multi-purpose back support and customization to individual seating needs	Kit includes one seat and one back cushion each for added comfort & support
Item # / Dimensions (WxLxH)	M8063 - 16" x 16" x 3" M8064 - 18" x 16" x 3" M8065 - 20" x 16" x 3" M8070 - 20" x 18" x 3" M8077 - 22" X 18" x 3"	14880 - 16" x 16" x 2.5" 14887 - 18" x 16" x 2.5" 14908 - 18" x 18" x 2.5" 14881 - 20" x 16" x 2.5" 14909 - 20" x 18" x 2.5"	14891 - 16" x 16" x 3" 14886 - 18" x 16" x 3" 14878 - 18" x 18" x 3" 14894 - 20" x 16" x 3" 14902 - 22" x 18" x 3" 14888 - 18" x 16" x 2"	8040-1 - 16" x 16" x 2" 8040-2 - 18" x 16" x 2" 8040-3 - 20" x 16" x 2" 8040-4 - 20" x 18" x 2" 8040-5 - 18" x 18" x 2" 8040-7 - 22" x 18" x 2" 8040-9 - 24" x 18" x 2"	8046-1 - 16" x 16" x 3" 8046-2 - 18" x 16" x 3" 8046-3 - 20" x 16" x 3" 8046-4 - 20" x 18" x 3" 8046-5 - 22" x 18" x 3" 8046-7 - 24" x 18" x 3"	8047-16 - 16" x 16" x 4" 8047-16-2 - 16" x 16" x 2" 8047-18 - 18" x 16" x 4" 8047-1818 - 18" x 18" x 4" 8046-3 - 20" x 16" x 3" 8047-1818-2 - 18" x 18" x 2" 8047-18-2 - 18" x 16" x 2" 8047-20 - 20" x 16" x 4" 8047-2018 - 20" x 18" x 4" 8047-20-2 - 20" x 16" x 2" 8047-18- 18" x 16" x 4"	8033 - 22" x 2.5" x 17" 14920 - 20" x 2.5" x 17" 14906 - 16" x 2.5" x 17" 14889 - 18" x 2.5" x 17"	14300 - Fits 16"-21" Wide Wheelchairs 14301 - Fits 22"-26" Wide Wheelchairs	14880KIT - Fits 16" Wide Wheelchairs 14887KIT - Fits 18" Wide Wheelchairs 14881KIT - Fits 20" Wide Wheelchairs
Weight Capacity	M8063 - 250 lbs M8064 - 250 lbs M8065 - 250 lbs M8070 - 250 lbs M8077 - 400 lbs	14880 - 250 lbs 14887 - 250 lbs 14908 - 250 lbs 14881 - 250 lbs 14909 - 300 lbs 14891 - 250 lbs	14891 - 250 lbs 14886 - 250 lbs 14878 - 250 lbs 14894 - 250 lbs 14902 - 500 lbs 14888 - 250 lbs	8040-1 - 275 lbs 8040-2 - 275 lbs 8040-3 - 275 lbs 8040-4 - 275 lbs 8040-5 - 275 lbs 8040-7 - 500 lbs 8040-9 - 500 lbs	8046-1 - 275 lbs 8046-2 - 275 lbs 8046-3 - 275 lbs 8046-4 - 275 lbs 8046-5 - 500 lbs 8046-7 - 500 lbs	500 lbs	--	350 lbs	--
Anti-Slip Feature	12" Tie Straps	12" Tie Straps	Non Skid Bottom & Carry Handle	12" Tie Straps	12" Tie Straps	Velcro Straps	Velcro Straps	Adjustable Velcro Straps	12" Tie Straps
Composition	Polyurethane Foam	Molded Polyurethane Foam	Gel / Foam	Gel / Foam	Gel / Foam	Air	Polyurethane Foam	Polyurethane Foam	Polyurethane Foam
Top Fabric	Fluid-Resistant Stretch Top	Fluid-Resistant Stretch Top	Urethane-coated Polyester Top	Fluid-Resistant Stretch Top	Fluid-Resistant Stretch Top	Fluid-Resistant Stretch Polyester Top	Stretch Fabric, Water Resistant	Polyurethane Foam	Fluid-Resistant Stretch Top
Bottom Fabric	Waterproof Vinyl Bottom	Waterproof Vinyl Bottom	Waterproof Vinyl Bottom	Waterproof Vinyl Bottom	Waterproof Vinyl Bottom	Non-Skid Bottom	Vinyl, Water Proof		Waterproof Vinyl Bottom
Anti-microbial Treated	--	--	--	✓	✓	--	--	--	--
Warranty	1 Year Limited	18 months Limited	18 months Limited	18 months Limited	18 months Limited	2 Year Limited	18 Months Limited	18 Months Limited	18 Months Limited
HCPC Code	E2601	E2601	14891 - E2603 14886 - E2603 14878 - E2603 14894 - E2603 14902 - E2604 14888 - E2601	8040-1 - E2601 8040-2 - E2601 8040-3 - E2601 8040-4 - E2601 8040-5 - E2601 8040-7 - E2602 8040-9 - E2602	8046-1 - E2603 8046-2 - E2603 8046-3 - E2603 8046-4 - E2603 8046-5 - E2604 8046-7 - E2604	E2622	8033 - E2612 14920 - E2611 14906 - E2611 14889 - E2611	14300 - E2611 14301 - E2612	E2601 (seat cushion) & E2611 (back cushion)
Unit of Measure	1/EA or 4/CV	1/EA or 6/CV	1/EA or 4/CV	1/EA or 4/CV	1/EA or 4/CV	1/EA or 6/CV	1/EA or 6/CV	1/EA or 6/CV	1/EA