



	★ STANDARD ★			★★ ENHANCED ★★	★★★ PREMIUM ★★★	SPECIALTY		BARIATRIC
Item #	15019	15670	15770	15870	15970	6500-DE-2-FB	BA9600-NP	1500SC-42-FB
Product Name	Therapeutic 5 Zone	Gravity 6	Gravity 7	Gravity 8	Gravity 9	Dynamic Elite	Balanced Aire	Shearcare 1500
Sales Channel								
Key Differences	Single layer 5 zone pressure redistribution	Single layer 5 zone pressure redistribution with heel slope	Designated heel section	Channel cut heel slope	Visco foam top layer	Gel infused visco foam top layer	Hybrid self-adjusting air/foam	Durable bariatric foam
Dimensions (WxLxH)	36" x 80" x 6"	36" x 80" x 6"	36" x 80" x 6"	36" x 80" x 6"	36" x 80" x 6"	36" x 80" x 6"	35" x 80" x 7"	42" x 80" x 6"
Additional Sizes Available (Item#: WxLxH)	--	--	15876: 36" x 76" x 6" 15784: 36" x 84" x 6" 15777: 36" x 80" x 6" (w/raised rails) 15785: 36" x 84" x 6" (w/raised rails) 15877: 36" x 76" x 6" (w/raised rails) Custom sizes available	15886: 36" x 76" x 6" 15884: 36" x 84" x 6" 15888: 36" x 80" x 6" (w/raised rails) Custom sizes available	15996: 36" x 76" x 6" 15984: 36" x 84" x 6" G9-39WX80LX6H: 39" x 80" x 6" 15970-4280: 42" x 80" x 6" 15970-4284: 42" x 84" x 6" 15977: 36" x 80" x 6" (w/raised rails) Custom sizes available	Custom sizes available	BA9600-NP-42: 42" x 80" x 7" BA9600-NP-84: 36" x 84" x 7" Custom sizes available	1500SC-48-FB: 48" x 80" x 6" Custom sizes available
Layers of foam	1	1	3	3	4	4	1	2
Weight Capacity	350 lbs	350 lbs.	350 lbs	400 lbs	450 lbs	400 lbs	500 lbs	650 lbs (48" W-750 lbs)
Available w/ Raised Rail	--	--	✓ - Standard	✓ - Standard	✓ - Standard	✓ - Custom	✓ - Custom	✓ - Custom
Unique Heel section	--	--	✓	✓	✓	✓	--	✓
Heel Slope	--	✓	--	✓	✓	✓	✓	--
Top Cover Material	Durable nylon	Stretch Polyurethane Nylon	Stretch polyurethane nylon	Stretch polyurethane nylon	Stretch polyurethane nylon	Multi-stretch polyurethane	Multi-stretch polyurethane	Multi-stretch polyurethane
Warranty	3 year foam, 1 year cover	3 year foam, 3 year cover	7 year foam, 3 year cover	10 year foam, 3 year cover	15 year foam, 3 year cover	15 year foam, 3 year cover	5 year foam, 3 year cover	7 year foam, 3 year cover
HCPCS Code	E0184	E0184	E0184	E0184	E0184	E0184	E0373 Also available E0277	E0184

LEGEND



STANDARD ★
Entry level systems suitable to meet basic user needs for a good value.

ENHANCED ★★
Meets basic needs plus supplemental features to improve the user's experience.

PREMIUM ★★★
Includes multiple advanced features for the ultimate user experience.