



Pressure Prevention | Standard, Enhanced and Premium Options



Control Unit Type	ANALOG CONTROL UNITS					DIGITAL CONTROL UNITS				
	★ STANDARD ★					★★ ENHANCED ★★			★★★ PREMIUM ★★★	
Item #	14025NS	14026	14530	14027	14360-B	14029, 14029-84	14029DP	14360-P	14200	14200N
Product Name	Med-Aire 5"	Med-Aire Melody	Med-Aire Assure	Med-Aire	Med Aire Edge	Med-Aire Plus	Med-Aire Plus	Med Aire Edge	Harmony Tri-therapy	Harmony Tri-therapy
Image									TO BE DISCONTINUED 	
Sales Channel										
Key Differences	Air overlay 5" height	Economical nylon cover, visual alarm	3" foam base	Quick cell disconnects	Firm foam base & full foam perimeter	Cell-on-cell design for power outage protection	10" Defined Perimeter	Firm foam base & full foam perimeter	True Low Air Loss w/ Pulsation	True Low Air Loss w/ Pulsation
Dimensions (W x L x H)	36" x 80" x 5"	35" x 80" x 8"	35" x 80" x 8"	36" x 80" x 8"	36" x 80" x 8"	14029 - 36" x 80" x 8" 14029-84 - 36" x 84" x 8"	36" x 80" x 8-10"	36" x 80" x 8"	36" x 80" x 10"	36" x 80" x 10"
Weight Capacity	350 lbs	350 lbs	350 lbs	350 lbs	450 lbs	450 lbs	450 lbs	450 lbs	440 lbs	550 lbs
LPM	8	8	8	8	8	8	8	8	1300	1300
Alternating Pressure	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Low Air Loss (LAL)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
True LAL	--	--	--	--	--	--	--	--	✓	✓
Pillow Function	✓	✓	✓	✓	✓	✓	--	✓	--	--
Quick Cell Disconnects	--	--	--	✓	✓	✓	--	✓	✓	✓
Autofirm	--	--	--	--	--	✓	✓	✓	✓	✓
Anti-microbial Treated	--	--	--	--	✓	--	--	✓	✓	✓
Seat Inflate / Fowler	--	--	--	--	--	✓	✓	✓	✓	✓
Low Pressure Alarm	Visual	Visual	Audible & Visual	Audible & Visual	Audible & Visual	Audible & Visual	Audible & Visual	Audible & Visual	Audible & Visual	Audible & Visual
Power Failure Alarm	--	--	--	--	--	Audible & Visual	Audible & Visual	Audible & Visual	Audible & Visual	Audible & Visual
System Failure Alarm	--	--	--	--	--	Audible & Visual	Audible & Visual	Audible & Visual	Audible & Visual	Audible & Visual
Auto Return from Autofirm to Therapy	--	--	--	--	--	30 min	30 min	30 min	20 min	20 min
Auto Return from Static to Therapy	--	--	--	--	--	120 min	120 min	120 min	--	--
Power Outage Protection	--	--	✓	--	✓	✓	✓	✓	✓	✓
Raised Perimeters	--	--	--	--	--	Optional w/DP model	✓	--	--	--
Top Cover Material	Nylon/Polyurethane	Nylon/Polyurethane	Quilted Nylon/Polyurethane	Quilted Nylon/Polyurethane	Stretch Polyester/Polyurethane	Quilted Stretch Polyester/Polyurethane	Polyester/Polyurethane	Quilted Stretch Polyester/Polyurethane	Stretch Silver Ion Polyurethane	Stretch Silver Ion Polyurethane
Air Cell Count	17	20	17	20	20	20	20	20	20	20
Warranty	24 months	12 months	12 months	12 months	18 months	18 months	18 months	18 months	24 months	24 months
HCPCS	E0372	E0277	E0277	E0277	E0277	E0277	E0277	E0277	E0277	E0277

LEGEND



E-COMMERCE
Available to order online



PROVIDER
Available in independent DME provider locations



LONG TERM CARE
Available in long term care facilities

STANDARD ★
Entry level systems suitable to meet basic user needs for a good value.

ENHANCED ★★
Meets basic needs plus supplemental features to improve the user's experience.

PREMIUM ★★★
Includes multiple advanced features for the ultimate user experience.



Pressure Prevention | Standard, Enhanced and Premium Options



Control Unit Type	ANALOG CONTROL UNIT		DIGITAL CONTROL UNITS					
	SPECIALTY		BARIATRIC					
Item #	BA9600-P, BA9600-P-42, BA9600-P-KIT	LS9500N	14030	14048	14054	14060	14200-42	14200-48
Product Name	Balanced Aire	Lateral Rotation	Med-Aire Plus Bariatric	Med-Aire Plus Bariatric	Med-Aire Plus Bariatric	Med-Aire Plus Bariatric	Harmony Tri-therapy	Harmony Tri-therapy
Sales Channel								
Key Differences	Self-Adjusting Hybrid System w/ heel slope, available in 42"W	Specialty lateral rotation system	42" W sleep surface	48" W sleep surface	54" W sleep surface	60" W sleep surface	Feature rich True low air loss 42" W sleep surface	Feature rich True low air loss 48" W sleep surface
Dimensions (W x L x H)	35" x 80"L x 7"-5"	36" x 80" x 10"	42" x 80" x 10"	48" x 80" x 10"	54" x 80" x 10"	60" x 80" x 10"	42" x 80" x 10"	48" x 80" x 10"
Weight Capacity	500 lbs	500 lbs	600 lbs	750 lbs	1000 lbs	1000 lbs	700 lbs	990 lbs
LPM	8	50	12	12	12	12	1300	1300
Alternating Pressure	✓	--	✓	✓	✓	✓	✓	✓
Low Air Loss (LAL)	--	✓	✓	✓	✓	✓	✓	✓
True LAL	--	--	--	--	--	--	✓	✓
Pillow Function	✓	--	✓	✓	✓	✓	--	--
Quick Cell Disconnects	--	✓	--	--	--	--	✓	✓
Autofirm	--	✓	✓	✓	✓	✓	✓	✓
Anti-microbial Treated	--	✓	--	--	--	--	--	--
Seat Inflate / Fowler	--	--	✓	✓	✓	✓	✓	✓
Low Pressure Alarm	--	Audible & Visual	Audible & Visual	Audible & Visual	Audible & Visual	Audible & Visual	Audible & Visual	Audible & Visual
Power Failure Alarm	--	Audible & Visual	Audible & Visual	Audible & Visual	Audible & Visual	Audible & Visual	Audible & Visual	Audible & Visual
System Failure Alarm	--	--	Audible & Visual	Audible & Visual	Audible & Visual	Audible & Visual	Audible & Visual	Audible & Visual
Auto Return from Autofirm to Therapy	--	30 min	30 min	30 min	30 min	30 min	20 min	20 min
Auto Return from Static to Therapy	--	--	120 min	120 min	120 min	120 min	--	--
Power Outage Protection	✓	✓	✓	✓	✓	✓	✓	✓
Raised Perimeters	--	--	--	--	--	--	--	--
Top Cover Material	4 Way Stretch Dermaplush	Quilted Polyester/Polyurethane	Quilted Nylon/Polyurethane	Quilted Nylon/Polyurethane	Quilted Nylon/Polyurethane	Quilted Nylon/Polyurethane	Stretch Silver Ion Polyurethane	Stretch Silver Ion Polyurethane
Air Cell Count	8	20	20	20	20	20	22	26
Warranty	5 year mattress, 3 year cover, 12 month control unit	24 months	18 months	18 months	18 months	18 months	24 months	24 months
HPCPS	E0277	E0277	E0277	E0277	E0277	E0277	E0277	E0277

LEGEND



E-COMMERCE
Available to order online



PROVIDER
Available in independent DME provider locations



LONG TERM CARE
Available in long term care facilities

STANDARD ★
Entry level systems suitable to meet basic user needs for a good value.

ENHANCED ★★
Meets basic needs plus supplemental features to improve the user's experience.

PREMIUM ★★★
Includes multiple advanced features for the ultimate user experience.