



1 Contoured backstrap with viewing window for improved comfort & visibility



2 Frame design provides added strength & durability with a 350lb weight capacity
PATENTED



3 Integrated slow-down brake controls speed when navigating diverse terrain
PATENTED



4 Braking system reduces tire wear & allows for easy replacement of brake pad & cable
PATENTED



5 Improved tire tread for better stopping ability & large 10" front wheels provide optimal stability and rolling comfort



6 Built-in accessory mount easily attaches a cup, cane or smartphone holder*



7 Built-in brake cable for improved durability and safety



8 Folds side to side with one hand and remains locked for easy transport and storage



9 Padded seat provides added comfort while sitting



10 Fully assembled out of the box



11 Storage bag mounts in multiple locations for easy access



User Guide Video



ITEM #	PRODUCT DESCRIPTION	UPC	UOM
102662RD	Nitro Sprint Rollator, Red	822383018584	1/ea
102662BL	Nitro Sprint Rollator, Blue	822383018591	1/ea
102662BK	Nitro Sprint Rollator, Black	822383018607	1/ea

PRODUCT SPECS

Frame	Aluminum
Width Between Arms	18"
Product Dimensions	24" (W) x 27.5" (D)
Handle Height Min/Max	30" (min) 36.50" (max)
Seat Dimensions	18" (W) x 10" (D) x 21" (H)
Storage Bag Dimensions	14" (L) x 6" (W) x 9" (H)
Wheel Size Front/Rear	10" (front) 8" (rear)
Product Weight	19 lbs.
Color	Red ● Blue ● Black ●
Weight Capacity	350 lbs.
Warranty	Limited Lifetime, Frame**

The sleek Nitro Sprint Rollator is designed to provide comfort, convenience and control at every turn. The rollator sports a unique braking system that reduces tire wear and an integrated slow-down brake for safer maneuverability over diverse terrain. When not in use, the rollator can be folded side to side with one hand and locked in place for easy transport. Other features include multiple storage bag mounting points and a wider backrest with a window for improved comfort and visibility.

BOX DIMENSIONS

Shipping Box Dimension	29" (L) x 9.6" (W) x 33.5" (H)
Shipping Box Weight	23 lb.

*Accessories sold separately.

**Refer to instruction guide or website for warranty details



ITEM #	PRODUCT DESCRIPTION	UPC	UOM
102662RD-H, 102662BL-H, 102662BK-H	Nitro Sprint Rollator, Hemi, Red, Blue, Black	82238018959, 82238018966, 82238018973	1/ea
102662RD-T, 102662BL-T, 102662BK-T	Nitro Sprint Rollator, Tall, Red, Blue, Black	82238018980, 82238018997, 82238019000	1/ea
102662BK-HD	Nitro Sprint Rollator, Heavy Duty, Black	82238019017	1/ea

PRODUCT SPECS	Hemi	Tall	Heavy Duty
Frame	Aluminum	Aluminum	Aluminum
Width Between Arms	18"	18"	23"
Product Dimensions	24" (W) x 27.5" (D)	24" (W) x 27.5" (D)	29" (W) x 28" (D)
Handle Height Min/Max	29.1" (min) 33" (max)	32.0" (min) 39.4" (max)	32.7" (min) 37.8" (max)
Seat Dimensions	18" (W) x 10" (D) x 18" (H)	18" (W) x 10" (D) x 23" (H)	22.5" (W) x 10" (D) x 22" (H)
Storage Bag Dimensions	14" (L) x 6" (W) x 9" (H)	14" (L) x 6" (W) x 9" (H)	18" (L) x 5.5" (W) x 10" (H)
Wheel Size Front/Rear	10" (front) 8" (rear)	10" (front) 8" (rear)	10" (front) 8" (rear)
Product Weight	19 lbs.	19.5 lbs.	22 lbs.
Color	Red ● Blue ● Black ●	Red ● Blue ● Black ●	Black ●
Weight Capacity	350 lbs.	350 lbs.	500 lbs.
Warranty	Limited Lifetime, Frame**	Limited Lifetime, Frame**	Limited Lifetime, Frame**

NITRO™ SPRINT Accessories



ITEM #	PRODUCT DESCRIPTION	UPC	UOM
102662-CUP	Nitro Universal Cup Holder	822383019550	1/ea
102662-CELL	Nitro Universal Cell Phone Holder	822383019567	1/ea
102662-CANE	Nitro Cane Holder	822383019574	1/bx
102662-TRAY	Nitro Sprint Rollator Tray	822383019611	1/ea