

1 **PATENTED** WIDTH EXPANDABLE FRONT AXLES

Expanding front axle adjusts to easily navigate narrow spaces and tight turns or provide additional stability on diverse terrain.

2 Large 10" front wheels provide optimal stability and rolling comfort

3 Durable, and lightweight frame supports users up to 350 lbs.

4 Integrated disc brake technology for increased stopping power that is quick and controlled

5 Built-in brake cable for improved durability and safety

6 Padded leg support provides comfort, stability, adjustability and is easily compatible with left or right leg

7 Height-adjustable handlebar equipped with dual handbrakes and push-button parking brake

8 Removable and portable storage basket that can hold up to 13 lbs

9 Foldable design allows for convenient storage and portability

10 Quick and tool-free assembly



ITEM #	PRODUCT DESCRIPTION	UPC	UOM
791RD	Nitro Knee Walker, Red	822383019086	1/ea
791BL	Nitro Knee Walker, Blue	822383019093	1/ea

BOX DIMENSIONS

Shipping Box Dimension	32.3" (L) x 17" (W) x 17" (H)
Shipping Box Weight	28.2 lbs.

PRODUCT SPECS

Frame	Aluminum/Steel
Product Dimensions (with basket)	16" (W) x 35.5" (D)
Handle Height Min/Max	30" (min) 39" (max)
Knee Pad Height Min/Max	17" (min) 23" (max)
Knee Pad Dimensions	13.5" (L) x 7" (W) x 3" (H)
Wheel Size Front/Rear	10" (front) 7.5" (rear)
Product Weight	23 lb.
Color	Red ● Blue ●
Weight Capacity	350 lbs.
Warranty	Limited Lifetime, Brake Cable**

The Nitro™ Glide Knee Walker offers users recovering from a below-the-knee injury the ability to remain mobile and glide through their daily activities in an all-new sleek design. Offered in two vibrant colors, this patented design features an adjustable-width front axle providing optimized support in both indoor and outdoor terrains. The ergonomic design of the adjustable knee pad provides the necessary comfort and elevation to the injured leg while the adjustable handlebar is equipped with 8 adjustment holes. The lightweight frame safely encloses the brake cable inside to enhance its modern design and the foldable feature makes it convenient for storage and transport.

**Refer to instruction guide or website for warranty details
*Accessories sold separately.