

HURRYCANE®

The All-Terrain Cane

The redesigned HurryCane Freedom Edition features thoughtful upgrades from top to bottom to give you more confidence, stability and comfort than ever before. We've added a larger, ergonomically grooved handle for a more comfortable, secure grip and widened the Steadi-Grip base so your cane stands by your side wherever you go — even on carpet!

The pivoting Whisper-flex design features an increased range of motion, mimicking the motion of your ankle for improved balance and traction. As always, your HurryCane Freedom Edition easily folds and unfolds in seconds for unparalleled convenience. It now also features a stay-tight, snapping wrist strap to secure your cane when folded. Try the new Freedom Edition for yourself to experience the cutting-edge technology and premium craftsmanship that set HurryCane apart and have made it the #1 selling cane in America.



Drive PreserveTech™ products are manufactured with an antimicrobial* agent that helps to minimize the growth of stain and odor-causing bacteria and mold on treated surfaces.

*Antimicrobial properties are built-in to protect the products. These products do not protect users or others against bacteria, viruses, germs or other disease organisms.

WRIST STRAP INCLUDED

COMFORT-FIT, SHOCK-ABSORBING HANDLE

8 CONVENIENT HEIGHT ADJUSTMENTS FROM 31.3"–38.23"

FOLDS & UNFOLDS IN SECONDS

ANODIZED ALUMINUM SHAFT

WHISPERFLEX PIVOTING BASE DESIGN FOR STABILITY & BALANCE

STEADIGRIP BASE FOR ALL-TERRAIN TRACTION

ITEM #

PRODUCT DESCRIPTION

UPC

UOM

| | | | |
|-------------|----------------------------|--------------|------------|
| HCANE-BK-C2 | HurryCane, Freedom, Black | 822383571034 | 1/ea, 6/cv |
| HCANE-BL-C2 | HurryCane, Freedom, Blue | 822383571041 | 1/ea, 6/cv |
| HCANE-RD-C2 | HurryCane, Freedom, Red | 822383571065 | 1/ea, 6/cv |
| HCANE-PR-C2 | HurryCane, Freedom, Purple | 822383571058 | 1/ea, 6/cv |



PRODUCT SPECS

| | |
|--------------------|--|
| Product Dimensions | 5.04" (W) x 4.21" (L) x 31.3" - 38.23" (H) |
| Product Weight | 1.034 lbs. |
| Color | Blue, Purple, Red, Black |
| Weight Capacity | 350 lbs. |
| Warranty | Limited Lifetime** |

BOX SPECS

| | |
|----------------------------|------------------------------------|
| Retail Box Dimensions (EA) | 5.12" (W) x 3.94" (L) x 14.96" (H) |
| Retail Box Weight | 7.408 lbs. |

COLORS AVAILABLE



**Refer to instruction guide or website for warranty details.

HURRYCANE 1.0



THE "OTHER GUYS"

FEATURES & BENEFITS

ERGONOMIC HANDLE WITH COMFORT GRIPS

WRIST STRAP WITH STAY-TIGHT SNAP

FOLDABLE


REPLACEABLE RUBBER FEET WITH
CROSS HATCH TRACTION

"ANKLE MOTION" WITH INCREASED FLEX

HEIGHT ADJUSTABILITY

WEIGHT CAP

SLEEK MODERN DESIGN



31" - 38"

31" - 38"

30" - 37"

350 lbs.

350 lbs.

250 lbs.

EXPLORE OTHER HURRY ITEMS

HurryCycle
Pedal Exerciser



H19E008BK

HurryFlex
Cane Tip



HTIP-BK

HurryCane GO
Offset Cane



HCANE-BK-G2

HurryCane Feet
Replacement Feet



HCFEET-BK