

# NITRO DUET COMBINATION ROLLATOR AND TRANSPORT CHAIR

If you're seeking a safe, convenient aid to improving your daily mobility, a combination rollator/transport chair can be the ideal solution



The Nitro Duet's winning design feature is the ease by which the user can transform the unit from a rollator to a transport chair. This rollator/transport chair is constructed with exceptional detail for powerful control, effortless mobility, eye-catching style and innovative functionality.

- As a rollator, the wheels allow the user to walk with an easy, smooth gait
- As a transport chair, the user can be escorted with ease by a caregiver
- Reliable for traveling indoors or outdoors, it is great for your busy, active lifestyle

## SPECS

Color	Red
Handle Height	33"-37"
Seat Height	20.5"
Unit Dimensions	23.5" (W) x 40" (L)
Seat Dimensions	18" (W) x 14" (D)
Weight	22 lbs.
Weight Limit	300 lbs.
Front Wheels	10"
Warranty (Frame)	Limited Lifetime
Warranty (Brake and Handle Housing)	5 Year Limited
UPC	822383972497
Each Dimensions	28.5" x 10.6" x 34.6"

## FEATURES

- Euro-style design offers an attractive, modern appearance
- Easily transforms from rollator to transport chair
- The seat provides a comfortable, convenient place to rest
- Handle height easily adjusts with unique push button and minimal effort
- Back support height easily adjusts with tool-free thumb screw
- Brake cable located inside frame for added safety, avoiding snags & tears
- Caster fork design enhances turning radius for added security during movement
- 10" front casters allow optimal steering and rolling comfort
- Easy, one-hand folding allows ultra-compact size for minimal storage space
- Cross-brace design allows for side-to-side folding and added stability while folded
- Included clip holds the Nitro Duet closed for added protection while frame is folded and during storage
- Removable zippered storage bag with unique attachment keeps it securely in place when Nitro Duet is open or folded
- Aluminum frame is lightweight for ease of picking up or carrying

For more information or to place an order, please contact customer service.  
877-224-0946 [www.drivemedical.com](http://www.drivemedical.com) or Visit [DriveMedical.com](http://DriveMedical.com) TODAY!

**drive**

**DeVilbiss**  
HEALTHCARE