



1 ERGONOMIC HANDGRIPS & PUSH DOWN BRAKE



2 HANDLE HEIGHT ADJUSTMENT 32.5-37.5"



3 LARGE REMOVABLE STORAGE POUCH



4 LIGHTWEIGHT ALUMINUM FRAME WITH BUILT-IN BRAKE CABLE



5 DUAL FRONT FORK DESIGN & LARGE 9" FRONT WHEEL



ITEM #	PRODUCT DESCRIPTION	UPC	UOM
RTL10266TWHL	PreserveTech™ Nitro 3-Wheel Rollator, Red	822383006314	1/ea

PRODUCT SPECS

Frame	Aluminum
Width Between Arms	16"
Width / Depth	22" (W) x 29" (D)
Handle Height Min/Max	32.50" (min) 37.50" (max)
Wheel size	9.8" (front) 7.5" (back)
Product Weight	13.5 lb.
Color	Red
Weight Capacity	300 lb.
Warranty	Limited Lifetime, Brake and Handle Housing 5 Years**

Featuring a hidden inner brake cable, the functional and modern Nitro 3 is designed with safety in mind, both indoors and outdoors. The Nitro 3 easily navigates narrow spaces and folds flat while in the standing position so it can roll, park and store anywhere. Also included is a large, removable pouch to keep personal items close. The Nitro 3 comes pre-assembled, so you can unbox, unfold and get your customers going!

BOX SPECS

Ship Carton Dimensions	32" (L) x 11" (W) x 34" (H)
Ship Carton Weight	17.9 lb.

**Refer to instruction guide or website for warranty details.