

Splash Defense™ Transfer Bench

PATENTED SHOWER CURTAIN GUARD



DURABLE ALUMIUM FRAME SUPPORTS UP TO 400 LB.

REVERSIBLE TO ACCOMMODATE ANY BATHROOM

SEAT HEIGHT ADJUSTMENTS 17.5–22.5"



Splash Defense™ Transfer Bench is unique and the first of its kind allowing your shower curtain liner to remain inside the bath tub. This product enhancement meets the needs of many by helping eliminate water from spilling out onto the bathroom floor, creating a danger zone for not only the user but also for the care taker. The Splash Defense™ transfer bench provides security for the care taker, while giving an individual their independence back and allowing them to bathe safely and privately again.



HOW TO CHOOSE A SHOWER CHAIR



Scan me



| ITEM # | PRODUCT DESCRIPTION | UPC | UOM |
|-------------|--------------------------------|--------------|------|
| RTL12032KDR | Splash Defense™ Transfer Bench | 822384004142 | 1/ea |

PRODUCT SPECS

| | |
|-------------------------|-------------------------|
| Outside Leg Width/Depth | 33.50" (W) x 23.00" (D) |
| Seat Width / Depth | 28.00" (W) x 16.00" (D) |
| Product Weight | 12.70 lb. |
| Color | White |
| Weight Capacity | 400 lb. |
| Warranty | Limited Lifetime* |

BOX SPECS

| | |
|-----------------------|-------------------------------------|
| Retail Box Dimensions | 9.30" (W) x 21.10" (L) x 17.50" (H) |
| Retail Box Weight | 17 lb. |

*Refer to instruction guide or website for warranty details.