




PreserveTech™ 360° Swivel Bath Chair

Stays Cleaner • Lasts Longer

Compact and customizable, the PreserveTech™ Swivel Bath Chair provides a safe transfer in and out of the tub for those with limited mobility. The seat rotates a full 360° in either direction and automatically locks at 90° intervals to help prevent skin irritation and friction wounds during transfers. The seat height can be adjusted from 15"–22", and the padded armrests can be raised out of the way.

Drive PreserveTech™ products are manufactured with an antimicrobial* agent that helps to minimize the growth of stain and odor-causing bacteria and mold on treated surfaces.

*Antimicrobial properties are built-in to protect the products. These products do not protect users or others against bacteria, viruses, germs or other disease organisms.



Scan me

ITEM #	PRODUCT DESCRIPTION	UPC	UOM
RTL12A001-GR	Swivel Bath Chair	822383005638/ea / 50822383005633/cv	1/ea, 1/cv

PRODUCT SPECS

Product Dimensions	18" (W) x 16" (L) x 15"-21" (H)
Product Weight	8.80 lbs.
Color	Grey
Weight Capacity	350 lbs.
Warranty	Limited Lifetime**

BOX SPECS

Retail Box Dimensions	7.09" (W) x 20.08" (L) x 16.14" (H)
Retail Box Weight	10.30 lbs.

**Refer to instruction guide or website for warranty details.