

2025

Rollators Catalog

Enhancing quality
of life is our
driving force



drive

DeVilbiss
HEALTHCARE

Dependable Independent Mobility

When it comes to independent mobility, Drive DeVilbiss Healthcare is proud to offer an array of dependable, high-quality rollators to meet your needs—with features you'll love, in styles you'll be proud to showcase. Whether your lifestyle calls for light and portable travel-friendly options, or a sturdy and durable selection for everyday use, you will appreciate the thought and care we put into the design and production of every Drive rollator.

Table of Contents

Rollators	3-17
Accessories	18-21



Rollators

Deluxe 3 Wheel Aluminum Rollator with 7.5" Casters

Item #	Color
10289RD	Red
10289BL	Blue

- Comes standard with basket, tray and pouch
- Basket conveniently folds for easy storage
- Ships in ultra-compact carton Comes with frame folded and compact Pull lightweight, aluminum trombone-style frame to lock in place and secure with lock knob
- 7.5" casters for indoor or outdoor use
- Handles easily adjust in height with self-threading knob Special loop lock made of internal aluminum casting operates easily and ensures safety
- Easy, one-hand folding

[Shop Now](#)



Winnie Lite Supreme/Go Lite 3 Wheel Aluminum Rollator

Item #	Color
RTL199	Tan Plaid

- Special loop lock made of internal aluminum casting operates easily and ensures safety
- 7.5" casters with soft-grip tires are ideal for indoor and outdoor use
- Comes standard with carry pouch
- Handles easily adjust in height with self-threading knob
- Easy, one-hand folding

[Shop Now](#)



Rollators

Nitro™ 3-Wheel Rollator

Item #	Color
RTL10266TWHL	Red



- Stylish and modern design features lightweight aluminum frame with built-in brake cable
- Ergonomic hand grips provide comfort and control, push down on handbrakes to engage parking brake feature
- Quickly and easily adjust the height of the handles with the push-button height adjustment
- Large, removable storage pouch
- Easy side to side fold allows the rollator to stay in the standing position, making it easy to roll, park and store
- Dual front fork design enhances turning radius while providing added safety
- Large 9" front wheel allows for optimal steering
- Out of box assembled

[Shop Now](#)



Steel Rollator with 6" Wheels, Knockdown

Item #	Color
R800KD-RD	Red
R800KD-BL	Blue



- Strong, durable steel frame
- Height adjustable arms (31"-36")
- Padded seat allows user to remain comfortable during periods of extended sitting
- Padded backrest adds additional comfort and support
- Convenient storage pouch under seat for easy access to your personal items
- Ergonomic hand grip with easy to use loop-lock brakes provides comfort and control
- 6" wheels with soft grip tires for indoor and outdoor use
- Compact and foldable design
- Tool-free quick and easy assembly
- Weighs just 17 lb

[Shop Now](#)



Rollators

Durable 4 Wheel Rollator with 7.5" Casters



Item #	Color
10257BL-1	Blue
10257RD-1	Red

- Durable steel frame
- Convenient mobility with a comfortable, restful seat
- Pouch under seat for added privacy and security of personal items
- Easy-turn lever for quick height adjustment
- Large 7.5 inch caster wheels are great for indoor and outdoor use
- Deluxe hand brake can be pressed when needed or pushed down to lock rollator into place
- Removable, hinged, padded backrest can be folded up or down
- Easily folds for transport or storage

[Shop Now](#)



Aluminum Rollator with 6" Casters



Item #	Color
R726RD	Red
R726BK	Black
R726BL	Blue

- Lightweight aluminum frame
- Comes with a comfortable, durable padded seat
- Zippered pouch under seat for added privacy and security of personal items
- Removable, hinged, padded backrest can be folded up or down
- 6" casters with soft grip tires ideal for indoor and outdoor use and make the ride smooth and stable
- Ergonomic handles are easy-to-grip and are adjustable

[Shop Now](#)



Rollators

Adjustable Height Rollator with 6" Casters



Item #	Color
RTL10261RD	Red
RTL10261BL	Blue

- Universal height adjustment on frame allows seat height to be adjusted from 18" to 22" in 1" increments
- Comes with a comfortable, durable padded seat
- Removable, hinged, padded backrest can be folded up or down
- 6" casters with soft grip tires ideal for indoor and outdoor use and make the ride smooth and stable
- Ergonomic handles are easy-to-grip and are adjustable in height with self-threading knob

[Shop Now](#)



Aluminum Rollator with 7.5" Casters



Item #	Color
R728RD	Red
R728BL	Blue

- Comes with a comfortable, durable padded seat
- Zippered pouch under seat for added privacy and security of personal items
- Removable, hinged, padded backrest can be folded up or down
- 7.5" casters with soft grip tires ideal for indoor and outdoor use and make the ride smooth and stable
- Ergonomic handles are easy-to-grip and are adjustable

[Shop Now](#)



Rollators

Compact Folding Rollator



Item #	Color
RTL10288BK	Black

- Easy side-to-side fold for daily use and compact fold for storage and travel
- This lightweight rollator only weighs 15.5 lbs for effortless travel and transport
- Strong, durable aluminum frame
- Adjustable handle height from 31"-36"
- Padded backrest provides additional support when in a seated position
- Seat allows you to take a break when needed and remain comfortable
- 8" wheels are ideal for maneuvering both indoors and outdoors
- Removeable storage bag for your personal items
- Cane holder included for added convenience
- Tool-free assembly

[Shop Now](#)



Nitro Aluminum Rollator, 10" Casters



Item #	Color
RTL10266	Red
RTL10266BK	Black

- This rolling walker is cleverly engineered with a comfortable seat, adjustable height handles, and large 10" front casters for seamless rolling.
- Brake cables are located inside the frame to provide a safe experience
- Adjustable contoured backrest
- Cross-brace design for added safety
- Clip hold and side-by-side folding
- Large and convenient storage bag included
- Weighs only 18 lb., easy-to-transport

[Shop Now](#)



Rollators

Nitro Aluminum Rollator, Hemi Height, 10" Casters

Item #
RTL10266BK-H

Color
Black

- Hemi-height ideal for users 5'2" and below
- Brake cable inside frame for added safety
- Handle height easily adjusts with unique push button
- Back support height adjusts with tool-free thumb screw
- Caster fork design enhances turning radius
- Large 10" front casters allow optimal steering and rolling comfort
- Easily folds with one hand to ultra-compact size for storage
- Plastic clip holds the Nitro closed when frame is folded
- Cross-brace design allows for side-to-side folding and added stability
- Lightweight aluminum frame



[Shop Now](#)

Nitro Aluminum Rollator, Tall Height, 10" Casters

Item #
RTL10266-T
RTL10266BK-T

Color
Red
Black

- Tall height ideal for users 6'0" and above
- Brake cable inside frame for added safety
- Handle height easily adjusts with unique push button
- Back support height adjusts with tool-free thumb screw
- Caster fork design enhances turning radius
- Large 10" front casters allow optimal steering and rolling comfort
- Easily folds with one hand for storage
- Plastic clip holds the Nitro closed when frame is folded
- Cross-brace design allows for side-to-side folding and added stability
- Lightweight aluminum frame



[Shop Now](#)

Rollators

Nitro DLX Rollator

Item

RTL10266CH-HS

RTL10266BL-HS

Color

Champagne

Blue



- Handle height easily adjusts with unique push button
- Back support height adjusts with tool-free thumb screw
- Caster fork design enhances turning radius
- 10" front casters allow excellent steering and rolling comfort
- Easily folds with one hand to ultra-compact size
- Cross-brace design allows for side-to-side folding and added stability
- Handsome and removable deluxe zippered storage bag with unique attachment to keep bag securely in place when Nitro is open or folded
- Lightweight frame

[Shop Now](#)



Nitro Sprint Rollator

Item

102662RD

102662BK

102662BL

Color

Red

Black

Blue



- Contoured backstrap for improved comfort and visibility
- The patented frame design provides added strength and durability with a 350lb weight capacity.
- Integrated slow-down brake controls speed when navigating diverse terrain.
- The enhanced braking system reduces tire wear and allows for easy replacement of brake pads and cables
- Built-in accessory mount easily attaches a cup, cane or smartphone holder.
- Built-in brake cable for improved durability and safety
- Large 10" front wheels provide optimal stability'
- Folds side to side with one hand
- A padded seat provides added comfort while sitting
- Storage bag mounts in multiple locations for easy access

[Shop Now](#)



Rollators

Nitro Sprint Rollator with Accessories

Item #	Color
102662RD-ACC	Red
102662BK-ACC	Black
102662BL-ACC	Blue

- Built-in accessory mount easily attaches a cup, cane or smartphone holder. ***ALL 3 ACCESSORIES INCLUDED!**
- Contoured backstrap for improved comfort and visibility
- The patented frame design provides added strength
- Integrated slow-down brake controls speed when navigating diverse terrain.
- The enhanced braking system reduces tire wear and allows for easy replacement of brake pads and cables
- Built-in brake cable for improved durability and safety
- Large 10" front wheels provide optimal stability'
- Folds side to side with one hand
- A padded seat provides added comfort while sitting
- Storage bag mounts in multiple locations for easy access.

[Shop Now](#)



Nitro Sprint Rollator, Tall

Item #	Color
102662RD-T	Red
102662BK-T	Black
102662BL-T	Blue

- Contoured backstrap for improved comfort and visibility
- The patented frame design provides added strength and durability with a 350lb weight capacity.
- Integrated slow-down brake controls speed when navigating diverse terrain.
- The enhanced braking system reduces tire wear and allows for easy replacement of brake pads and cables
- Built-in accessory mount easily attaches a cup, cane or smartphone holder.
- Built-in brake cable for improved durability and safety
- Large 10" front wheels provide optimal stability'
- Folds side to side with one hand
- A padded seat provides added comfort while sitting
- Storage bag mounts in multiple locations for easy access.

[Shop Now](#)



Rollators

Nitro Sprint Rollator, Tall with Accessories

Item #	Color
102662RD-T-ACC	Red
102662BK-T-ACC	Black
102662BL-T-ACC	Blue

- Built-in accessory mount easily attaches a cup, cane or smartphone holder. ***ALL 3 ACCESSORIES INCLUDED!**
- Contoured backstrap for improved comfort and visibility
- The patented frame design provides added strength
- Integrated slow-down brake controls speed when navigating diverse terrain.
- The enhanced braking system reduces tire wear and allows for easy replacement of brake pads and cables
- Built-in brake cable for improved durability and safety
- Large 10" front wheels provide optimal stability'
- Folds side to side with one hand
- A padded seat provides added comfort while sitting
- Storage bag mounts in multiple locations for easy access.

[Shop Now](#)



Nitro Sprint Rollator, Hemi

Item #	Color
102662RD-H	Red
102662BK-H	Black
102662BL-H	Blue

- Contoured backstrap for improved comfort and visibility
- The patented frame design provides added strength and durability with a 350lb weight capacity.
- Integrated slow-down brake controls speed when navigating diverse terrain.
- The enhanced braking system reduces tire wear and allows for easy replacement of brake pads and cables
- Built-in accessory mount easily attaches a cup, cane or smartphone holder.
- Built-in brake cable for improved durability and safety
- Large 10" front wheels provide optimal stability'
- Folds side to side with one hand
- A padded seat provides added comfort while sitting
- Storage bag mounts in multiple locations for easy access.

[Shop Now](#)



Rollators

Nitro Sprint Rollator, Hemi with Accessories



Item #	Color
102662RD-H-ACC	Red
102662BK-H-ACC	Black
102662BL-H-ACC	Blue

- Built-in accessory mount easily attaches a cup, cane or smartphone holder. ***ALL 3 ACCESSORIES INCLUDED!** Contoured backstrap for improved comfort and visibility
- The patented frame design provides added strength
- Integrated slow-down brake controls speed when navigating diverse terrain.
- The enhanced braking system reduces tire wear and allows for easy replacement of brake pads and cables
- Built-in brake cable for improved durability and safety
- Large 10" front wheels provide optimal stability'
- Folds side to side with one hand
- A padded seat provides added comfort while sitting
- Storage bag mounts in multiple locations for easy access.

[Shop Now](#)



AeroWalk™ Rollator

Item #	Color
10277GY	Grey
10277WT	White
10277BL	Blue

Coming Soon!



Rollators

Nitro Elite CF, Carbon Fiber Rollator

Item #

RTL10266CF

Color

Black Carbon Fiber

- Luxurious super lightweight carbon fiber frame
- Frame weighs only 12 lb for powerful performance with swift and effortless mobility
- Opulent and removable zippered storage bag
- Unique bag attachment keeps bag securely in place when rollator is open or folded
- Brake cable inside frame for added safety
- Handle height easily adjusts with unique push button
- Caster fork design enhances turning radius
- Large 10" front casters allow optimal steering and rolling comfort
- Easily fold with one hand to ultra-compact size for storage
- Cross-brace design allows for side-to-side folding
- EZ clip cane holder included

[Shop Now](#)

Aluminum Junior Rollator with 6" Casters

Item #

301PSBN

Color

Blue

301PSRN

Red

- 6" casters with soft-grip are ideal for indoor and outdoor use
- Tool-free removable padded back rest with cover for comfort
- 18" floor to seat height
- Comfortable padded seat
- Handles are adjustable height
- Special loop lock made of internal aluminum casting operates easily and ensures safety
- Easy, one-hand folding

[Shop Now](#)

Rollators

Go-Lite Bariatric Steel Rollator



Item #	Color
10215RD-1	Red
10215BL-1	Blue

- 8" casters with soft-grip tires are ideal for indoor and outdoor use
- Soft padded oversized seat
- Comes with large 20.5" (W) x 4.5" (H) x 9" (D) basket that can be mounted under seat
- Tool-free removable padded backrest for comfort
- Special loop lock made of internal aluminum casting operates easily and ensures safety
- Handles easily adjustable in height with self-threading knob

[Shop Now](#)



Nitro HD Rollator



Item #	Color
10266HD-BK	Black

- Brake cable inside frame for added safety
- Handle height easily adjusts with unique push button
- Unique suspension system for support and comfort
- Back support height easily adjusts tool-free
- Caster fork design enhances turning radius
- Large 10" front casters allow optimal steering
- Easily folds with one hand to ultra-compact size
- Under seat locking mechanism holds the Nitro closed when frame is folded
- Cross-brace design allows for side-to-side folding and added stability
- Handsome and removable zippered storage bag with unique attachment to keep bag securely in place when Nitro is open or folded
- Lightweight, aluminum frame

[Shop Now](#)



Rollators

Nitro® Sprint HD Rollator

Item #

102662BK-HD

Color

Black



- Contoured backstrap with viewing window
- The patented frame design provides added strength and durability with a 500 lb weight capacity
- Integrated slow-down brake controls speed
- Enhanced braking system reduces tire wear and allows for easy replacement of brake pads and cables
- Improved tire tread for better stopping ability
- Built-in accessory mount easily attaches a cup, cane or smartphone holder
- Built-in brake cable for improved durability and safety
- Large 10" front wheels provide optimal stability and rolling comfort
- Folds side to side with one hand and remains locked for easy transport and storage
- Storage bag mounts in multiple locations for easy access

[Shop Now](#)

Duet Rollator/Transport Chair with 6" Casters

Item #

795BU

795BK

795B

Color

Burgundy

Black

Blue



- Combines the features of a rollator and transport chair in one unit
- An individual can ambulate independently or be pushed safely by a caregiver
- The tool-free two-position contoured padded backrest is reversible
- Comes with padded seat with carry pouch
- The folding frame has a lock-out feature to ensure the unit stays open while in use and ensures the unit will stay folded during transport
- 8" casters are ideal for indoor and outdoor use
- Handles are adjustable in height

[Shop Now](#)

Rollators

Nitro Duet Rollator and Transport Chair

Item #

RTL10266DT

Color

Red

- Combines the features of a rollator and transport chair into one unit. The user can ambulate independently or be pushed safely by a caregiver
- The tool-free two-position contoured backrest is reversible depending on whether you are using the Nitro Duet as a rollator or transport chair
- Back support height easily adjusts with tool-free buttons
- Caster fork design and large 10" front casters enhances turning radius and rolling comfort
- Cross-brace design allows for side-to-side folding and added stability
- Brake cables routed inside the frame add safety and convenience

[Shop Now](#)

Elevate Upright Walker

Item #

UPWKRBL10

Color

Blue

UPWKRGY10

Grey

- Natural line of sight than traditional mobility aids.
- Height-adjustable upholstered forearm pads keep arms supported and level for comfortable use, while large, ten-inch wheels make it perfect for indoor and outdoor use
- Padded armrests and upright handle adjust to 44"
- Durable and comfortable padded seat
- Simple one-handed, side-to-side folding
- Includes removable zipper storage bag with strap
- Dual front fork for optimal stability and
- Loop lock brakes for added safety and control
- Adjustable length dual handbrake system
- Large 10" multi-terrain front wheels for
- Height adjustable backrest

[Shop Now](#)

Nitro O2 Tank Holder

Item # 10266-4OX

- Oxygen tank holder for exclusive use with Drive's Nitro rollator (Not compatible with Nitro Sprint, Nitro Glide, Nitro Duet, or Nitro 3 Wheel)
- Accommodates the most commonly used tank sizes (ML6/A tanks, M4/A tanks, M9/C tanks, M6/B tanks, M15/D tanks)
- Conveniently attaches to the Nitro's front end without affecting any of the Nitro's functionality

[Shop Now](#)



Nitro O2 Tank Holder

Item # 10266HD-4OX

- Oxygen tank holder for exclusive use with Drive's Nitro Heavy Duty rollator (Not compatible with Nitro Sprint, Nitro Glide, Nitro Duet, or Nitro 3 Wheel)
- Accommodates the most commonly used tank sizes (ML6/A tanks, M4/A tanks, M9/C tanks, M6/B tanks, M15/D tanks)
- Conveniently attaches to the Nitro's front end without affecting any of the Nitro's functionality

[Shop Now](#)



Nitro Slow Down Brake

Item # 10266-BRAKE ENH

- Compatible with Nitro Rollator

[Shop Now](#)



Slow-Down Brakes for Side-Fold

Item # 730-942

- Fits Airgo eXcursion and Fusion models
- Not compatible with eXcursion6 due to 6" wheels
- Retro-compatible with all revisions, but works best on new low-profile wheels
- Can apply pressure to both sides, or one side to correct gait

[Shop Now](#)



Nitro Cup Holder

Item # 10266-CH

- Securely holds cup or bottle
- Fits all Nitro models
- Tool free installation

[Shop Now](#)

Nitro Sprint and Glide Cup Holder

Item # 102662-CUP

- The Nitro Sprint and Glide Cup Holder is designed to fit on all Nitro Sprint Rollators and Nitro Glide Knee Walkers with built-in accessory mount to provide convenience at every turn
- Fits on all Nitro Sprint Rollators and Nitro Glide Knee Walkers with built-in accessory mount
- Securely holds cups, water bottles, and travel mugs all at an arms reach

[Shop Now](#)



Accessories

Nitro Sprint Rollator Tray

Item # 102662-TRAY

- The Nitro Sprint Tray is designed to fit on the Nitro Sprint Rollator to provide convenience at every turn
- Fits on all Nitro Sprint Rollators
- Tray sits securely on top of the rollator seat
- Molded cup holder provides added stability for beverages

[Shop Now](#)



Nitro Sprint and Glide Cell Phone Holder

Item # 102662-CELL

- The Nitro Sprint and Glide Cell Phone Holder is designed to fit on all Nitro Sprint Rollators and Nitro Glide Knee Walkers with built-in accessory mount to provide convenience at every turn
- Fits on all Nitro Sprint Rollators and Nitro Glide Knee Walkers with built-in accessory mount
- Securely holds most smart phones with ability to view and use cell phone while in cell phone holder
- Locking mechanism ensures cell phone remains in its holder

[Shop Now](#)



Nitro Sprint Cane Holder

Item # 102662-CANE

- The Nitro Cane Holder is designed to fit on the Nitro Sprint Rollator with built-in accessory mount to provide convenience at every turn
- Fits on all Nitro Sprint Rollators with built-in accessory mount
- Securely holds most standard canes on left side of rollator
- Easily mount cane or remove cane from rollator for seamless transition from cane to rollator

[Shop Now](#)



Nitro Cane Holder

Item # 892

- Fits all Nitro products
- Holds most manufactures canes

[Shop Now](#)





CORPORATE OFFICE 🏠 99 Seaview Blvd., Port Washington, NY 11050 📞 877-374-8763 🌐 www.drivemedical.com
CANADA OFFICE 🏠 256 Aviva Park Drive, Vaughan, Ontario Canada L4L 9C7 📞 866-752-1363

©2025 Medical Depot, Inc. dba Drive DeVilbiss Healthcare, 99 Seaview Boulevard, Port Washington, NY 11050. All trademarks used in association with the sale of products of Drive DeVilbiss are trademarks owned by Medical Depot, Inc. All other trademarks, trade names, service marks and logos referenced herein belong to their respective companies. Refer to product instruction guide for warranty details.